**STERİL PRİMER KALÇA PROTEZLERİ / STERILE PRIMARY  
HIP PROSTHESIS**

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*Total hip arthroplasty is designed to improve patient mobility and support components. By replacing the hip joint, it usually relieves pain associated with arthritis and other hip diseases or injuries and increases patient mobility.*

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*Stem and acetabular component types were designed according to bone deformation.*

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*The hip joint is a joint that plays an important role in the realization of many movements such as standing, sitting, running, squatting.*

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*Tıpmad Primary Total Hip Prosthesis are the designs of metal and polyethylene components on both sides of the hip joint (head and pelvis in the pelvis).*

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*In order to manage the primary operations within Tıpmad, different stem designs have been formed due to different patient anatomy and deformation types.*

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*The product lists of the systems with the model names above are available in the Design Files. (Annex 1)*

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*TCK Plus is the most up-to-date product in stem designs. It is designed to provide cemented and uncemented use.*

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*The stem produced from titanium and steel materials is used*

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*The TCK Plus hip prosthesis stem unit has a uniform design and is designed to treat the pain and the indications that prevent daily activities of patients.*

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*TCK 10 Hip prosthesis is designed for cemented use. There are 5 different types of coating by combining sand blasting, Hydroxapatite Coating, Plasma Spray Coating, Porous Coating and coating types on the hydroxapatite coating.*

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*In order to provide appropriate treatment for the patient's anatomy, it is designed in 14 different sizes. Thanks to the 12/14 tapered connection, the hip hip prosthesis is used in combination with acetabular components.*

### TMS 30 Kalça Protezi, TMS 30 Kalça Protezi – Yakalıklık / TMS 30 Hip Prosthesis, TMS 30 Hip Prosthesis Collared

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*The shaft provides wedge-shaped rotational stability with 3-*

*dimensional tapered round edges. The femur channel has a structure suitable for self-centering and creating an effective cement cover. The proximal portion of the shaft has a conical structure 12/14, which is attached to the necked orthopedic femoral heads which accept optimal proximal positioning. It is used as a component in distal centering and cemented use. There is also a model for supporting the Kalkar region. The product is designed in 5 different stem sizes in order to adapt to the differences caused by patient anatomy.*

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*Tipmed Straight Stem Hip Prosthesis is a monolithic hip prosthesis design that can replace the proximal femur region. It has a slim design to protect bone reserve in patients with good bone quality. 9 different stem sizes are used with combined hip prosthesis acetabular components.*

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*Tipmed Uncemented Straight Stem Hip Prosthesis is designed for cementless, press-fit applications. It has a medial curvature that is shaped in such a way that it is proportional to its size in order to provide optimum adaptation to the pressure applied to the patient's anatomy according to the anatomy of the patient.*

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*It is produced from 10 different titanium and titanium materials. Thanks to its 12/14 conical structure, Hekimmed hip prosthesis works in combination with acetabular*

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Tipmed S2 Monoblock Hip Prosthesis is designed to have different varieties in the neck region to support the weak spongiosis in the medulla canal. The prosthesis is designed to be used for femoral fractures of patients undergoing primary hip surgery. It is intended to support the calcareous region with bone loss in the medial by filling the metaphyseal region. It is designed for completely cementless application.

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*Because of the bone loss in the calcar area, it is preferred to use it when the bone quality is poor or the deformation occurs in the limestone area. Serkilaj wire is designed to be used in combination with. It has 5 different stem lengths due to patient anatomy and deformation differences. 5 different stem sizes are designed to include 4 different calcareous zone dimensions. There are 20 different size sizes in total.*

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*Head for Hemi Prosthesis are intended for cemented and uncemented use, such as hemi hip replacements. It is designed for use in femoral head / neck fractures, femoral head / neck aseptic necrosis, osteo, rheumatoid and post-traumatic arthritis of the hip with minimum acetabular.*

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*Unipolar Head for Hemi Prosthesis are for cemented and uncemented use, such as hemi hip replacements. It is designed for use in femoral head / neck fractures, femoral head / neck aseptic necrosis, osteo, rheumatoid and post-traumatic arthritis of the hip with minimum acetabular.*

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*Modular head for total prostheses for use in uncemented applications. Osteoarthritis including avascular necrosis, rheumatoid arthritis, noninflammatory joint diseases, functional deformity correction, nonunion treatments, femoral neck fractures, proximal femur trochanteric fractures are used. The head, which is placed on the stoma, provides part of the prosthesis to the prosthesis.*

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*The Polyethylene Cemented Acetabular Cup is placed into the degenerated anatomic acetabulum.*

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*Bipolar Cup is designed for uncemented applications. The insert locks in the femoral head and can move freely within the metallic shell. This double mobility provides the advantage of low friction with the head. This design has low dislocation rates in both primary hip replacement and recurrent revision cases. The optimized cup design prevents soft tissue and femoral neck from hitting each other.*

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*Acetabular Titanium metal Shell, part of the total hip replacement, is used for structuring severe hip joint injuries caused by osteoarthritis, rheumatoid arthritis, and trauma arthritis. Unopened screw holes can be closed with screw plug. In this case there is no contact between metal and polyethylene.*

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*UHMPWPE Insert is inserted into metal acetabular cup. The liner can be removed with two slots. Polyethylene insert is produced as 0° as in standard and 10°- 20° to avoid dislocation.*

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*Special bone screw with sphere-shaped head and Ø 6.5 mm spongios groove is used for fixation. The screw is in good contact with the place where the insertion angle is different, even if it is different. This prevents wear and provides screw stability.*

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*Acetabular Reinforcement Cage is used for revision of the bony defects or protrusion. Acetabular Reinforcement Cage bone grafts ensure that the combination of. The acetabular ring is implanted with cement. The dorsal part can provide better rotation. Holes on the edge can provide better fixation in the pelvic wall. The distal nose is inserted through the container in the ischium. The container must be cemented in the correct position.*

#### **Asetabular Halka, Anatomik Halka, Kancalı Halka / Acetabular Ring, Anatomic Ring, Hook Ring**

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*In Total Hip Replacements; osteoporosis, rheumatoid arthritis, acetabular cysts, concurrent bone grafting, post arthrodesis and socket revisions are used in the acetabular ring. Acetabular It is recommended to use the rings in combination with polyethylene containers with low wall thickness. Acetabular Ring use is recommended in cases of severe bone loss and protrusion revision. Acetabular Ring is fixed with 6 spongios screws.*

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*Distal Centralizer is manufactured from UHMWPE. It is intended to be used with Stem. The centerer has 4 winged exterior views for better centering. The design of the wing helps to provide better cortical diaphyseal and stem trajectory.*

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*The plug is designed as a femoral cement restrictor and limits the working area. Stops acrylic cement process when printing starts*

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*Products do not have any accessories.*

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*The products are only component with Tıpmed products.*

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*It is designed for use by doctors who have completed their specialist training in Orthopedics and Traumatology.*

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*It is suitable for patients who have completed their skeletal development.*

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*In order to minimize the side effects in the patient population, the following criteria should be taken into consideration in patient selection.*

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*Active infection in the organism; skeletal immaturity; clinical conditions presented by patients who may compromise the stability of the implant (muscle atrophy, neuromuscular iris and vascular insufficiency); women of childbearing age (age 45 years); kidney failure; severe obesity (body mass index greater than 40 kg / m<sup>2</sup>); patients with immunological depression (acquired immunodeficiency syndrome, use of corticoid at immunosuppressive doses or use of other immunosuppressive agents); and precision known to the metal. Radiographic parameters: osteoporosis;*

*osteonecrosis affecting more than 50% of the femoral head; Multiple cysts or cysts larger than 1 cm; Bone mineral density less than 0.65 g / cm<sup>2</sup>*

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Implants are completely implanted in the patient. Surgical invasive contact time is more than 30 days. (Long Term)

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*The physical conditions of the patient, which prevents adequate implant support or prevents the use of the appropriate size of the implant, and the inadequate bone quality resulting from these conditions. For example: cancer or congenital dislocation, upper femoral metabolic bone diseases*

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*inadequate fixation of the implant to the bone*

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Ürünler MR ortamında değerlendirilmemiştir. MR ortamında ısınma veya göçme için test edilmemiştir. *Tıpmed products are combined only with their components and surgical instruments. The products are matched according to the information you have on the product. For this reason, it is inconvenient for the appliance to be combined with different products or different hand tools. The excessive force to surgical instruments can cause to break. It is recommended that surgical instruments of Tıpmed Primary Hip Prosthesis be checked before each use.*

*The products were not evaluated on MR. It has not been tested for heating or dislocation in the MR environment*

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*Mobilization of the patient with an assistive walking device should be made as early as possible. Postoperative rehabilitation should be completed in accordance with each hospital's own practices.*

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The first component of the prosthesis is removed from the bone and then taken out of the patient's body with the help of the extractor.

The excreted product must be disposed of in accordance with national legislation and hospital procedures.

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*TIPMED implants are only accepted by the hospital or surgeon when the original package and label are delivered without tearing.*

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*The implants are exposed to a minimum of 25 kgy gamma irradiation. UHMWPE implants are sterilized by Ethylene Oxide. Surgical Instruments are non-sterile. Must be sterilized before use. Before the operation, the sterile package should be checked from the angle of the hole, tears etc.*

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