**STERİL REVİZYON KALÇA PROTEZLERİ / STERILE
REVISION HIP PROSTHESIS**

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Total hip arthroplasty is designed to improve patient mobility and support components. By replacing the hip joint, it usually relieves pain associated with arthritis and other hip diseases or injuries and increases patient mobility.

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Stem and acetabular component types were designed according to bone deformation.

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The hip joint is a joint that plays an important role in the realization of many movements such as standing, sitting, running, squatting.

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Tıpmad Revision Hip Prosthesis are the designs of metal and polyethylene components on both sides of the hip joint (head and pelvis in the pelvis).

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In order to manage the revision operations within Tıpmad, different stem designs have been formed due to different patient anatomy and deformation types.

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The product lists of the systems with the model names above are available in the Design Files. (Annex 1)

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TCK Plus Hip Prosthesis combines successful designs with cutting edge technology. TCK Plus is intended for use in noninflammatory degenerative joint disease including osteoarthritis and avascular necrosis, rheumatoid arthritis, correction of functional deformity, and revision procedures where other treatments have failed. TCK Plus is used for femoral neck fractures. TCK Plus can be used both cemented and cementless (coating). The proximal part of the stem can be coated with titanium plasma spray, porous, HA and HA + titanium plasma spray (Double coated). For revision patients, TCK Plus has both long and X long sizes. TCK Plus

instrumentation has been designed to make surgical procedures efficient.

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TCK-10 Hip Prosthesis is intended for uncemented use in the primary application as treatment of the followings; patient conditions of noninflammatory degenerative joint disease, avascular necrosis, osteoarthritis and inflammatorily joint disease, rheumatoid arthritis. TCK-10 Hip Prosthesis is designed to be durable for revision hip arthroplasty when there is an intact metaphyseal-diaphyseal junction for adequate press-fit stability at surgery. The interlock grit- blast all over on the stem surface provides increased fixation at the prosthesis/cement interface compared with that of traditional smooth surface implants. Optional Ti. Plasma spray (PL), HA coated or porous coated stems uses the osteo-conductive qualities of this efficient biomaterial or porous which produce the mineral fraction of natural bone. For revision patients, TCK-10 has both long and X long sizes. TCK-10 instrumentation has been designed to make surgical procedures efficient.

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Head for Hemi Prosthesis are intended for cemented and uncemented use, such as hemi hip replacements. It is designed for use in femoral head / neck fractures, femoral head / neck aseptic necrosis, osteo, rheumatoid and post-traumatic arthritis of the hip with minimum acetabular.

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Unipolar Head for Hemi Prosthesis are for cemented and uncemented use, such as hemi hip replacements. It is designed for use in femoral head / neck fractures, femoral head / neck aseptic necrosis, osteo, rheumatoid and post-traumatic arthritis of the hip with minimum acetabular.

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Modular head for total prostheses for use in uncemented applications. Osteoarthritis including avascular necrosis, rheumatoid arthritis, noninflammatory joint diseases, functional deformity correction, nonunion treatments, femoral neck fractures, proximal femur trochanteric fractures are used. The head, which is placed on the stoma, provides part of the prosthesis to the prosthesis.

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The Polyethylene Cemented Acetabular Cup is placed into the degenerated anatomical acetabulum.

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Bipolar Cup is designed for uncemented applications. The insert locks in the femoral head and can move freely within the metallic shell. This double mobility provides the advantage of low friction with the head. This design has low dislocation rates in both primary hip replacement and recurrent revision cases. The optimized cup design prevents soft tissue and femoral neck from hitting each other.

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Acetabular Titanium metal Shell, part of the total hip replacement, is used for structuring severe hip joint injuries caused by osteoarthritis, rheumatoid arthritis, and trauma arthritis. Unopened screw holes can be closed with screw plug. In this case there is no contact between metal and polyethylene.

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UHMPWPE Insert is inserted into metal acetabular cup. The liner can be removed with two slots. Polyethylene insert is produced as 0° as in standard and 10°- 20° to avoid dislocation.

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Special bone screw with sphere-shaped head and Ø 6.5 mm spongios groove is used for fixation. The screw is in good contact with the place where the insertion angle is different, even if it is different. This prevents wear and provides screw stability.

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Acetabular Reinforcement Cage is used for revision of the bony defects or protrusion. Acetabular Reinforcement Cage bone grafts ensure that the combination of. The acetabular ring is implanted with cement. The dorsal part can provide better rotation. Holes on the edge can provide better fixation in the pelvic wall. The distal nose is inserted through the container in the ischium. The container must be cemented in the correct position.

Asetabular Halka, Anatomik Halka, Kancalı Halka / Acetabular Ring, Anatomic Ring, Hook Ring

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In Total Hip Replacements; osteoporosis, rheumatoid arthritis, acetabular cysts, concurrent bone grafting, post arthrodesis and socket revisions are used in the acetabular ring. Acetabular It is recommended to use the rings in combination with polyethylene containers with low wall thickness. Acetabular Ring use is recommended in cases of severe bone loss and protrusion revision. Acetabular Ring is fixed with 6 spongios screws.

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Distal Centralizer is manufactured from UHMWPE. It is intended to be used with Stem. The centerer has 4 winged exterior views for better centering. The design of the wing helps to provide better cortical diaphyseal and stem trajectory.

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The plug is designed as a femoral cement restrictor and limits the working area. Stops acrylic cement process when printing starts

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Products do not have any accessories.

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The products are only component with Tıpmad products.

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The products are verified by quality control of combinations with all components. Other companies do not have combinations with their components.

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It is designed for use by doctors who have completed their specialist training in Orthopedics and Traumatology.

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vascular diseases that compromise the relevant limb

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joint area. Before finalizing cleaning is insufficient (bone cement waste, surgical removal of residues) can result in excessive wear

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inadequate fixation of the implant to the bone

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Tıpmed products are combined only with their components and surgical instruments. The products are matched according to the information you have on the product. For this reason, it is inconvenient for the appliance to be combined with different products or different hand tools. The excessive force to surgical instruments can cause to break. It is recommended that

surgical instruments of Tıpmed Revision Hip Prosthesis be checked before each use.

The products were not evaluated on MR. It has not been tested for heating or dislocation in the MR environment

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Mobilization of the patient with an assistive walking device should be made as early as possible. Postoperative rehabilitation should be completed in accordance with each hospital's own practices.

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The first component of the prosthesis is removed from the bone and then taken out of the patient's body with the help of the extractor.

The excreted product must be disposed of in accordance with national legislation and hospital procedures.

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TIPMED implants are only accepted by the hospital or surgeon when the original package and label are delivered without tearing.

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The implants are exposed to a minimum of 25 kgy gamma irradiation. UHMWPE implants are sterilized by Ethylene Oxide. Surgical Instruments are non-sterile. Must be sterilized before use. Before the operation, the sterile package should be checked from the angle of the hole, tears etc.

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