**STERİL TPM REVİZYON DİZ PROTEZLERİ / STERILE TPM
REVISION KNEE PROSTHESIS**

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Tıpm TPM8 Revision Knee Prosthesis is designed as a revision of the TPM8 Total Knee Prosthesis system. It consists of Femoral Komponent, Tibia Metal Back, and components (Augment, Offset Adapter, Stem and screws), Insert and Patella components.

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Knee Revision surgery due to the damaged structure of the patients due to the femoral component is designed as a Cruiate Retaining

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The Femoral Component is made of CoCrMo material (ASTM F75 / ISO 5832-4). There are 2 different anatomy designs. The component has been designed as 8 different to meet the differences due to bone structure. Through the cement channels, it supports cemented use. Working surfaces have a design that minimizes wear.

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Similar to the TPM8 Total Knee Prosthesis, it contains 8 different types and it is different from the measurement structure used in the primary operation. The main reason for this difference is that poor bone reserve of patient for revision surgery

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ISO 5832-3 material was designed with the help of the stem in the revision operations due to poor bone quality instability situation has been minimized.

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Tibia Metal Back is made of CoCrMo material (ASTM F75 / ISO 5832-4). A cement channel is designed for cemented uses. The patient is designed to provide you with an anatomical design suitable for bone structure.

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The material of Insert component is subject to the standard

(ASTM F638 / ISO 5834-2). In order to tolerate differences due to patient anatomy, 2 different types of insert have been designed. Of these, the Plus group has restricted the valgus-varus movements due to the damaged structures of the patient.

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The Patella component is a component that is preferred to use or not to use according to the doctor's determination in patients with patellar injury. There are 2 different design, 4 different kinds of you are available. It is subject to the raw material of the Patella component produced from polyethylene material (ASTM F638 / ISO 5834-2).

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Revision is the wedge-shaped implants used to support deformational bone loss in the tibia and / or femur region and to fill the gaps that occur. It is raw material is ISO 5832-3 It is used because of poor bone reserve of the patients who need knee revision surgery and because of the high incision of the bone in revision surgery.

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It is fixation augment fixation screw to femoral or tibia component. It is raw material is ISO 5832-3

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The eccentric shaft used for axial deflection and displacement of the femoral and tibial components in the tibial and femoral region. It is raw material is ISO 5832-3

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It is fixation offset adapter fixation screw to femoral or tibia component. It is raw material is ISO 5832-3

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A summary of the properties of the components is as follows:

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Products do not have any accessories.

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The products are only component with Tipmed products.

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The products are verified by quality control of combinations with all components. Other companies do not have combinations with their components.

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- İmplant kaybı veya aşınma / Implant loosening or wear
- Enfeksiyon / Infection
- Instabilite / Instability
- Sertlik / Stiffness
- Protezde kırılma / Fracture of prosthesis
- Operasyon yapılan dizde artan ağrı / Increased pain in the affected knee

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before and during incorrect operation may lead to cracking, abrasion, stress fracture and / or excessive wear

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Yan Etkiler / Possible Adverse Effects

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Tıpmed products are combined only with their components and surgical instruments. The products are matched according to the information you have on the product. For this reason, it is inconvenient for the appliance to be combined with different products or different hand tools. The excessive force to surgical instruments can cause to break. It is recommended that surgical instruments of Tıpmed TPM 8 Revision Knee Prosthesis be checked before each use.

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Recovery after revision knee surgery, is slightly slower than the primary knee surgery. It is recommended to give physical therapy support to harvest and to perform walking exercises with walker or walking stick in the first period after the operation.

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Survival of Product: As a result of the literature studies, the product has been targeted to perform on the patient for 10 years.

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TIPMED implants are only accepted by the hospital or surgeon when the original package and label are delivered without tearing.

Sterilizasyon ve Tekrar Sterilizasyon & Bertaraf / Sterilisation & Resterilisation & Disposal

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Most implants are delivered to the market with sterile and protective packaging. The sterilization method is indicated on the package label.

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