*Regain  
your  
mobility*

**STERİL TPM-8  
TOTAL DİZ PROTEZLERİ**

**STERILE TPM-8  
TOTAL KNEE PROSTHESIS**

**STERİLTPM-8 TOTAL DİZ PROTEZLERİ / STERILE TPM-8  
TOTAL KNEE PROSTHESIS**

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## TIPMED KNEE PROSTHESIS (For Sterile Products)

### Description

Tipmed produces a variety of knee joint prostheses designed to be cemented or uncemented. Knee joint prosthesis components include femoral, tibial, and patellar components. The components are available in various designs and sizes for primer applications. Some special components are available. Femoral components, tibial metal back, insert for cruciate retaining and posterior stabilized.

### Material

Femoral Components	CoCrMo Alloyed
Tibial Metal Back	CoCrMo Alloyed
Tibial Metal Back Stem	Titanium Ti6AL4V Alloyed
Patellar Component	UHMWPE/High Cross Bağlı
UHMWPE (ISO 5832-4)	
Posterior Stabilized Component	CoCrMo Alloyed
Insert Holder	Titanium Alloyed
Screws for Posterior Stabilized or Insert Holder (ISO 5832-3)	Titanium Alloyed

### Indications

1. Painful or crippled knee joints in which one or more compartments are affected, resulting from osteoarthritis, rheumatoid arthritis or traumatic arthritis
2. Correction of varus, valgus or posttraumatic deformity
3. Failed osteotomy, arthrodesis or correction or revision of the previous failed joint replacement procedure.
4. Patients with severe knee pain
5. Femur Avascular Necrosis
6. Moderate Valgus, Varus or Stretching Deformations

Note: Patella with deep bowl appearance is used to increase painful movement. And it provides anterior — posterior stability.

### Contraindications

Absolute contraindications include: infection, sepsis or osteomyelitis.

Relative contraindications are:

1. A patient who does not cooperate or a patient with a neurological disorder that does not comply with the instructions.

3. Metabolic disorders that can damage bone formation
4. Osteomalasia
5. Distant infection that can spread to the implant site focus
6. Rapid joint destruction, determined bone loss or marked bone resorption on x-ray.
7. Vascular insufficiency, mucous atrophy or neuromuscular disease
8. The soft tissue surrounding the knee is missing or inadequate
9. Systemic infection
10. Affected joint
11. Skeletal immature (immaturity)
12. Severe imbalance in the joint

### Warnings

Incorrect selection, placement, positioning, alignment, and fixation of the implant components can lead to unexpected stress situations that can lead to shortened lifetime of the knee prosthetic parts. Misalignment of components or incorrect implantation can cause excessive wear and / or malfunction of the implant or procedure. Insufficient cleaning (removal of surgical debris) before closing can cause excessive wear. Use clean gloves when handling implants. Laboratory tests have shown that implants exposed to body fluids, surgical waste, or oily tissue show less adherence to cement than implants held with clean gloves. Improper use or damage (scratches, nicks) of implants before and during the operation can lead to cracked choroxia, wear, stress fractures and / or wear. Do not make changes to the implants. Do not reuse (used for only one patient.) Do not use damaged components or knee components from different companies. The surgeon must be familiar with the implants and devices before performing the operation.

Misalignment of components or soft tissue imbalance can create excessive pressure on the components, causing excessive wear on patellar or tibial bed articulation surfaces.

Revision operation can be done to prevent damage to the components. It should be ensured that all its components in bone cement are fully supported to prevent stress concentrations that can lead to failure of the procedure. In order to minimize the wear of the implant moving surfaces, it is very important to completely remove and remove bone cement residues, metallic and other

the surgeon's responsibility to determine whether there is sufficient fixation and stability. The Tipmed Knee Prosthesis provides surgeons with a tool to reduce pain and restore function for patients. These devices are generally successful in achieving these goals. However, these devices cannot be expected to reach a normal healthy bone and tissue activity level and load. Practices accepted in post-operative care are important.

The patient's negligence about following post-operative care instructions, including rehabilitation, may jeopardize the success of the procedure. The patient needs to be protected for full load bearing until the patient has sufficient fixation and healing. You should be informed about the limit of this reconstruction. Excessive activity, trauma and excessive weight can cause relaxation and / or wear, leading to an early failure of the implant. Loosening of the implants can lead to further production of wear particles. It can speed up damage to the bone and complicate successful revision surgery. The patient should be warned and raised to comply with the doctor's instructions, including general surgical risks and possible side effects and follow-up visits.

### Precautions

The safety and efficacy of this device is contraindicated for patients who need rheumatoid arthritis, polyarthrosis, pseudogout, or revision. A special tool has been designed for joint replacement systems to assist error-free implantation of knee prosthetic components.

Using tools and implant components from other systems can lead to misplacement, extent of error, excessive wear, and malfunction of the device. Devices that are exposed to excessive use and excessive power may occur in the fracture. Surgical devices should be used for their intended purpose.

Tipmed should be warned and raised to follow the directions of the doctor, including general surgical risks and possible side effects listed, as well as follow-up visits.

Reuse may cause infection problems or cause deformation and breakage due to fatigue of metal that cannot be detected visually. Do not reuse implants. Even if the implant appears undamaged, the stress it has previously exposed may have caused defects that could reduce the implant's lifetime.

The products have not been evaluated in MR environment, there is no information about MR compatibility.

the specialist physician who will apply the revision to the patient should know the surgical operation technique and inform the patient about the contraindications. Disposal of the extracted Implant must be ensured in accordance with biological waste regulations.

- If there is an unexpected change in product performance after the operation, the product should be removed with the appropriate surgical protocol and the patient should be transmitted to the company with x-ray films.

- Disposed product should be disposed of according to national legislation and hospital procedures.

#### **Possible Adverse Effects**

1. Material sensitivity reactions. Foreign material implantation can cause histological reactions in tissues, including macrophages and fibroblasts of various sizes.

The clinical significance of this effect is uncertain. Similar changes may occur during the healing process or earlier. The surrounding tissue or fluid may contain discoloration due to metallic and polyethylene components of joint implants and increased particle wear. Waste from abrasion

It has been reported that it may cause a cellular response resulting in osteolysis or may be the result of loosening of the implant.

2. Infection and allergic reaction early or late post operation.

3. Poor bone stock from osteoporosis, previous bone damage caused by operations, bone bone perforation or fracture may occur during operation due to resorption or wearing the prosthesis.

4. Fixation loss, trauma, misalignment bone due to excessive activity

5. Periarticular calcification or ossification, with or without an impediment to joint mobility.

6. Insufficient range of motion due to incorrect selection or positioning of components

7. Undesirable shortening of the limb

8. Dislocation or subluxation caused by inadequate fixing or incorrect positioning. Muscle, fibrous tissue laxity may contribute to these conditions.

9. Stress fracture may occur in the component as a result of loss of fixation, forceful activity, misalignment, trauma, non-union, or excessive weight.

10. Crack corrosion and friction can occur in units between components.

11. Wear and / or deformation of moving surfaces

12. Valgus-varus deformity

13. Knee in patients with severe flexion and valgus deformity

Following arthroplasty, temporary peroneal palsy has been reported secondary to surgical manipulation and increased joint motion.

14. Patellar tendon rupture and ligament laxity

15. Intraoperative or postoperative bone fracture and / or postoperative pain.

#### **Packaging & Labelling**

TIPMED implants are accepted only when received by the hospital or surgeon in their original packaging and without tearing the label.

#### **METAL COMPONENTS**

Metal components (hand tools) are provided non-sterile. The device must be sterilized before use. Steam sterilization is recommended at 121C, 1 bar and 30 minutes.

sterilized. If necessary, a steam autoclave can be applied with the appropriate protective pouch. Protect medical equipment from metal or other hard objects, which can damage the product and especially interworking parts.

#### **Sterilisation/Resterilisation**

Most implants are supplied to the market in sterile packaging and protective packaging. The sterilization method is indicated on the package label. Implants are exposed to a minimum of 25 kGy gamma irradiation. Implants and hand tools are provided non-sterile if not labeled sterile on the label. It must be sterilized before use. Check the sterile package for holes, tears, etc. before the operation. UHMWPE Implants are not allowed to be re-sterilized. Only High Crosslink UHMWPE Materials Ethylene Oxide is sterilized.

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TIPMED Medical Device Manufacturing Ltd. Co.