



SHOULDER REPLACEMENT

*Regain
your
mobility*



SHOULDER

SHOULDER REPLACEMENT

TPM REVERSE SHOULDER COMPONENTS	PAGE NO
TPM REVERSE SHOULDER STANDARD HUMERAL STEM	C3
TPM REVERSE SHOULDER MINI HUMERAL STEM	C4
TPM REVERSE SHOULDER MICRO HUMERAL STEM	C5
TPM REVERSE SHOULDER FRACTURE HUMERAL STEM	C6
TPM REVERSE SHOULDER REVISION HUMERAL STEM	C7
TPM REVERSE SHOULDER HUMERAL TRAY	C8
TPM REVERSE SHOULDER HUMERAL CUP	C8
TPM REVERSE SHOULDER BASEPLATE	C9
TPM REVERSE SHOULDER BASEPLATE WITH EXTENSION (ECCENTRIC)	C9
TPM REVERSE SHOULDER STANDARD BASEPLATE WITH EXTENSION	C9
TPM REVERSE SHOULDER GLENOSPHERE (ECCENTRIC)	C10
TPM REVERSE SHOULDER STANDARD GLENOSPHERE	C10
TPM REVERSE SHOULDER SCREW OPTIONS	C11
TPM REVERSE SHOULDER HUMERAL HEAD (ECCENTRIC)	C12
TPM REVERSE SHOULDER STANDARD HUMERAL HEAD	C13

TPM REVERSE SHOULDER INSTRUMENTS	PAGE NO
TPM REVERSE SHOULDER INSTRUMENT SET	C14
TPM REVERSE SHOULDER REVISION HUMERAL STEM INSTRUMENT SET	C25

TPM REVERSE SHOULDER STANDARD HUMERAL STEM



ORDER INFORMATION 8698791 ...

Stem Diameter	Stem Length	Cemented Porous Coated Titanium	Uncemented Porous Coated Titanium
7 mm	122 mm	385608	385462
8 mm	122 mm	385615	385479
9 mm	122 mm	385622	385486
10 mm	122 mm	385639	385493
11 mm	122 mm	385646	385509
12 mm	122 mm	385653	385516
13 mm	122 mm	385660	385523
14 mm	122 mm	385677	385530
15 mm	122 mm	385684	385547
16 mm *	122 mm	385691	385554
17 mm *	122 mm	385707	385561
18 mm *	122 mm	385714	385578
19 mm *	122 mm	385721	385585
20 mm *	122 mm	385738	385592

* Optional size

TPM REVERSE SHOULDER MINI HUMERAL STEM



ORDER INFORMATION 8698791 ...

Stem Diameter	Stem Length	Cemented Porous Coated Titanium	Uncemented Porous Coated Titanium
7 mm	83 mm	385882	385745
8 mm	83 mm	385899	385752
9 mm	83 mm	385905	385769
10 mm	83 mm	385912	385776
11 mm	83 mm	385929	385783
12 mm	83 mm	385936	385790
13 mm	83 mm	385943	385806
14 mm	83 mm	385950	385813
15 mm	83 mm	385967	385820
16 mm *	83 mm	385974	385837
17 mm *	83 mm	385981	385844
18 mm *	83 mm	385998	385851
19 mm *	83 mm	386001	385868
20 mm *	83 mm	386018	385875

TPM REVERSE SHOULDER MICRO HUMERAL STEM



ORDER INFORMATION 8698791 ...

Stem Diameter	Stem Length	Cemented Porous Coated Titanium	Uncemented Porous Coated Titanium
7 mm *	55 mm	386162	386025
8 mm *	55 mm	386179	386032
9 mm *	55 mm	386186	386049
10 mm *	55 mm	386193	386056
11 mm *	55 mm	386209	386063
12 mm *	55 mm	386216	386070
13 mm *	55 mm	386223	386087
14 mm *	55 mm	386230	386094
15 mm *	55 mm	386247	386100
16 mm *	55 mm	386254	386117
17 mm *	55 mm	386261	386124
18 mm *	55 mm	386278	386131
19 mm *	55 mm	386285	386148
20 mm *	55 mm	386292	386155

* Optional for all size

TPM REVERSE SHOULDER FRACTURE HUMERAL STEM



ORDER INFORMATION 8698791 ...

Stem Diameter	Stem Length	Cemented Plasma Spray CoCrMo	Uncemented Plasma Spray CoCrMo
6 mm	122 mm	386308	386353
8 mm	122 mm	386315	386360
10 mm	122 mm	386322	386377
12 mm	122 mm	386339	386384
14 mm	122 mm	386346	386391

TPM REVERSE SHOULDER REVISION HUMERAL STEM



ORDER INFORMATION 8698791 ...

Stem Diameter	Stem Length	Cemented Porous Coated Titanium	Uncemented Porous Coated Titanium
6 mm	160 mm	455783	455745
6 mm	194 mm	106982	107033
8 mm	160 mm	455790	455752
8 mm	194 mm	106999	107040
10 mm	160 mm	455806	455769
10 mm	194 mm	106951	107002
12 mm	160 mm	455813	455776
12 mm	194 mm	106968	107019

TPM REVERSE SHOULDER HUMERAL TRAY



ORDER INFORMATION 8698791 ...

Size	Thickness	HUMERAL TRAY Titanium	Size	Thickness	MINI HUMERAL TRAY Titanium
44 mm	0 mm	389194	40 mm	0 mm	106944
44 mm	+5 mm	389187	40 mm	+5 mm	106937
44 mm	+10 mm	389170	40 mm	+10 mm	106920


TPM REVERSE SHOULDER HUMERAL CUP



ORDER INFORMATION 8698791 ...

Size	Thickness	HUMERAL CUPS UHMWPE	Size	Thickness	MINI HUMERAL CUPS UHMWPE
36 mm	0 mm	386438	36 mm	0 mm	106876
36 mm	+3 mm	386445	36 mm	+3 mm	106869
38 mm	0 mm	386452	38 mm	0 mm	106890
38 mm	+3 mm	386469	38 mm	+3 mm	106883
41 mm	0 mm	386476	41 mm	0 mm	106913
41 mm	+3 mm	386483	41 mm	+3 mm	106906


TPM REVERSE SHOULDER BASEPLATE



ORDER INFORMATION 8698791 ...

Size	Eccentric Uncemented Porous + HA Coated Titanium	Standard Uncemented Porous + HA Coated Titanium
Ø 25 mm	389200	453161


TPM REVERSE SHOULDER BASEPLATE WITH EXTENSION (ECCENTRIC)



ORDER INFORMATION 8698791 ...

Size	Short Uncemented Porous + HA Coated Titanium	Long Uncemented Porous + HA Coated Titanium
Ø 25 mm	473596	473602

TPM REVERSE SHOULDER STANDARD BASEPLATE WITH EXTENSION



ORDER INFORMATION 8698791 ...

Size	Short Uncemented Porous + HA Coated Titanium	Long Uncemented Porous + HA Coated Titanium
Ø 25 mm	473619	473626

TPM REVERSE SHOULDER GLENOSPHERE (ECCENTRIC)



ORDER INFORMATION 8698791 ...

Size	Offset	CoCrMo	Titanium
36 mm	0 mm	409137	409045
36 mm	+3 mm	409144	409052
36 mm	+6 mm	409151	409069
38 mm	0 mm	409168	409076
38 mm	+3 mm	409175	409083
38 mm	+6 mm	409182	409090
41 mm	0 mm	409199	409106
41 mm	+3 mm	409205	409113
41 mm	+6 mm	409212	409120


TPM REVERSE SHOULDER STANDARD GLENOSPHERE



ORDER INFORMATION 8698791 ...

Size	Offset	CoCrMo	Titanium
36 mm	0 mm	453079	452980
36 mm	+3 mm	453086	452997
36 mm	+6 mm	453093	453000
38 mm	0 mm	453109	453017
38 mm	+3 mm	453116	453024
38 mm	+6 mm	453123	453031
41 mm	0 mm	453130	453048
41 mm	+3 mm	453147	453055
41 mm	+6 mm	453154	453062

TPM REVERSE SHOULDER SCREW OPTIONS



ORDER INFORMATION 8698791 ...

Length	Ø 4.75 mm Titanium Screw Self-Tapping	Ø 4.75 mm Locking Titanium Screw Self-Tapping	Length	Ø 6.5mm Titanium Screw Self-Tapping
15 mm *	386681	386759	20 mm *	386827
20 mm	386698	386766	25 mm *	386834
25 mm	386704	386773	30 mm *	386841
30 mm	386711	386780	35 mm	386858
35 mm	386728	386797	40 mm	386865
40 mm	386735	386803	45 mm	386872
45 mm	386742	386810	50 mm	386889
50 mm	388104	388098	55 mm	388111
			60 mm	388128
			65 mm	388135

* Optional size

**TPM REVERSE SHOULDER
HUMERAL HEAD (ECCENTRIC)**

ORDER INFORMATION 8698791 ...

Size	Thickness	CoCrMo	Titanium
38 mm *	19 mm	388838	389002
38 mm *	21 mm	388845	389019
40 mm	15 mm	387350	387251
40 mm	20 mm	387367	387268
40 mm	22 mm	387374	387275
42 mm	18 mm	388852	389026
42 mm	21 mm	388869	389033
42 mm	24 mm	388876	389040
44 mm *	15 mm	387381	387282
44 mm *	17 mm	387398	387299
44 mm *	22 mm	387404	387305
44 mm *	27 mm	387411	387312
46 mm	18 mm	388883	389057
46 mm	21 mm	388890	389064
46 mm	24 mm	388906	389071
46 mm *	27 mm	388913	389088
48 mm *	19 mm	387428	387329
48 mm *	24 mm	387435	387336
48 mm *	27 mm	387442	387343
50 mm	21 mm	388920	389095
50 mm	24 mm	388937	389101
50 mm	27 mm	388944	389118
54 mm *	21 mm	388951	389125
54 mm *	24 mm	388968	389132
54 mm *	27 mm	388975	389149
58 mm *	24 mm	388982	389156
58 mm *	27 mm	388999	389163

* Optional size

TPM REVERSE SHOULDER STANDARD HUMERAL HEAD



ORDER INFORMATION 8698791 ...

Size	Thickness	CoCrMo	Titanium
38 mm *	19 mm	390022	390299
38 mm *	21 mm	390039	390305
40 mm	15 mm	390046	390312
40 mm	20 mm	390053	390329
40 mm	22 mm	390060	390336
42 mm	18 mm	390077	390343
42 mm	21 mm	390084	390350
42 mm	24 mm	390091	390367
44 mm *	15 mm	390107	390374
44 mm *	17 mm	390114	390381
44 mm *	22 mm	390121	390398
44 mm *	27 mm	390138	390404
46 mm	18 mm	390145	390411
46 mm	21 mm	390152	390428
46 mm	24 mm	390169	390435
46 mm *	27 mm	390176	390442
48 mm *	19 mm	390183	390459
48 mm *	24 mm	390190	390466
48 mm *	27 mm	390206	390473
50 mm	21 mm	390213	390480
50 mm	24 mm	390220	390497
50 mm	27 mm	390237	390503
54 mm *	21 mm	390244	390510
54 mm *	24 mm	390251	390527
54 mm *	27 mm	390268	390534
58 mm *	24 mm	390275	390541
58 mm *	27 mm	390282	390558

* Optional size

TPM REVERSE SHOULDER INSTRUMENT SET 90117021000


TRAY NO 1 - 90117021001



Item No.	Description	Unit
390848	TPM Reverse Shoulder System, Trial Fracture Humeral Stem, CoCrMo, 6*122	1
390855	TPM Reverse Shoulder System, Trial Fracture Humeral Stem, CoCrMo, 8*122	1
390862	TPM Reverse Shoulder System, Trial Fracture Humeral Stem, CoCrMo, 10*122	1
390879	TPM Reverse Shoulder System, Trial Fracture Humeral Stem, CoCrMo, 12*122	1
390886	TPM Reverse Shoulder System, Trial Fracture Humeral Stem, CoCrMo, 14*122	1
387459	TPM Reverse Shoulder System, Standard-Mini Humeral Stem Rasp, S Steel, 7 mm	1
387466	TPM Reverse Shoulder System, Standard-Mini Humeral Stem Rasp, S Steel, 8 mm	1
387473	TPM Reverse Shoulder System, Standard-Mini Humeral Stem Rasp, S Steel, 9 mm	1
387480	TPM Reverse Shoulder System, Standard-Mini Humeral Stem Rasp, S Steel, 10 mm	1
387497	TPM Reverse Shoulder System, Standard-Mini Humeral Stem Rasp, S Steel, 11 mm	1
387503	TPM Reverse Shoulder System, Standard-Mini Humeral Stem Rasp, S Steel, 12 mm	1
387510	TPM Reverse Shoulder System, Standard-Mini Humeral Stem Rasp, S Steel, 13 mm	1
387527	TPM Reverse Shoulder System, Standard-Mini Humeral Stem Rasp, S Steel, 14 mm	1
387534	TPM Reverse Shoulder System, Standard-Mini Humeral Stem Rasp, S Steel, 15 mm	1



290430055530240
St. Steel Container
530 x 240 x 55 mm



Item No.	Description	Unit
308898	TPM Reverse Shoulder System, Humeral Stem Rasp & Humeral Stem Extractor, S Steel	1
262572	Hammer (Silicon Handle), S Steel	1
308867	TPM Reverse Shoulder System, 2-Prong Glenosphere Inserter / Impactor	1
308911	TPM Reverse Shoulder System, Humeral Stem Rasp & Humeral Stem Impactor, S Steel	2



TPM REVERSE SHOULDER INSTRUMENT SET 90117021000

TRAY NO 2 - 90117021002



Item No.	Description	Unit
387732	TPM Reverse Shoulder System, Trial Glenosphere, Plastic, 36, 0 (Eccentric)	1
387749	TPM Reverse Shoulder System, Trial Glenosphere, Plastic, 36, +3 (Eccentric)	1
387756	TPM Reverse Shoulder System, Trial Glenosphere, Plastic, 36, +6 (Eccentric)	1
387763	TPM Reverse Shoulder System, Trial Glenosphere, Plastic, 38, 0 (Eccentric)	1
387770	TPM Reverse Shoulder System, Trial Glenosphere, Plastic, 38, +3 (Eccentric)	1
387787	TPM Reverse Shoulder System, Trial Glenosphere, Plastic, 38, +6 (Eccentric)	1
387794	TPM Reverse Shoulder System, Trial Glenosphere, Plastic, 41, 0 (Eccentric)	1
387800	TPM Reverse Shoulder System, Trial Glenosphere, Plastic, 41, +3 (Eccentric)	1
387817	TPM Reverse Shoulder System, Trial Glenosphere, Plastic, 41, +6 (Eccentric)	1
453178	TPM Reverse Shoulder System, Trial Standard Glenosphere, Plastic, 36, 0	1
453185	TPM Reverse Shoulder System, Trial Standard Glenosphere, Plastic, 36, +3	1
453192	TPM Reverse Shoulder System, Trial Standard Glenosphere, Plastic, 36, +6	1
453208	TPM Reverse Shoulder System, Trial Standard Glenosphere, Plastic, 38, 0	1
453215	TPM Reverse Shoulder System, Trial Standard Glenosphere, Plastic, 38, +3	1
453222	TPM Reverse Shoulder System, Trial Standard Glenosphere, Plastic, 38, +6	1
453239	TPM Reverse Shoulder System, Trial Standard Glenosphere, Plastic, 41, 0	1
453246	TPM Reverse Shoulder System, Trial Standard Glenosphere, Plastic, 41, +3	1
453253	TPM Reverse Shoulder System, Trial Standard Glenosphere, Plastic, 41, +6	1



Item No.	Description	Unit
308430	TPM Reverse Shoulder System, Trial Humeral Tray & Cup, Plastic, 36, 0, 0	1
308447	TPM Reverse Shoulder System, Trial Humeral Tray & Cup, Plastic, 36, 5, 0	1
308454	TPM Reverse Shoulder System, Trial Humeral Tray & Cup, Plastic, 36, 10, 0	1
308461	TPM Reverse Shoulder System, Trial Humeral Tray & Cup, Plastic, 36, 0, +3	1
308478	TPM Reverse Shoulder System, Trial Humeral Tray & Cup, Plastic, 36, 5, +3	1
308485	TPM Reverse Shoulder System, Trial Humeral Tray & Cup, Plastic, 36, 10, +3	1
308492	TPM Reverse Shoulder System, Trial Humeral Tray & Cup, Plastic, 38, 0, 0	1
308508	TPM Reverse Shoulder System, Trial Humeral Tray & Cup, Plastic, 38, 5, 0	1
308515	TPM Reverse Shoulder System, Trial Humeral Tray & Cup, Plastic, 38, 10, 0	1
308522	TPM Reverse Shoulder System, Trial Humeral Tray & Cup, Plastic, 38, 0, +3	1
308539	TPM Reverse Shoulder System, Trial Humeral Tray & Cup, Plastic, 38, 5, +3	1
308546	TPM Reverse Shoulder System, Trial Humeral Tray & Cup, Plastic, 38, 10, +3	1
308553	TPM Reverse Shoulder System, Trial Humeral Tray & Cup, Plastic, 41, 0, 0	1
308560	TPM Reverse Shoulder System, Trial Humeral Tray & Cup, Plastic, 41, 5, 0	1
308577	TPM Reverse Shoulder System, Trial Humeral Tray & Cup, Plastic, 41, 10, 0	1
308584	TPM Reverse Shoulder System, Trial Humeral Tray & Cup, Plastic, 41, 0, +3	1
308591	TPM Reverse Shoulder System, Trial Humeral Tray & Cup, Plastic, 41, 5, +3	1
308607	TPM Reverse Shoulder System, Trial Humeral Tray & Cup, Plastic, 41, 10, +3	1
453260	TPM Reverse Shoulder System, Trial Mini Humeral Tray & Mini Cup, Plastic, 36, 0, 0	1
453277	TPM Reverse Shoulder System, Trial Mini Humeral Tray & Mini Cup, Plastic, 36, 5, 0	1
453284	TPM Reverse Shoulder System, Trial Mini Humeral Tray & Mini Cup, Plastic, 36, 10, 0	1
453291	TPM Reverse Shoulder System, Trial Mini Humeral Tray & Mini Cup, Plastic, 36, 0, +3	1
453307	TPM Reverse Shoulder System, Trial Mini Humeral Tray & Mini Cup, Plastic, 36, 5, +3	1
453314	TPM Reverse Shoulder System, Trial Mini Humeral Tray & Mini Cup, Plastic, 36, 10, +3	1
453321	TPM Reverse Shoulder System, Trial Mini Humeral Tray & Mini Cup, Plastic, 38, 0, 0	1
453338	TPM Reverse Shoulder System, Trial Mini Humeral Tray & Mini Cup, Plastic, 38, 5, 0	1
453345	TPM Reverse Shoulder System, Trial Mini Humeral Tray & Mini Cup, Plastic, 38, 10, 0	1
453352	TPM Reverse Shoulder System, Trial Mini Humeral Tray & Mini Cup, Plastic, 38, 0, +3	1
453369	TPM Reverse Shoulder System, Trial Mini Humeral Tray & Mini Cup, Plastic, 38, 5, +3	1
453376	TPM Reverse Shoulder System, Trial Mini Humeral Tray & Mini Cup, Plastic, 38, 10, +3	1
453383	TPM Reverse Shoulder System, Trial Mini Humeral Tray & Mini Cup, Plastic, 41, 0, 0	1
453390	TPM Reverse Shoulder System, Trial Mini Humeral Tray & Mini Cup, Plastic, 41, 5, 0	1
453406	TPM Reverse Shoulder System, Trial Mini Humeral Tray & Mini Cup, Plastic, 41, 10, 0	1
453413	TPM Reverse Shoulder System, Trial Mini Humeral Tray & Mini Cup, Plastic, 41, 0, +3	1
453420	TPM Reverse Shoulder System, Trial Mini Humeral Tray & Mini Cup, Plastic, 41, 5, +3	1
453437	TPM Reverse Shoulder System, Trial Mini Humeral Tray & Mini Cup, Plastic, 41, 10, +3	1



TPM REVERSE SHOULDER INSTRUMENT SET 90117021000

TRAY NO 2 - 90117021002



Item No.	Description	Unit
387848	TPM Shoulder System, Trial Humeral Head, Plastic, 40*15 (Eccentric)	1
387855	TPM Shoulder System, Trial Humeral Head, Plastic, 40*20 (Eccentric)	1
387862	TPM Shoulder System, Trial Humeral Head, Plastic, 40*22 (Eccentric)	1
387879	TPM Shoulder System, Trial Humeral Head, Plastic, 42*18 (Eccentric)	1
387886	TPM Shoulder System, Trial Humeral Head, Plastic, 42*21 (Eccentric)	1
387893	TPM Shoulder System, Trial Humeral Head, Plastic, 42*24 (Eccentric)	1
387947	TPM Shoulder System, Trial Humeral Head, Plastic, 46*18 (Eccentric)	1
387954	TPM Shoulder System, Trial Humeral Head, Plastic, 46*21 (Eccentric)	1
387961	TPM Shoulder System, Trial Humeral Head, Plastic, 46*24 (Eccentric)	1
388012	TPM Shoulder System, Trial Humeral Head, Plastic, 50*21 (Eccentric)	1
388029	TPM Shoulder System, Trial Humeral Head, Plastic, 50*24 (Eccentric)	1
388036	TPM Shoulder System, Trial Humeral Head, Plastic, 50*27 (Eccentric)	1



290430100530240
St. Steel Container
530 x 240 x 100 mm



Item No.	Description	Unit
390589	TPM Shoulder System, Trial Standard Humeral Head, Plastic, 40*15	1
390596	TPM Shoulder System, Trial Standard Humeral Head, Plastic, 40*20	1
390602	TPM Shoulder System, Trial Standard Humeral Head, Plastic, 40*22	1
390619	TPM Shoulder System, Trial Standard Humeral Head, Plastic, 42*18	1
390626	TPM Shoulder System, Trial Standard Humeral Head, Plastic, 42*21	1
390633	TPM Shoulder System, Trial Standard Humeral Head, Plastic, 42*24	1
390688	TPM Shoulder System, Trial Standard Humeral Head, Plastic, 46*18	1
390695	TPM Shoulder System, Trial Standard Humeral Head, Plastic, 46*21	1
390701	TPM Shoulder System, Trial Standard Humeral Head, Plastic, 46*24	1
390756	TPM Shoulder System, Trial Standard Humeral Head, Plastic, 50*21	1
390763	TPM Shoulder System, Trial Standard Humeral Head, Plastic, 50*24	1
390770	TPM Shoulder System, Trial Standard Humeral Head, Plastic, 50*27	1



TPM REVERSE SHOULDER INSTRUMENT SET 90117021000

TRAY NO 3 - 90117021003




Item No.	Description	Unit
448099	TPM Reverse Shoulder System, Cutting Guide Handle, S Steel	1
448297	TPM Reverse Shoulder System, Cutting Guide , S Steel	1
448259	TPM Reverse Shoulder System, Cutting Guide Free Style, S Steel	1
214991	Feeler Blade, Long, S Steel	1
214984	Feeler Blade, Small, S Steel	1
448105	TPM Reverse Shoulder System, T Handle	1
448112	TPM Reverse Shoulder System, I Handle	1
308928	TPM Reverse Shoulder System, Glenoid Planer with Blades, S Steel	1
473640	TPM Reverse Shoulder System, Reamer For Baseplate With Extension, S Steel	1
473657	TPM Reverse Shoulder System, Stopper For Baseplate With Extension Reamer, S Steel	1



290430055530240

St. Steel Container
530 x 240 x 55 mm



Item No.	Description	Unit
309024	TPM Reverse Shoulder System, Humeral Tray Adaptor, 0, S Steel	1
309031	TPM Reverse Shoulder System, Humeral Tray Adaptor, 5, S Steel	1
309017	TPM Reverse Shoulder System, Humeral Cup & Tray Compression, S Steel	1
390909	TPM Reverse Shoulder System, Glenoid Reamer, S Steel	1
309055	TPM Reverse Shoulder System, Cannulated Starting Reamer For Glenoid, S Steel	1
387596	TPM Reverse Shoulder System, Humeral Reamer, S Steel, 7	1
387602	TPM Reverse Shoulder System, Humeral Reamer, S Steel, 8	1
387619	TPM Reverse Shoulder System, Humeral Reamer, S Steel, 9	1
387626	TPM Reverse Shoulder System, Humeral Reamer, S Steel, 10	1
387633	TPM Reverse Shoulder System, Humeral Reamer, S Steel, 11	1
387640	TPM Reverse Shoulder System, Humeral Reamer, S Steel, 12	1
387657	TPM Reverse Shoulder System, Humeral Reamer, S Steel, 13	1
387664	TPM Reverse Shoulder System, Humeral Reamer, S Steel, 14	1
387671	TPM Reverse Shoulder System, Humeral Reamer, S Steel, 15	1



TPM REVERSE SHOULDER INSTRUMENT SET 90117021000

TRAY NO 5 - 90117021005



Item No.	Description	Unit
448143	TPM Reverse Shoulder System, Head Impactor	1
448129	TPM Reverse Shoulder System, Humeral Cup Impactor	1
448136	TPM Reverse Shoulder System, Glenoid Baseplate Impactor	1
308935	TPM Reverse Shoulder System, Baseplate Reamer, S Steel	1
308997	TPM Reverse Shoulder System, AO Shaft Screw Driver, S Steel, 2.5	1
448303	TPM Reverse Shoulder System, AO Shaft Screw Driver, S Steel, 3.5	1
390831	TPM Reverse Shoulder System, Trial Base Plate, Plastic	1
457046	Screw Driver for Screw Holder AA2.5	1
219491	Screw Driver for Screw Holder AA3.5	1
215127	Pin Extractor, S Steel	1
215219	Pin With Head, 55, S Steel	2
215226	Pin With Head, 65, S Steel	2
359876	TPM Rev. Knee Sys., Pin with Drill, 80 mm, S Steel	2
308942	TPM Reverse Shoulder System, Pin with Drill, 100 mm, S Steel	2



290430055530240

St. Steel Container
530 x 240 x 55 mm



Item No.	Description	Unit
244103	Pin Adapter, S Steel	1
448082	TPM Reverse Shoulder System, Sizer with 10 Degree Inferior Tilt Handle, S Steel	2
448167	TPM Reverse Shoulder System, Guide Handle, S Steel, Long	2
448174	TPM Reverse Shoulder System, Guide Handle, S Steel, Small	2
448266	TPM Reverse Shoulder System, Baseplate Central Screw Drill Guide	1
448280	TPM Reverse Shoulder System, Peripheral Drill and Screw Guide, S Steel	1
309048	TPM Reverse Shoulder System, Baseplate Drill Guide, S Steel	2
308980	TPM Reverse Shoulder System, Depth Gauge for 6.5 Shoulder Screw, S Steel	1
308973	TPM Reverse Shoulder System, Depth Gauge for 4.75 Shoulder Screw, S Steel	1
448273	TPM Reverse Shoulder System, Sizer with 10 Degree Inferior Tilt, S Steel	1
473664	TPM Reverse Shoulder System, Drill Guide For Baseplate With Extension, S Steel	1
308959	TPM Reverse Shoulder System, Drill with Scale for 6.5 Shoulder Screw, S Steel	1
308966	TPM Reverse Shoulder System, Drill with Scale for 4.75 Shoulder Screw, S Steel	1
085072	Drill with Round Shaft, S Steel, 3.2*200	2
448150	TPM Reverse Shoulder System, Guide Pin, S Steel, 3,2*230	2
086505	3.2 Trocar Point Kirschner Wire, S Steel, 230	2
309000	TPM Reverse Shoulder System, Humeral Cup Extractor, S Steel	1



TPM REVERSE SHOULDER INSTRUMENT SET 90117021000

TRAY NO 4 - 90117021004

29043005530240

St. Steel Container
530 x 240 x 55 mm



Item No.	Description	Unit
390893	TPM Reverse Shoulder System, Humeral Resection Plate, S Steel	1
448327	TPM Reverse Shoulder System, Humeral Head Removal Fork, S Steel	1
448310	TPM Reverse Shoulder System, Humeral Tray Removal Fork, S Steel	1
308850	TPM Reverse Shoulder System, Shoehorn	1
448235	TPM Reverse Shoulder System, Wide Glenoid Retractor, S Steel	1
448242	TPM Reverse Shoulder System, Thin Glenoid Retractor, S Steel	1
448228	TPM Reverse Shoulder System, 2-Prong Capsular Retractor, S Steel	1
448211	TPM Reverse Shoulder System, Modified Darrach Retractor, S Steel	1
308904	TPM Reverse Shoulder System, Slide Hammer, S Steel	1
448204	TPM Reverse Shoulder System, Golf Club Retractor, S Steel	1
448198	TPM Reverse Shoulder System, Bent Ring Fukuda	1
448181	TPM Reverse Shoulder System, Large Ring Fukuda	1



TPM REVERSE SHOULDER REVISION HUMERAL STEM INSTRUMENT SET 90125001101

TRAY NO 1

290430055530240

St. Steel Container
530 x 240 x 55 mm



Item No.	Description	Unit
478768	TPM Reverse Shoulder System, Revision Humeral Stem Rasp, S Steel, 6 mm	1
478775	TPM Reverse Shoulder System, Revision Humeral Stem Rasp, S Steel, 8 mm	1
478782	TPM Reverse Shoulder System, Revision Humeral Stem Rasp, S Steel, 10 mm	1
478799	TPM Reverse Shoulder System, Revision Humeral Stem Rasp, S Steel, 12 mm	1
478683	TPM Reverse Shoulder System, Revision Humeral Reamer, S Steel, 6*160	1
478706	TPM Reverse Shoulder System, Revision Humeral Reamer, S Steel, 8*160	1
478720	TPM Reverse Shoulder System, Revision Humeral Reamer, S Steel, 10*160	1
478744	TPM Reverse Shoulder System, Revision Humeral Reamer, S Steel, 12*160	1
478690	TPM Reverse Shoulder System, Revision Humeral Reamer, S Steel, 6*194	1
478713	TPM Reverse Shoulder System, Revision Humeral Reamer, S Steel, 8*194	1
478737	TPM Reverse Shoulder System, Revision Humeral Reamer, S Steel, 10*194	1
478751	TPM Reverse Shoulder System, Revision Humeral Reamer, S Steel, 12*194	1
478805	TPM Reverse Shoulder System, Revision Humeral Stem Rasp & Revision Humeral Stem Impactor, S Steel	2
478812	TPM Reverse Shoulder System, Revision Humeral Stem Glenoid Planer with Blades, S Steel	1



Regain
your
mobility



RIISING SUN OF YOUR LIFE

You can feel the
experience, confidence
and quality in field
of arthroplasty,
traumatology and
spine.



Company Profile:

Tipmed is a prominent medical products manufacturer, specializing in implants and instruments in the field of orthopedics, traumatology and neurochirurgie. It was established in 2001 in Izmir by the industry veteran Sedat Onel to provide a path-breaking solution to the most demanding medical device manufacturing challenges.

Products:

Orthopedic instruments and implant systems such as Total Hip Replacement, Knee and Shoulder Replacement, Spine Systems, External Fixators and Trauma Systems like Locking Plating Systems, Intramedullary Nailing Systems, Screw Systems.

www.tipmed.com

*Regain
your
mobility*



TIPMED Medical Device Manufacturing Ltd. Co.

Kemalpaşa Mahallesi 7401/1 Sokak No:16
Bornova / İzmir - TÜRKİYE

T : +90 232 479 33 22

M : +90 543 479 33 22

F : +90 232 479 38 25

info@tipmed.com

marketing@tipmed.com

www.tipmed.com

