

Regain
your
mobility

S
S
R
M
T
P



TPM
Reverse
SHOULDER
SYSTEM

STANDARD HUMERAL STEM



Titanium Alloyed,
Cemented &
Uncemented,
Porous Coated

Length : 122 mm
Diameter: 7 mm
up to 15 mm with
1 mm increment

Optional Size:
16-17-18-19-20 mm

MINI HUMERAL STEM



Titanium Alloyed,
Cemented &
Uncemented,
Porous Coated

Length : 83 mm
Diameter: 7 mm
up to 15 mm with
1 mm increment

Optional Size:
16-17-18-19-20 mm

MICRO HUMERAL STEM



Titanium Alloyed,
Cemented &
Uncemented,
Porous Coated

Length : 55 mm
Diameter: 7 mm
up to 15 mm with
1 mm increment

Optional Size:
16-17-18-19-20 mm


FRACTURE HUMERAL STEM



CoCrMo Alloyed,
Cemented &
Uncemented,
Plasma Spray
Coated

Length : 122 mm
Diameter: 6 mm
up to 14 mm with
2 mm increment

REVISION HUMERAL STEM



Titanium Alloyed,
Cemented &
Uncemented,
Porous Coated

Length :
160 mm, 194 mm
Diameter: 6 mm
up to 12 mm with
2 mm increment

HUMERAL TRAY

Titanium Alloyed
Thickness : 0 - 5 - 10 mm, Size : 44 mm

MINI HUMERAL TRAY

Titanium Alloyed
Thickness : 0 - 5 - 10 mm, Size : 40 mm



TPM Reverse SHOULDER SYSTEM

BASEPLATE (ECCENTRIC) & STANDARD BASEPLATE

Titanium Alloyed
Uncemented,
Porous + H.A. Coated
Size : Ø25 mm



HUMERAL CUP

UHMWPE Alloyed
Thickness : 0 - 3 mm, Size : 36 - 38 - 41 mm

MINI HUMERAL CUP

UHMWPE Alloyed
Thickness : 0 - 3 mm, Size : 36 - 38 - 41 mm



HUMERAL HEAD (ECCENTRIC)

Titanium & CoCrMo Alloyed
Size : 40 mm / 15 - 20 - 22 mm
42 - 46 mm / 18 - 21 - 24 mm
50 mm / 21 - 24 - 27 mm
Optional Size : 38 - 44 - 48 - 54 - 58



BASEPLATE (ECCENTRIC) WITH EXT. Short/Long & STANDARD BASEPLATE WITH EXT. Short/Long

Titanium Alloyed
Uncemented,
Porous + H.A. Coated
Size : Ø25 mm



GLENOSPHERE (ECCENTRIC) & STANDARD GLENOSPHERE

Titanium & CoCrMo Alloyed
Offset : 0 - 3 - 6 mm
Size : 36 - 38 - 41 mm



STANDARD HUMERAL HEAD

Titanium & CoCrMo Alloyed
Size : 40 mm / 15 - 20 - 22 mm
42 - 46 mm / 18 - 21 - 24 mm
50 mm / 21 - 24 - 27 mm
Optional Size : 38 - 44 - 48 - 54 - 58



4,75 MM SHOULDER SCREW

Titanium Alloyed
Locking & Non-Locking
Length : 20 - 25 - 30 -
35 - 40 - 45 - 50 mm



6,5 MM SHOULDER SCREW

Titanium Alloyed
Length : 20 - 25 - 30 - 35 -
40 - 45 - 50 55 - 60 - 65 mm

