



KNEE REPLACEMENT

*Regain
your
mobility*

B

KNEE REPLACEMENT

TPM8 PRIMARY KNEE COMPONENTS **PAGE NO**

FEMORAL COMPONENTS CRUCIATE RETAINING (CR) B6

FEMORAL COMPONENTS POSTERIOR STABILIZED (PS) B7

TPM8 TIBIAL METAL BACK COMPONENTS (TOP FITTINGS) B8

TIBIAL METAL BACK COMPONENTS (SIDE SHIFTING) B9

INSERT CRUCIATE RETAINING (CR) SIDE SHIFTING, UHMWPE B10

INSERT POSTERIOR STABILIZED (PS) SIDE SHIFTING, UHMWPE B12

TPM8 INSERT CRUCIATE RETAINING (CR) TOP FITTINGS, UHMWPE B14

TPM8 INSERT POSTERIOR STABILIZED (PS) TOP FITTINGS, UHMWPE B16

PATELLA COMPONENT B18

TIBIAL METAL BACK STEM B18

TPM8 PRIMARY KNEE INSTRUMENTS B19**TPM REVISION KNEE COMPONENTS** **PAGE NO**

FEMORAL PS COMPONENTS B30

POSTERIOR FEMORAL AUGMENT B31

DISTAL FEMORAL AUGMENT B32

EXTENSION STEM B33

TIBIAL AUGMENT B34

TIBIAL METAL BACK B35

UHMWPE POSTERIOR STABILIZED (PS) TIBIAL INSERT B36

UHMWPE PLUS TIBIAL INSERT B37

OFFSET ADAPTER B38

PATELLA COMPONENT B38

TPM8 REVISION KNEE INSTRUMENTS B39



TPM8
Primary
KNEE

FEMORAL COMPONENTS

CoCrMo Alloyed

Cruciate Retaining (CR)



Posterior Stabilized (PS)



Cemented

TIBIAL METAL BACK

CoCrMo Alloyed

TPM8 Tibial Metal Back
Top Fittings



Tibial Metal Back
Side Shifting



Fixed Bearing
Mobile Bearing

Cemented

Fixed Bearing
Mobile Bearing

TIBIAL PE INSERT

UHMWPE - High Crosslink

Cruciate Retaining (CR)
Top Fittings



Posterior Stabilized (PS)
Top Fittings



Fixed Bearing
Mobile Bearing

Cruciate Retaining (CR)
Side Shifting



Posterior Stabilized (PS)
Side Shifting



Fixed Bearing
Mobile Bearing

Deep Dish Insert
Cruciate Retaining (CR)
Side Shifting



Mobile Bearing

PATELLA COMPONENT

UHMWPE Alloyed



TIBIAL METAL BACK STEM

Titanium Alloyed




**FEMORAL COMPONENTS
CRUCIATE RETAINING (CR)**

ORDER INFORMATION 8698791 ...

Optional *

Size	Femoral Component (CR) Right Side	Femoral Component (CR) Left Side
*54 mm	179344	179429
*56 mm	179351	179436
60 mm	179368	179443
64 mm	179375	179450
*66 mm	398981	399001
68 mm	179382	179467
*70 mm	398998	399018
72 mm	179399	179474
76 mm	179405	179481
80 mm	179412	179498

FEMORAL COMPONENTS POSTERIOR STABILIZED (PS)



ORDER INFORMATION 8698791 ...

Optional *

Size	Femoral Component (PS) Right Side	Femoral Component (PS) Left Side
*54 mm	179559	179634
*56 mm	179566	179641
60 mm	179573	179658
64 mm	179580	179665
*66 mm	399148	399162
68 mm	179597	179672
*70 mm	399155	399179
72 mm	179603	179689
76 mm	179610	179696
80 mm	179627	179702

Ancorge Peg for
Femoral Component
Titanium

194033

TPM8 TIBIAL METAL BACK COMPONENTS (TOP FITTINGS)



ORDER INFORMATION 8698791 ...

Optional *

Size	Fixed Tibial Metal Back Cemented	Mobile Tibial Metal Back Cemented
*54 mm	190332	245391
*56 mm	190349	245407
60 mm	190356	245414
64 mm	190363	245421
*66 mm	399568	399469
68 mm	190370	245438
*70 mm	399575	399476
72 mm	190387	245445
76 mm	190394	245452
80 mm	190400	245469

TPM8 Insert Holder Titanium

Mobile	Fixed
408925	408932



Ancorge Peg for Tibia Metal Back, Titanium

194040

TIBIAL METAL BACK COMPONENTS (SIDE SHIFTING)



Ancorge Peg for
Tibia Metal Back



ORDER INFORMATION 8698791 ...

Optional *

Size	Fixed Tibial Metal Back Cemented	Mobile Tibial Metal Back Cemented
*54 mm	180142	187592
*56 mm	180159	187608
60 mm	180166	187615
64 mm	180173	187622
*66 mm	399384	399308
68 mm	180180	187639
*70 mm	399391	399315
72 mm	213406	187646
76 mm	180197	187653
80 mm	180203	187660

Insert Holder Titanium


Mobile	Fixed
179719	179726



Ancorage Peg for Tibia Metal Back, Titanium

194040

INSERT CRUCIATE RETAINING (CR) SIDE SHIFTING, UHMWPE



ORDER INFORMATION 8698791 ...

Optional *

Size	CR Mobile Bearing UHMWPE	CR Mobile Bearing High Crosslink	CR Fixed Bearing UHMWPE	CR Fixed Bearing High Crosslink	CR Mobile UHMWPE Deep Dish	CR Mobile High Crosslink UHMWPE Deep Dish
*8 x 54 mm	187837	359418	187677	473992	180319	185291
*8 x 56 mm	187844	190912	187684	474005	180326	185307
8 x 60 mm	187851	356066	187691	355700	180333	185314
8 x 64 mm	187868	356073	187707	355717	180340	185321
*8 x 66 mm	399667	399902	404163	474012	404682	404842
8 x 68 mm	187875	356080	187714	355724	180357	185338
*8 x 70 mm	399674	399919	404170	474029	404699	404859
8 x 72 mm	187882	356097	187721	355731	180364	185345
8 x 76 mm	187899	356103	187738	355748	180371	185352
8 x 80 mm	187905	356110	187745	355755	180388	185369
*10 x 54 mm	191117	359364	179733	474036	180395	185376
*10 x 56 mm	191124	190929	179740	474043	180401	185383
10 x 60 mm	191131	355762	179757	355403	180418	185390
10 x 64 mm	191148	355779	179764	355410	180425	185406
*10 x 66 mm	399681	399926	404187	474050	404705	404866
10 x 68 mm	191155	355786	179771	355427	180432	185413
*10 x 70 mm	399698	399933	404194	474067	404712	404873
10 x 72 mm	191162	355793	179788	355434	180449	185420
10 x 76 mm	191179	355809	179795	355441	180456	185437
10 x 80 mm	191186	355816	179801	355458	180463	185444

INSERT CRUCIATE RETAINING (CR) SIDE SHIFTING, UHMWPE

ORDER INFORMATION 8698791 ...

Optional *

Size	CR Mobile Bearing UHMWPE	CR Mobile Bearing High Crosslink	CR Fixed Bearing UHMWPE	CR Fixed Bearing High Crosslink	CR Mobile UHMWPE Deep Dish	CR Mobile High Crosslink UHMWPE Deep Dish
*12 x 54 mm	191193	359371	179818	474074	180470	185451
*12 x 56 mm	191209	190936	179825	474081	180487	185468
12 x 60 mm	191216	355823	179832	355465	180494	185475
12 x 64 mm	191223	355830	179849	355472	180500	185482
*12 x 66 mm	399704	399940	404200	474098	404729	404880
12 x 68 mm	191230	355847	179856	355489	180517	185499
*12 x 70 mm	399711	399957	404217	474104	404736	404897
12 x 72 mm	191247	355854	179863	355496	180524	185505
12 x 76 mm	191254	355861	179870	355502	180531	185512
12 x 80 mm	191261	355878	179887	355519	180548	185529
*14 x 54 mm	191278	359388	179894	474111	180555	185536
*14 x 56 mm	191285	190943	179900	474128	180562	185543
14 x 60 mm	191292	355885	179917	355526	180579	185550
14 x 64 mm	191308	355892	179924	355533	180586	185567
*14 x 66 mm	399728	399964	404224	474135	404743	404903
14 x 68 mm	191315	355908	179931	355540	180593	185574
*14 x 70 mm	399735	399971	404231	474142	404750	404910
14 x 72 mm	191322	355915	179948	355557	180609	185581
14 x 76 mm	191339	355922	179955	355564	180616	185598
14 x 80 mm	191346	355939	179962	355571	180623	185604
*16 x 54 mm	191353	359395	179979	474159	180630	185611
*16 x 56 mm	191360	190950	179986	474166	180647	185628
16 x 60 mm	191377	355946	179993	355588	180654	185635
16 x 64 mm	191384	355953	180005	355595	180661	185642
*16 x 66 mm	399742	399988	404248	474173	404767	404927
16 x 68 mm	191391	355960	180012	355601	180678	185659
*16 x 70 mm	399759	399995	404255	474180	404774	404934
16 x 72 mm	191407	355977	180029	355618	180685	185666
16 x 76 mm	191414	355984	180036	355625	180692	185673
16 x 80 mm	191421	355991	180043	355632	180708	185680
*18 x 54 mm	191438	359401	180050	474197	180715	185697
*18 x 56 mm	191445	190967	180067	474203	180722	185703
18 x 60 mm	191452	356004	180074	355649	180739	185710
18 x 64 mm	191469	356011	180081	355656	180746	185727
*18 x 66 mm	399766	403906	404262	474210	404781	404941
18 x 68 mm	191476	356028	180098	355663	180753	185734
*18 x 70 mm	399773	403913	404279	474227	404798	404958
18 x 72 mm	191483	356035	180104	355670	180760	185741
18 x 76 mm	191490	356042	180111	355687	180777	185758
18 x 80 mm	191506	356059	180128	355694	180784	185765

INSERT POSTERIOR STABILIZED (PS) SIDE SHIFTING, UHMWPE



ORDER INFORMATION 8698791 ...

Optional *

Size	PS Mobile Bearing UHMWPE	PS Mobile Bearing High Crosslink	PS Fixed Bearing UHMWPE	PS Fixed Bearing High Crosslink	PS Mobile UHMWPE Deep Dish	PS Mobile High Crosslink UHMWPE Deep Dish
*8 x 54 mm	187912	407447	187752	474234	184096	186496
*8 x 56 mm	187929	407454	187769	474241	184102	186502
8 x 60 mm	187936	407461	187776	356783	184119	186519
8 x 64 mm	187943	407478	187783	356790	184126	186526
*8 x 66 mm	403920	407485	404446	474258	405009	405160
8 x 68 mm	187950	407492	187790	356806	184133	186533
*8 x 70 mm	403937	407508	404453	474265	405016	405177
8 x 72 mm	187967	407515	187806	356813	184140	186540
8 x 76 mm	187974	407522	187813	356820	184157	186557
8 x 80 mm	187981	407539	187820	356837	184164	186564
*10 x 54 mm	191513	407546	187196	474272	184171	186571
*10 x 56 mm	191520	407553	187202	474289	184188	186588
10 x 60 mm	191537	407560	187219	356486	184195	186595
10 x 64 mm	191544	407577	187226	356493	184201	186601
*10 x 66 mm	403944	407584	404460	474296	405023	405184
10 x 68 mm	191551	407591	187233	356509	184218	186618
*10 x 70 mm	403951	407607	404477	474302	405030	405191
10 x 72 mm	191568	407614	187240	356516	184225	186625
10 x 76 mm	191575	407621	187257	356523	184232	186632
10 x 80 mm	191582	407638	187264	356530	184249	186649

INSERT POSTERIOR STABILIZED (PS) SIDE SHIFTING, UHMWPE

ORDER INFORMATION 8698791 ...
Optional *

Size	PS Mobile Bearing UHMWPE	PS Mobile Bearing High Crosslink	PS Fixed Bearing UHMWPE	PS Fixed Bearing High Crosslink	PS Mobile UHMWPE Deep Dish	PS Mobile High Crosslink UHMWPE Deep Dish
*12 x 54 mm	191599	407645	187271	474319	184256	186656
*12 x 56 mm	191605	407652	187288	474326	184263	186663
12 x 60 mm	191612	407669	187295	356547	184270	186670
12 x 64 mm	191629	407676	187301	356554	184287	186687
*12 x 66 mm	403968	407683	404484	474333	405047	405207
12 x 68 mm	191636	407690	187318	356561	184294	186694
*12 x 70 mm	403975	407706	404491	474340	405054	405214
12 x 72 mm	191643	407713	187325	356578	184300	186700
12 x 76 mm	191650	407720	187332	356585	184317	186717
12 x 80 mm	191667	407737	187349	356592	184324	186724
*14 x 54 mm	191674	407744	187356	474357	184331	186731
*14 x 56 mm	191681	407751	187363	474364	184348	186748
14 x 60 mm	191698	407768	187370	356608	184355	186755
14 x 64 mm	191704	407775	187387	356615	184362	186762
*14 x 66 mm	403982	407782	404507	474371	405061	405221
14 x 68 mm	191711	407799	187394	356622	184379	186779
*14 x 70 mm	403999	407805	404514	474388	405078	405238
14 x 72 mm	191728	407812	187400	356639	184386	186786
14 x 76 mm	191735	407829	187417	356646	184393	186793
14 x 80 mm	191742	407836	187424	356653	184409	186809
*16 x 54 mm	191759	407843	187431	474395	184416	186816
*16 x 56 mm	191766	407850	187448	474401	184423	186823
16 x 60 mm	191773	407867	187455	356660	184430	186830
16 x 64 mm	191780	407874	187462	356677	184447	186847
*16 x 66 mm	404002	407881	404521	474418	405085	405245
16 x 68 mm	191797	407898	187479	356684	184454	186854
*16 x 70 mm	404019	407904	404538	474425	405092	405252
16 x 72 mm	191803	407911	187486	356691	184461	186861
16 x 76 mm	191810	407928	187493	356707	184478	186878
16 x 80 mm	191827	407935	187509	356714	184485	186885
*18 x 54 mm	191834	408222	187516	474432	184492	186892
*18 x 56 mm	191841	408239	187523	474449	184508	186908
18 x 60 mm	191858	408246	187530	356721	184515	186915
18 x 64 mm	191865	408253	187547	356738	184522	186922
*18 x 66 mm	404026	408260	404545	474456	405108	405269
18 x 68 mm	191872	408277	187554	356745	184539	186939
*18 x 70 mm	404033	408284	404552	474463	405115	405276
18 x 72 mm	191889	408291	187561	356752	184546	186946
18 x 76 mm	191896	408307	187578	356769	184553	186953
18 x 80 mm	191902	408314	187585	356776	184560	186960

TPM8 INSERT CRUCIATE RETAINING (CR) TOP FITTINGS, UHMWPE



ORDER INFORMATION 8698791 ...

Optional *

Size	CR Mobile Bearing UHMWPE	CR Mobile Bearing High Crosslink	CR Fixed Bearing UHMWPE	CR Fixed Bearing High Crosslink	CR Mobile UHMWPE Deep Dish	CR Mobile High Crosslink UHMWPE Deep Dish
* 8 x 54 mm	269359	359531	268871	359470	270799	271437
* 8 x 56 mm	269366	190974	268888	474470	270805	271444
8 x 60 mm	269373	358220	268895	357865	270812	271451
8 x 64 mm	269380	358237	268901	357872	270829	271468
* 8 x 66 mm	405320	405566	405924	474487	406402	406563
8 x 68 mm	269397	358244	268918	357889	270836	271475
* 8 x 70 mm	405337	405573	405931	474494	406419	406570
8 x 72 mm	269403	358251	268925	357896	270843	271482
8 x 76 mm	269410	358268	268932	357902	270850	271499
8 x 80 mm	269427	358275	268949	357919	270867	271505
* 10 x 54 mm	269434	359487	268956	359425	270874	271512
* 10 x 56 mm	269441	190981	268963	474500	270881	271529
10 x 60 mm	269458	181507	268970	357568	270898	271536
10 x 64 mm	269465	181514	268987	357575	270904	271543
* 10 x 66 mm	405344	405580	405948	474517	406426	406587
10 x 68 mm	269472	181521	268994	357582	270911	271550
* 10 x 70 mm	405351	405597	405955	474524	406433	406594
10 x 72 mm	269489	181538	269007	357599	270928	271567
10 x 76 mm	269496	357964	269014	357605	270935	271574
10 x 80 mm	269502	357971	269021	357612	270942	271581

TPM8 INSERT CRUCIATE RETAINING (CR) TOP FITTINGS, UHMWPE

ORDER INFORMATION 8698791 ...
Optional *

Size	CR Mobile Bearing UHMWPE	CR Mobile Bearing High Crosslink	CR Fixed Bearing UHMWPE	CR Fixed Bearing High Crosslink	CR Mobile UHMWPE Deep Dish	CR Mobile High Crosslink UHMWPE Deep Dish
*12 x 54 mm	269519	359494	269038	359432	270959	271598
*12 x 56 mm	269526	190998	269045	474531	270966	271604
12 x 60 mm	269533	181545	269052	357629	270973	271611
12 x 64 mm	269540	181552	269069	357636	270980	271628
*12 x 66 mm	405368	405603	405962	474548	406440	406600
12 x 68 mm	269557	358008	269076	357643	270997	271635
*12 x 70 mm	405375	405610	405979	474555	406457	406617
12 x 72 mm	269564	358015	269083	357650	271000	271642
12 x 76 mm	269571	358022	269090	357667	271017	271659
12 x 80 mm	269588	358039	269106	357674	271024	271666
*14 x 54 mm	269595	359500	269113	359449	271031	271673
*14 x 56 mm	269601	191001	269120	474562	271048	271680
14 x 60 mm	269618	358046	269137	357681	271055	271697
14 x 64 mm	269625	358053	269144	357698	271062	271703
*14 x 66 mm	405382	405627	405986	474579	406464	406624
14 x 68 mm	269632	358060	269151	357704	271079	271710
*14 x 70 mm	405399	405634	405993	474586	406471	406631
14 x 72 mm	269649	358077	269168	357711	271086	271727
14 x 76 mm	269656	358084	269175	357728	271093	271734
14 x 80 mm	269663	358091	269182	357735	271109	271741
*16 x 54 mm	269670	359517	269199	359456	271116	271758
*16 x 56 mm	269687	191018	269205	474593	271123	271765
16 x 60 mm	269694	358107	269212	357742	271130	271772
16 x 64 mm	269700	358114	269229	357759	271147	271789
*16 x 66 mm	405405	405641	406006	474609	406488	406648
16 x 68 mm	269717	358121	269236	357766	271154	271796
*16 x 70 mm	405412	405658	406013	474616	406495	406655
16 x 72 mm	269724	358138	269243	357773	271161	271802
16 x 76 mm	269731	358145	269250	357780	271178	271819
16 x 80 mm	269748	358152	269267	357797	271185	271826
*18 x 54 mm	269755	359524	269274	359463	271192	271833
*18 x 56 mm	269762	191025	269281	474623	271208	271840
18 x 60 mm	269779	358169	269298	357803	271215	271857
18 x 64 mm	269786	358176	269304	357810	271222	271864
*18 x 66 mm	405429	405665	406020	474630	406501	406662
18 x 68 mm	269793	358183	269311	357827	271239	271871
*18 x 70 mm	405436	405672	406037	474647	406518	406679
18 x 72 mm	269809	358190	269328	357834	271246	271888
18 x 76 mm	269816	358206	269335	357841	271253	271895
18 x 80 mm	269823	358213	269342	357858	271260	271901

TPM8 INSERT POSTERIOR STABILIZED (PS) TOP FITTINGS, UHMWPE



ORDER INFORMATION 8698791 ...

Optional *

Size	PS Mobile Bearing UHMWPE	PS Mobile Bearing High Crosslink	PS Fixed Bearing UHMWPE	PS Fixed Bearing High Crosslink	PS Mobile UHMWPE Deep Dish	PS Mobile High Crosslink UHMWPE Deep Dish
* 8 x 54 mm	272557	408321	272076	474654	274001	274643
* 8 x 56 mm	272564	408338	272083	474661	274018	274650
8 x 60 mm	272571	408345	272090	358947	274025	274667
8 x 64 mm	272588	408352	272106	358954	274032	274674
* 8 x 66 mm	405689	408802	406167	474678	406723	406884
8 x 68 mm	272595	408369	272113	358961	274049	274681
* 8 x 70 mm	405696	408819	406174	474685	406730	406891
8 x 72 mm	272601	408376	272120	358978	274056	274698
8 x 76 mm	272618	408383	272137	358985	274063	274704
8 x 80 mm	272625	408390	272144	358992	274070	274711
* 10 x 54 mm	272632	408406	272151	474692	274087	274728
* 10 x 56 mm	272649	408413	272168	474708	274094	274735
10 x 60 mm	272656	408420	272175	358640	274100	274742
10 x 64 mm	272663	408437	272182	358657	274117	274759
* 10 x 66 mm	405702	408826	406181	474715	406747	406907
10 x 68 mm	272670	408444	272199	358664	274124	274766
* 10 x 70 mm	405719	408833	406198	474722	406754	406914
10 x 72 mm	272687	408451	272205	358671	274131	274773
10 x 76 mm	272694	408468	272212	358688	274148	274780
10 x 80 mm	272700	408475	272229	358695	274155	274797

TPM8 INSERT POSTERIOR STABILIZED (PS) TOP FITTINGS, UHMWPE

ORDER INFORMATION 8698791 ...
Optional *

Size	PS Mobile Bearing UHMWPE	PS Mobile Bearing High Crosslink	PS Fixed Bearing UHMWPE	PS Fixed Bearing High Crosslink	PS Mobile UHMWPE Deep Dish	PS Mobile High Crosslink UHMWPE Deep Dish
*12 x 54 mm	272717	408482	272236	474739	274162	274803
*12 x 56 mm	272724	408499	272243	474746	274179	274810
12 x 60 mm	272731	408505	272250	358701	274186	274827
12 x 64 mm	272748	408512	272267	358718	274193	274834
*12 x 66 mm	405726	408840	406204	474753	406761	406921
12 x 68 mm	272755	408529	272274	358725	274209	274841
*12 x 70 mm	405733	408857	406211	474760	406778	406938
12 x 72 mm	272762	408536	272281	358732	274216	274858
12 x 76 mm	272779	408543	272298	358749	274223	274865
12 x 80 mm	272786	408550	272304	358756	274230	274872
*14 x 54 mm	272793	408567	272311	474777	274247	274889
*14 x 56 mm	272809	408574	272328	474784	274254	274896
14 x 60 mm	272816	408581	272335	358763	274261	274902
14 x 64 mm	272823	408598	272342	358770	274278	274919
*14 x 66 mm	405740	408864	406228	474791	406785	406945
14 x 68 mm	272830	408604	272359	358787	274285	274926
*14 x 70 mm	405757	408871	406235	474807	406792	406952
14 x 72 mm	272847	408611	272366	358794	274292	274933
14 x 76 mm	272854	408628	272373	358800	274308	274940
14 x 80 mm	272878	408635	272380	358817	274315	274957
*16 x 54 mm	272885	408642	272397	474814	274322	274964
*16 x 56 mm	272892	408659	272403	474821	274339	274971
16 x 60 mm	272908	408666	272410	358824	274346	274988
16 x 64 mm	272915	408673	272427	358831	274353	274995
*16 x 66 mm	405764	408888	406242	474838	406808	406969
16 x 68 mm	272922	408680	272434	358848	274360	275008
*16 x 70 mm	405771	408895	406259	474845	406815	406976
16 x 72 mm	272939	408697	272441	358855	274377	275015
16 x 76 mm	272946	408703	272458	358862	274384	275022
16 x 80 mm	272953	408710	272465	358879	274391	275039
*18 x 54 mm	272960	408727	272472	474852	274407	275046
*18 x 56 mm	272977	408734	272489	474869	274414	275053
18 x 60 mm	272984	408741	272496	358886	274421	275060
18 x 64 mm	272991	408758	272502	358893	274438	275077
*18 x 66 mm	405788	408901	406266	474876	406822	406983
18 x 68 mm	273004	408765	272519	358909	274445	275084
*18 x 70 mm	405795	408918	406273	474883	406839	406990
18 x 72 mm	273011	408772	272526	358916	274452	275091
18 x 76 mm	273028	408789	272533	358923	274469	275107
18 x 80 mm	273035	408796	272540	358930	274476	275114

PATELLA COMPONENT



ORDER INFORMATION 8698791 ...

Size	Type	Patella
28 mm	3 Peg Flat	188353
31 mm	3 Peg Flat	188360
34 mm	3 Peg Flat	188377
37 mm	3 Peg Flat	188384

TIBIAL METAL BACK STEM



ORDER INFORMATION 8698791 ...

Size	Tibial Metal Back Stem Grid Blast
10 mm	180210
20 mm	180227
30 mm	180234
40 mm	180241
50 mm	180258
60 mm	180265
70 mm	180272
80 mm	180289
90 mm	180296
100 mm	180302

TPM8 KNEE INSTRUMENTATION

— Trial Tibial
CR Insert B23

— Trial Tibial Stem
Punch Guide B22

— Trial Femoral
Component B20

— Trial Tibial Punch B22

TPM8
Primary
KNEE

FEMORAL CRUCIATE RETAINING INSTRUMENT SET 90118011001

FEMUR CONTAINER NO 1 - 46453024075

530 x 240 x 75 mm




Item No.	Description	Unit
*Not installed in tray. Only available on request.		
214557	*Femoral Trial Component Right 54 mm	1
214564	*Femoral Trial Component Right 56 mm	1
214571	Femoral Trial Component Right 60 mm	1
214588	Femoral Trial Component Right 64 mm	1
456841	*Femoral Trial Component Right 66 mm	1
214595	Femoral Trial Component Right 68 mm	1
456858	*Femoral Trial Component Right 70 mm	1
214601	Femoral Trial Component Right 72 mm	1
214618	Femoral Trial Component Right 76 mm	1
214625	Femoral Trial Component Right 80 mm	1
214632	*Femoral Trial Component Left 54 mm	1
214649	*Femoral Trial Component Left 56 mm	1
214656	Femoral Trial Component Left 60 mm	1
214663	Femoral Trial Component Left 64 mm	1
456865	*Femoral Trial Component Left 66 mm	1
214670	Femoral Trial Component Left 68 mm	1
456872	*Femoral Trial Component Left 70 mm	1
214687	Femoral Trial Component Left 72 mm	1
214694	Femoral Trial Component Left 76 mm	1
214700	Femoral Trial Component Left 80 mm	1
214793	Femoral Impactor, Large	1
214809	Femoral Size Caliper Pin, S Steel	1
214786	Femoral Impactor, Small	1
214779	Femur Size Calliper Pin Guide, 0°, S Steel	1
214717	Femur Size Calliper Pin Guide, Right, 3°, S Steel	1
214724	Femur Size Calliper Pin Guide, Right, 5°, S Steel	1
214748	Femur Size Calliper Pin Guide, Left, 3°, S Steel	1
214755	Femur Size Calliper Pin Guide, Left, 5°, S Steel	1
214816	Femoral Size Caliper, S Steel	1



FEMORAL CRUCIATE RETAINING INSTRUMENT SET 90118011001

FEMUR CONTAINER NO 2 - 46453024065

530 x 240 x 65 mm




Item No.	Description	Unit
*Not installed in tray. Only available on request.		
214823	*Femoral Distal Chamfer Cutting Guide 54mm	1
214830	*Femoral Distal Chamfer Cutting Guide 56mm	1
214847	Femoral Distal Chamfer Cutting Guide 60mm	1
214854	Femoral Distal Chamfer Cutting Guide 64mm	1
456889	*Femoral Distal Chamfer Cutting Guide 66mm	1
214861	Femoral Distal Chamfer Cutting Guide 68mm	1
456896	*Femoral Distal Chamfer Cutting Guide 70mm	1
214878	Femoral Distal Chamfer Cutting Guide 72mm	1
214885	Femoral Distal Chamfer Cutting Guide 76mm	1
214892	Femoral Distal Chamfer Cutting Guide 80mm	1
215257	Patella Holder Cutting Guide, S Steel	1
215271	Patella Trial, 28, Plastic	1
215288	Patella Trial, 31, Plastic	1
215295	Patella Trial, 34, Plastic	1
215301	Patella Trial, 37, Plastic	1
215325	Patella Drill Guide, 28-31, S Steel	1
215332	Patella Drill Guide, 34-37, S Steel	1
215042	Screw Driver, AA5	1
214922	Femur Distal Angle Cutting Equipment, 0,3,5,7,9, S Steel	1
214915	Femur Distal Cutting Guide, S Steel	1
309062	Femur Distal Angle Cutting Equipment Attachment	1
215318	Patella Drill, 6, S Steel	1
214977	Tighting Pin	1
214960	Femur Distal Chamfer Cutting Guide Pin, S Steel	2
214946	Femoral Medullary Drill 9 mm, S Steel	1
214939	Anchorage Peg Drill, 8, S Steel	1
215424	Patellar Clamp	1



FEMORAL CRUCIATE RETAINING INSTRUMENT SET 90118011001


TIBIA CONTAINER NO 3 - 46453024055



Item No.	Description	Unit
*Not installed in tray. Only available on request.		
239918	Gap Gauge, 8, Aluminum	1
214502	Gap Gauge, 10, Aluminum	1
214519	Gap Gauge, 12, Aluminum	1
214526	Gap Gauge, 14, Aluminum	1
214533	Gap Gauge, 16, Aluminum	1
214540	Gap Gauge, 18, Aluminum	1
215141	Trial Tibia Insert Holder, Titanium	2
214380	*Tibia Stem Punch Guide, 54, S Steel	1
214397	*Tibia Stem Punch Guide, 56, S Steel	1
214403	Tibia Stem Punch Guide, 60, S Steel	1
214410	Tibia Stem Punch Guide, 64, S Steel	1
456971	*Tibia Stem Punch Guide, 66, S Steel	1
214427	Tibia Stem Punch Guide, 68, S Steel	1
456988	*Tibia Stem Punch Guide, 70, S Steel	1
214434	Tibia Stem Punch Guide, 72, S Steel	1
214441	Tibia Stem Punch Guide, 76, S Steel	1
214458	Tibia Stem Punch Guide, 80, S Steel	1
215110	Tibial Impactor	1
214465	*Tibia Punch, 54 - 56, CoCrMo	1
214472	Tibia Punch, 60 - 64, CoCrMo	1
214489	Tibia Punch, 68 - 72, CoCrMo	1
214496	Tibia Punch, 76 - 80, CoCrMo	1
215028	Tibial Medullary Drill Guide	1
215233	Tibial Medullary Drill Ø12, S Steel	1
215134	Tibia Punch Impactor, S Steel	1
085065	Drill with Round Shaft, S Steel, 3.2*130	2
244110	T - Handle Type1, Short, S Steel	1



FEMORAL CRUCIATE RETAINING INSTRUMENT SET 90118011001



St. Steel Container
Bottom & Lid Perforated
530 x 240 x 55 mm

Item No.	Description	Unit
----------	-------------	------

218487	*Trial Tibial CR Insert, 54*8 mm, Plastic	1
215714	*Trial Tibial CR Insert, 54*10 mm, Plastic	1
215721	*Trial Tibial CR Insert, 54*12 mm, Plastic	1
215738	*Trial Tibial CR Insert, 54*14 mm, Plastic	1
215745	*Trial Tibial CR Insert, 54*16 mm, Plastic	1
215752	*Trial Tibial CR Insert, 54*18 mm, Plastic	1

*Not installed in tray. Only available on request.

218494	*Trial Tibial CR Insert, 56*8 mm, Plastic	1
215769	*Trial Tibial CR Insert, 56*10 mm, Plastic	1
215776	*Trial Tibial CR Insert, 56*12 mm, Plastic	1
215783	*Trial Tibial CR Insert, 56*14 mm, Plastic	1
215790	*Trial Tibial CR Insert, 56*16 mm, Plastic	1
215806	*Trial Tibial CR Insert, 56*18 mm, Plastic	1

*Not installed in tray. Only available on request.

218500	Trial Tibial CR Insert, 60*8 mm, Plastic	1
215813	Trial Tibial CR Insert, 60*10 mm, Plastic	1
215820	Trial Tibial CR Insert, 60*12 mm, Plastic	1
215837	Trial Tibial CR Insert, 60*14 mm, Plastic	1
215844	Trial Tibial CR Insert, 60*16 mm, Plastic	1
215851	Trial Tibial CR Insert, 60*18 mm, Plastic	1

218517	*Trial Tibial CR Insert, 64*8 mm, Plastic	1
215868	*Trial Tibial CR Insert, 64*10 mm, Plastic	1
215875	*Trial Tibial CR Insert, 64*12 mm, Plastic	1
215882	*Trial Tibial CR Insert, 64*14 mm, Plastic	1
215899	*Trial Tibial CR Insert, 64*16 mm, Plastic	1
215905	*Trial Tibial CR Insert, 64*18 mm, Plastic	1

*Not installed in tray. Only available on request.

Item No.	Description	Unit
----------	-------------	------

218524	Trial Tibial CR Insert, 68*8 mm, Plastic	1
215912	Trial Tibial CR Insert, 68*10 mm, Plastic	1
215929	Trial Tibial CR Insert, 68*12 mm, Plastic	1
215936	Trial Tibial CR Insert, 68*14 mm, Plastic	1
215943	Trial Tibial CR Insert, 68*16 mm, Plastic	1
215950	Trial Tibial CR Insert, 68*18 mm, Plastic	1

218531	*Trial Tibial CR Insert, 72*8 mm, Plastic	1
215967	*Trial Tibial CR Insert, 72*10 mm, Plastic	1
215974	*Trial Tibial CR Insert, 72*12 mm, Plastic	1
215981	*Trial Tibial CR Insert, 72*14 mm, Plastic	1
215998	*Trial Tibial CR Insert, 72*16 mm, Plastic	1
216001	*Trial Tibial CR Insert, 72*18 mm, Plastic	1

*Not installed in tray. Only available on request.


218548	Trial Tibial CR Insert, 76*8 mm, Plastic	1
216018	Trial Tibial CR Insert, 76*10 mm, Plastic	1
216025	Trial Tibial CR Insert, 76*12 mm, Plastic	1
216032	Trial Tibial CR Insert, 76*14 mm, Plastic	1
216049	Trial Tibial CR Insert, 76*16 mm, Plastic	1
216056	Trial Tibial CR Insert, 76*18 mm, Plastic	1

218555	*Trial Tibial CR Insert, 80*8 mm, Plastic	1
216063	*Trial Tibial CR Insert, 80*10 mm, Plastic	1
216070	*Trial Tibial CR Insert, 80*12 mm, Plastic	1
216087	*Trial Tibial CR Insert, 80*14 mm, Plastic	1
216094	*Trial Tibial CR Insert, 80*16 mm, Plastic	1
216100	*Trial Tibial CR Insert, 80*18 mm, Plastic	1

*Not installed in tray. Only available on request.

FEMORAL CRUCIATE RETAINING INSTRUMENT SET 90118011001


TIBIA CONTAINER NO 4 - 46453024155



Item No.	Description	Unit
214984	Feeler Blade, Small, S Steel	1
214991	Feeler Blade, Long, S Steel	1
215011	Tibia Stylus, S Steel	1
215073	External Guide, 6, Titanium	1
215059	Medula Guide, 9, Short, Titanium	1
215066	Medula Guide, 9, Long, Titanium	1
215240	Pin-Punch Impactor-Extractor	1
215158	Tibia Stopper Screw Driver	1
215189	Pin With Head, 25, S Steel	4
215196	Pin With Head, 35, S Steel	4
215202	Pin With Head, 45, S Steel	4
215219	Pin With Head, 55, S Steel	4
215226	Pin With Head, 65, S Steel	4
215165	Pin Straight, 55, S Steel	4
215172	Pin Straight, 65, S Steel	4
215127	Pin Extractor, S Steel	1



FEMORAL CRUCIATE RETAINING INSTRUMENT SET 90118011001



St. Steel Container Bot-
tom & Lid Perforated
530 x 240 x 55 mm

Item No.	Description	Unit
215080	Tibia Proximal Cutting Guide, Right, S Steel	1
215103	Tibia Proximal Cutting Guide, Left, S Steel	1
215004	Intramedullary Tibia Proximal Cutting Guide Holder, S Steel	1
215097	Tibia Proximal Cutting Guide, Straight, S Steel	1
215035	Extramedullary Tibia Proximal Cutting Guide Holder, S Steel	1
246626	Extramedullary Proximal Tibial Cutting Guide Component, S Steel	1



* Assembled

FEMORAL POSTERIOR STABILIZED INSTRUMENT SET 90120101001

FEMUR TRAY NO 1 - 80120101001



Item No.	Description	Unit
*Not installed in tray. Only available on request.		
215554	*Femoral PS Trial Component, 54 mm, Right, CoCrMo	1
215561	*Femoral PS Trial Component, 56 mm, Right, CoCrMo	1
215578	Femoral PS Trial Component, 60 mm, Right, CoCrMo	1
215585	Femoral PS Trial Component, 64 mm, Right, CoCrMo	1
456902	*Femoral PS Trial Component, 66 mm, Right, CoCrMo	1
215592	Femoral PS Trial Component, 68 mm, Right, CoCrMo	1
456919	*Femoral PS Trial Component, 70 mm, Right, CoCrMo	1
215608	Femoral PS Trial Component, 72 mm, Right, CoCrMo	1
215615	Femoral PS Trial Component, 76 mm, Right, CoCrMo	1
215622	Femoral PS Trial Component, 80 mm, Right, CoCrMo	1
215639	*Femoral PS Trial Component, 54 mm, Left, CoCrMo	1
215646	*Femoral PS Trial Component, 56 mm, Left, CoCrMo	1
215653	Femoral PS Trial Component, 60 mm, Left, CoCrMo	1
215660	Femoral PS Trial Component, 64 mm, Left, CoCrMo	1
456926	*Femoral PS Trial Component, 66 mm, Left, CoCrMo	1
215677	Femoral PS Trial Component, 68 mm, Left, CoCrMo	1
456933	*Femoral PS Trial Component, 70 mm, Left, CoCrMo	1
215684	Femoral PS Trial Component, 72 mm, Left, CoCrMo	1
215691	Femoral PS Trial Component, 76 mm, Left, CoCrMo	1
215707	Femoral PS Trial Component, 80 mm, Left, CoCrMo	1
215431	*Femoral PS Distal Guide, 54 mm, S Steel	1
215448	*Femoral PS Distal Guide, 56 mm, S Steel	1
215455	Femoral PS Distal Guide, 60 mm, S Steel	1
215462	Femoral PS Distal Guide, 64 mm, S Steel	1
456940	*Femoral PS Distal Guide, 66 mm, S Steel	1
215479	Femoral PS Distal Guide, 68 mm, S Steel	1
456957	*Femoral PS Distal Guide, 70 mm, S Steel	1
215486	Femoral PS Distal Guide, 72 mm, S Steel	1
215493	Femoral PS Distal Guide, 76 mm, S Steel	1
215509	Femoral PS Distal Guide, 80 mm, S Steel	1



FEMORAL POSTERIOR STABILIZED INSTRUMENT SET 90120101001

FEMUR TRAY NO 2 - 80120101002



46426551160
St. Steel Container
Bottom & Lid Perforated
550 x 265 x 160 mm



Item No.	Description	Unit
215516	*Femoral PS Distal Chisel, 54-56 mm, S Steel	1
	*Not installed in tray. Only available on request.	
215523	Femoral PS Distal Chisel, 60-64 mm, S Steel	1
215530	Femoral PS Distal Chisel, 68-72 mm, S Steel	1
215547	Femoral PS Distal Chisel, 76-80 mm, S Steel	1



FEMORAL POSTERIOR STABILIZED INSTRUMENT SET 90120101001



Item No.	Description	Unit
218562	*Trial Tibial Insert PS, 54*8 mm, Plastic	1
213581	*Trial Tibial Insert PS, 54*10 mm, Plastic	1
213598	*Trial Tibial Insert PS, 54*12 mm, Plastic	1
213604	*Trial Tibial Insert PS, 54*14 mm, Plastic	1
213611	*Trial Tibial Insert PS, 54*16 mm, Plastic	1
213628	*Trial Tibial Insert PS, 54*18 mm, Plastic	1
*Not installed in tray. Only available on request.		
218579	*Trial Tibial Insert PS, 56*8 mm, Plastic	1
213635	*Trial Tibial Insert PS, 56*10 mm, Plastic	1
213642	*Trial Tibial Insert PS, 56*12 mm, Plastic	1
213659	*Trial Tibial Insert PS, 56*14 mm, Plastic	1
213666	*Trial Tibial Insert PS, 56*16 mm, Plastic	1
213673	*Trial Tibial Insert PS, 56*18 mm, Plastic	1
*Not installed in tray. Only available on request.		
218586	Trial Tibial Insert PS, 60*8 mm, Plastic	1
213680	Trial Tibial Insert PS, 60*10 mm, Plastic	1
213697	Trial Tibial Insert PS, 60*12 mm, Plastic	1
213703	Trial Tibial Insert PS, 60*14 mm, Plastic	1
213710	Trial Tibial Insert PS, 60*16 mm, Plastic	1
213727	Trial Tibial Insert PS, 60*18 mm, Plastic	1
218593	Trial Tibial Insert PS, 64*8 mm, Plastic	1
213734	Trial Tibial Insert PS, 64*10 mm, Plastic	1
213741	Trial Tibial Insert PS, 64*12 mm, Plastic	1
213758	Trial Tibial Insert PS, 64*14 mm, Plastic	1
213765	Trial Tibial Insert PS, 64*16 mm, Plastic	1
213772	Trial Tibial Insert PS, 64*18 mm, Plastic	1

Item No.	Description	Unit
218609	Trial Tibial Insert PS, 68*8 mm, Plastic	1
213789	Trial Tibial Insert PS, 68*10 mm, Plastic	1
213796	Trial Tibial Insert PS, 68*12 mm, Plastic	1
213802	Trial Tibial Insert PS, 68*14 mm, Plastic	1
213819	Trial Tibial Insert PS, 68*16 mm, Plastic	1
213826	Trial Tibial Insert PS, 68*18 mm, Plastic	1
218616	Trial Tibial Insert PS, 72*8 mm, Plastic	1
213833	Trial Tibial Insert PS, 72*10 mm, Plastic	1
213840	Trial Tibial Insert PS, 72*12 mm, Plastic	1
213857	Trial Tibial Insert PS, 72*14 mm, Plastic	1
213864	Trial Tibial Insert PS, 72*16 mm, Plastic	1
213871	Trial Tibial Insert PS, 72*18 mm, Plastic	1
218623	Trial Tibial Insert PS, 76*8 mm, Plastic	1
213888	Trial Tibial Insert PS, 76*10 mm, Plastic	1
213895	Trial Tibial Insert PS, 76*12 mm, Plastic	1
213901	Trial Tibial Insert PS, 76*14 mm, Plastic	1
213918	Trial Tibial Insert PS, 76*16 mm, Plastic	1
213925	Trial Tibial Insert PS, 76*18 mm, Plastic	1
218630	Trial Tibial Insert PS, 80*8 mm, Plastic	1
213932	Trial Tibial Insert PS, 80*10 mm, Plastic	1
213949	Trial Tibial Insert PS, 80*12 mm, Plastic	1
213956	Trial Tibial Insert PS, 80*14 mm, Plastic	1
213963	Trial Tibial Insert PS, 80*16 mm, Plastic	1
213970	Trial Tibial Insert PS, 80*18 mm, Plastic	1



TPM
Revision
KNEE

FEMORAL PS COMPONENTS



ORDER INFORMATION 8698791 ...

Optional *

Side	Size	Femoral Component Co.Cr.Mo.
Left	60 mm	230069
Left	64 mm	230076
Left	*66 mm	347521
Left	68 mm	230083
Left	*70 mm	347514
Left	72 mm	230090
Left	76 mm	230106
Left	80 mm	230113

Side	Size	Femoral Component Co.Cr.Mo.
Right	60 mm	230007
Right	64 mm	230014
Right	*66 mm	347507
Right	68 mm	230021
Right	*70 mm	347491
Right	72 mm	230038
Right	76 mm	230045
Right	80 mm	230052

POSTERIOR FEMORAL AUGMENT



ORDER INFORMATION 8698791 ...

Optional *


Size	Titanium	Size	Titanium
5 x 60 mm	353201	10 x 60 mm	353195
5 x 64 mm	353225	10 x 64 mm	353218
*5 x 66 mm	353249	*10 x 66 mm	353232
5 x 68 mm	353263	10 x 68 mm	353256
*5 x 70 mm	353287	*10 x 70 mm	353270
5 x 72 mm	353300	10 x 72 mm	353294
5 x 76 mm	353324	10 x 76 mm	353317
5 x 80 mm	353348	10 x 80 mm	353331



Fixation screw of femoral
posterior block augment
(M8)

350149

DISTAL FEMORAL AUGMENT



ORDER INFORMATION 8698791 ...

Optional *

Size	Right Lateral / Left Medial	Right Medial / Left Lateral	Size	Right Lateral / Left Medial	Right Medial / Left Lateral
5 x 60 mm	352242	352488	10 x 60 mm	352228	352464
5 x 64 mm	352273	352518	10 x 64 mm	352259	352495
*5 x 66 mm	352303	352549	*10 x 66 mm	352280	352525
5 x 68 mm	352334	352570	10 x 68 mm	352310	352556
*5 x 70 mm	352365	352600	*10 x 70 mm	352341	352587
5 x 72 mm	352396	352631	10 x 72 mm	352372	352617
5 x 76 mm	352426	352662	10 x 76 mm	352402	352648
5 x 80 mm	352457	352693	10 x 80 mm	352433	352679

Size	Right Lateral / Left Medial	Right Medial / Left Lateral
15 x 60 mm	352235	352471
15 x 64 mm	352266	352501
*15 x 66 mm	352297	352532
15 x 68 mm	352327	352563
*15 x 70 mm	352358	352594
15 x 72 mm	352389	352624
15 x 76 mm	352419	352655
15 x 80 mm	352440	352686



Fixation screw of femoral
posterior block augment
(M8)

350149

EXTENSION STEM



ORDER INFORMATION 8698791 ...

Stem Size	Stem Length	Titanium Grit - Blast	Stem Size	Stem Length	Titanium Grit - Blast	Stem Size	Stem Length	Titanium Grit - Blast
10 mm	40 mm	354642	14 mm	40 mm	354680	18 mm	40 mm	354727
10 mm	80 mm	354338	14 mm	80 mm	354468	18 mm	80 mm	354581
10 mm	120 mm	354291	14 mm	120 mm	354413	18 mm	120 mm	354536
10 mm	160 mm	354628	14 mm	160 mm	354420	18 mm	160 mm	354543

Stem Size	Stem Length	Titanium Grit - Blast	Stem Size	Stem Length	Titanium Grit - Blast	Stem Size	Stem Length	Titanium Grit - Blast
11 mm	40 mm	354659	15 mm	40 mm	354703	19 mm	40 mm	354307
11 mm	80 mm	354345	15 mm	80 mm	354499	19 mm	80 mm	354598
11 mm	120 mm	354314	15 mm	120 mm	354444	19 mm	120 mm	354550
11 mm	160 mm	354321	15 mm	160 mm	354451	19 mm	160 mm	354567

Stem Size	Stem Length	Titanium Grit - Blast	Stem Size	Stem Length	Titanium Grit - Blast	Stem Size	Stem Length	Titanium Grit - Blast
12 mm	40 mm	354666	16 mm	40 mm	354697	20 mm	40 mm	354390
12 mm	80 mm	354369	16 mm	80 mm	354505	20 mm	80 mm	354635
12 mm	120 mm	354383	16 mm	120 mm	354482	20 mm	120 mm	354604
12 mm	160 mm	354376	16 mm	160 mm	354475	20 mm	160 mm	354611

Stem Size	Stem Length	Titanium Grit - Blast	Stem Size	Stem Length	Titanium Grit - Blast
13 mm	40 mm	354673	17 mm	40 mm	354710
13 mm	80 mm	354437	17 mm	80 mm	354574
13 mm	120 mm	354406	17 mm	120 mm	354512
13 mm	160 mm	354352	17 mm	160 mm	354529

TIBIAL AUGMENT



ORDER INFORMATION 8698791 ... Optional *

Size	Right Lateral / Left Medial	Right Medial / Left Lateral
5 x 60 mm	353843	354079
5 x 64 mm	353874	354109
5 x 66 mm*	353904	354130
5 x 68 mm	353935	354161
5 x 70 mm*	353966	354192
5 x 72 mm	353997	354222
5 x 76 mm	367925	354253
5 x 80 mm	354048	354284



Tibial Augment Screw

Titanium

	349730
20	084068

ORDER INFORMATION 8698791 ... Optional *

Size	Right Lateral / Left Medial	Right Medial / Left Lateral
10 x 60 mm	353829	354055
10 x 64 mm	353850	354086
10 x 66 mm*	353881	354116
10 x 68 mm	353911	354147
10 x 70 mm*	353942	354178
10 x 72 mm	353973	354208
10 x 76 mm	354000	354239
10 x 80 mm	354024	354260

ORDER INFORMATION 8698791 ... Optional *

Size	Right Lateral / Left Medial	Right Medial / Left Lateral
15 x 60 mm	353836	354062
15 x 64 mm	353867	354093
15 x 66 mm*	353898	354123
15 x 68 mm	353928	354154
15 x 70 mm*	353959	354185
15 x 72 mm	353980	354215
15 x 76 mm	354017	354246
15 x 80 mm	354031	354277

TIBIAL METAL BACK



ORDER INFORMATION 8698791 ...

Optional *

Size	Tibial Metal Back Cemented
60 mm	230120
64 mm	230137
66 mm*	347484
68 mm	230144
70 mm*	347477
72 mm	230151
76 mm	230168
80 mm	230175

UHMWPE POSTERIOR STABILIZED (PS) TIBIAL INSERT

ORDER INFORMATION 8698791 ... Optional *

Size	UHMWPE PS Insert	High Crosslink PS Insert
10 x 60 mm	230182	491361
10 x 64 mm	230298	491378
10 x 66 mm*	347682	491385
10 x 68 mm	230403	491392
10 x 70 mm*	347675	491408
10 x 72 mm	230519	491415
10 x 76 mm	230625	491422
10 x 80 mm	230731	491439
12 x 60 mm	230199	491446
12 x 64 mm	230304	491453
12 x 66 mm*	347668	491460
12 x 68 mm	230410	491477
12 x 70 mm*	347651	491484
12 x 72 mm	230526	491491
12 x 76 mm	230632	491507
12 x 80 mm	230748	491514
14 x 60 mm	230205	491521
14 x 64 mm	230311	491538
14 x 66 mm*	347644	491545
14 x 68 mm	230427	491552
14 x 70 mm*	347637	491569
14 x 72 mm	230533	491576
14 x 76 mm	230649	491583
14 x 80 mm	230755	491590
16 x 60 mm	230212	491606
16 x 64 mm	230328	491613
16 x 66 mm*	347620	491620
16 x 68 mm	230434	491637
16 x 70 mm*	347613	491644
16 x 72 mm	230540	491651
16 x 76 mm	230656	491668
16 x 80 mm	230762	491675
18 x 60 mm	230229	491682
18 x 64 mm	230335	491699
18 x 66 mm*	347606	491705
18 x 68 mm	230441	491712
18 x 70 mm*	347590	491729
18 x 72 mm	230557	491736



ORDER INFORMATION 8698791 ... Optional *

Size	UHMWPE PS Insert	High Crosslink PS Insert
18 x 76 mm	230663	491743
18 x 80 mm	230779	491750
20 x 60 mm	230236	491767
20 x 64 mm	230342	491774
20 x 66 mm*	347583	491781
20 x 68 mm	230458	491798
20 x 70 mm*	347576	491804
20 x 72 mm	230564	491811
20 x 76 mm	230670	491828
20 x 80 mm	230786	491835
22 x 60 mm	230243	491842
22 x 64 mm	230359	491859
22 x 66 mm*	347569	491866
22 x 68 mm	230465	491873
22 x 70 mm*	347552	491880
22 x 72 mm	230571	491897
22 x 76 mm	230687	491903
22 x 80 mm	230793	491910
24 x 60 mm	230250	491927
24 x 64 mm	230366	491934
24 x 66 mm*	347545	491941
24 x 68 mm	230472	491958
24 x 70 mm*	347538	491965
24 x 72 mm	230588	491972
24 x 76 mm	230694	491989
24 x 80 mm	230809	491996

UHMWPE PLUS TIBIAL INSERT

ORDER INFORMATION 8698791 ... Optional *

Size	UHMWPE PLUS Insert	High Crosslink PLUS Insert
10 x 60 mm	492009	492641
10 x 64 mm	492016	492658
10 x 66 mm*	492023	492665
10 x 68 mm	492030	492672
10 x 70 mm*	492047	492689
10 x 72 mm	492054	492696
10 x 76 mm	492061	492702
10 x 80 mm	492078	492719
12 x 60 mm	492085	492726
12 x 64 mm	492092	492733
12 x 66 mm*	492108	492740
12 x 68 mm	492115	492757
12 x 70 mm*	492122	492764
12 x 72 mm	492139	492771
12 x 76 mm	492146	492788
12 x 80 mm	492153	492795
14 x 60 mm	492160	492801
14 x 64 mm	492177	492818
14 x 66 mm*	492184	492825
14 x 68 mm	492191	492832
14 x 70 mm*	492207	492849
14 x 72 mm	492214	492856
14 x 76 mm	492221	492863
14 x 80 mm	492238	492870
16 x 60 mm	492245	492887
16 x 64 mm	492252	492894
16 x 66 mm*	492269	492900
16 x 68 mm	492276	492917
16 x 70 mm*	492283	492924
16 x 72 mm	492290	492931
16 x 76 mm	492306	492948
16 x 80 mm	492313	492955
18 x 60 mm	492320	492962
18 x 64 mm	492337	492979
18 x 66 mm*	492344	492986
18 x 68 mm	492351	492993
18 x 70 mm*	492368	493006
18 x 72 mm	492375	493013



ORDER INFORMATION 8698791 ... Optional *

Size	UHMWPE PLUS Insert	High Crosslink PLUS Insert
18 x 76 mm	492382	493020
18 x 80 mm	492399	493037
20 x 60 mm	492405	493044
20 x 64 mm	492412	493051
20 x 66 mm*	492429	493068
20 x 68 mm	492436	493075
20 x 70 mm*	492443	493082
20 x 72 mm	492450	493099
20 x 76 mm	492467	493105
20 x 80 mm	492474	493112
22 x 60 mm	492481	493129
22 x 64 mm	492498	493136
22 x 66 mm*	492504	493143
22 x 68 mm	492511	493150
22 x 70 mm*	492528	493167
22 x 72 mm	492535	493174
22 x 76 mm	492542	493181
22 x 80 mm	492559	493198
24 x 60 mm	492566	493204
24 x 64 mm	492573	493211
24 x 66 mm*	492580	493228
24 x 68 mm	492597	493235
24 x 70 mm*	492603	493242
24 x 72 mm	492610	493259
24 x 76 mm	492627	493266
24 x 80 mm	492634	493273

OFFSET ADAPTER



Offset Stem Adapter
fixation screw to
tibial metal back

Titanium

350156



Offset Stem Adapter
fixation screw to
femoral component

Titanium

350163

ORDER INFORMATION 8698791 ...

Offset	Offset Adapter Titanium
0 mm	232827
2 mm	232834
4 mm	232841
6 mm	232858

PATELLA COMPONENT



ORDER INFORMATION 8698791 ...

Size	Type	Patella
28 mm	3 Peg Flat	188353
31 mm	3 Peg Flat	188360
34 mm	3 Peg Flat	188377
37 mm	3 Peg Flat	188384

TPM REVISION KNEE INSTRUMENTATION




TPM
Revision
KNEE

TPM REVISION KNEE SYSTEM FEMORAL TEMPLATE TRAY

CONTAINER - TRAY - 90126101005



*
Femoral Distal
Cutting Guide




Item No.	Description	Unit
236283	Distal Cut Block, S Steel	1 *
351061	Femur Distal Cutting Device, 0°*6°, S Steel	1 *
351481	Distal Block Attachment Tower, S Steel	1 *
351092	Distal Bond Cut Chisel, 60*64, S Steel	1
351115	Distal Bond Cut Chisel, 68*72, S Steel	1
351122	Distal Bond Cut Chisel, 76*80, S Steel	1
236221	Femoral Template, 60, S Steel	1
236238	Femoral Template, 64, S Steel	1
236245	Femoral Template, 68, S Steel	1
236252	Femoral Template, 72, S Steel	1
236269	Femoral Template, 76, S Steel	1
236276	Femoral Template, 80, S Steel	1



TPM REVISION KNEE SYSTEM FEMORAL TEMPLATE TRAY

St. Steel Container
Lid Perforated
530 x 240 x 70 mm



Item No.	Description	Unit
351382	Femoral Box Resection Gauge	1
308423	Wrench, 15 mm	1
351085	Femoral Offset Coin Rotation Handle	1
214786	Femoral Impactor, Small	1
351146	Tighting Pin, S Steel	1
367826	Extractor Chisel, S Steel	1
368472	Femoral Impactor, Large	1



TPM REVISION KNEE SYSTEM FEMORAL OFFSET CAPABLE CUT BLOCK TRAY


CONTAINER - TRAY - 90126101002



Item No.	Description	Unit
214960	Femur Distal Chamfer Cutting Guide Pin, S Steel	2
351146	Tighting Pin, S Steel	1
236481	Femoral Box Resection Guide, Left, S Steel	1
236474	Femoral Box Resection Guide, Right, S Steel	1
351498	Femoral Box Chisel Guide, S Steel	1
236504	Femoral Reaming Housing, Left, S Steel	1
236498	Femoral Reaming Housing, Right, S Steel	1
236597	Femoral Boss Reamer Bushing, 0, S Steel	1
236603	Femoral Offset Reamer Bushing, 2, S Steel	1
236610	Femoral Offset Reamer Bushing, 4, S Steel	1
236627	Femoral Offset Reamer Bushing, 6, S Steel	1
236511	Femoral Offset Reamer, 0, S Steel	1
236528	Femoral Offset Reamer, 2, S Steel	1
236535	Femoral Offset Reamer, 4, S Steel	1
236542	Femoral Offset Reamer, 6, S Steel	1



TPM REVISION KNEE SYSTEM FEMORAL OFFSET CAPABLE CUT BLOCK TRAY



Item No.	Description	Unit
236351	Femoral Offset Capable Cut Block, Left, 60, S Steel	1
236368	Femoral Offset Capable Cut Block, Left, 64, S Steel	1
236375	Femoral Offset Capable Cut Block, Left, 68, S Steel	1
236382	Femoral Offset Capable Cut Block, Left, 72, S Steel	1
236399	Femoral Offset Capable Cut Block, Left, 76, S Steel	1
236405	Femoral Offset Capable Cut Block, Left, 80, S Steel	1
236290	Femoral Offset Capable Cut Block, Right, 60, S Steel	1
236306	Femoral Offset Capable Cut Block, Right, 64, S Steel	1
236313	Femoral Offset Capable Cut Block, Right, 68, S Steel	1
236320	Femoral Offset Capable Cut Block, Right, 72, S Steel	1
236337	Femoral Offset Capable Cut Block, Right, 76, S Steel	1
236344	Femoral Offset Capable Cut Block, Right, 80, S Steel	1
236559	Femoral Offset Coin, 0, S Steel	1
236566	Femoral Offset Coin, 2, S Steel	1
236573	Femoral Offset Coin, 4, S Steel	1
236580	Femoral Offset Coin, 6, S Steel	1
236412	Femoral Cut Block Insert, 60, S Steel	1
236429	Femoral Cut Block Insert, 64, S Steel	1
236436	Femoral Cut Block Insert, 68, S Steel	1
236443	Femoral Cut Block Insert, 72, S Steel	1
236450	Femoral Cut Block Insert, 76, S Steel	1
236467	Femoral Cut Block Insert, 80, S Steel	1
350620	Trial Offset Adapter, 0, S Steel	1
350637	Trial Offset Adapter, 2, S Steel	1
350644	Trial Offset Adapter, 4, S Steel	1
350651	Trial Offset Adapter, 6, S Steel	1



TPM REVISION KNEE SYSTEM FEMORAL TRIALS TRAY

CONTAINER - TRAY - 90126101001



*
Femoral Augment
Resection Guide




Item No.	Description	Unit
236047	Femoral Cut Through Trial, Left, 60, S Steel	1
236054	Femoral Cut Through Trial, Left, 64, S Steel	1
236061	Femoral Cut Through Trial, Left, 68, S Steel	1
236078	Femoral Cut Through Trial, Left, 72, S Steel	1
236085	Femoral Cut Through Trial, Left, 76, S Steel	1
236092	Femoral Cut Through Trial, Left, 80, S Steel	1
235989	Femoral Cut Through Trial, Right, 60, S Steel	1
235996	Femoral Cut Through Trial, Right, 64, S Steel	1
236009	Femoral Cut Through Trial, Right, 68, S Steel	1
236016	Femoral Cut Through Trial, Right, 72, S Steel	1
236023	Femoral Cut Through Trial, Right, 76, S Steel	1
236030	Femoral Cut Through Trial, Right, 80, S Steel	1
350743	Femoral Augment Resection, 60, S Steel	1
350750	Femoral Augment Resection, 64, S Steel	1
350774	Femoral Augment Resection, 68, S Steel	1
350798	Femoral Augment Resection, 72, S Steel	1
350804	Femoral Augment Resection, 76, S Steel	1
350811	Femoral Augment Resection, 80, S Steel	1



TPM REVISION KNEE SYSTEM FEMORAL TRIALS TRAY

St. Steel Container
Lid Perforated
530 x 240 x 80 mm



Item No.	Description	Unit
236160	Femoral Modular Box, Left, 60, S Steel	1
236177	Femoral Modular Box, Left, 64, S Steel	1
236184	Femoral Modular Box, Left, 68, S Steel	1
236191	Femoral Modular Box, Left, 72, S Steel	1
236207	Femoral Modular Box, Left, 76, S Steel	1
236214	Femoral Modular Box, Left, 80, S Steel	1
236108	Femoral Modular Box, Right, 60, S Steel	1
236115	Femoral Modular Box, Right, 64, S Steel	1
236122	Femoral Modular Box, Right, 68, S Steel	1
236139	Femoral Modular Box, Right, 72, S Steel	1
236146	Femoral Modular Box, Right, 76, S Steel	1
236153	Femoral Modular Box, Right, 80, S Steel	1



TPM REVISION KNEE SYSTEM TRIALS FEMORAL AUGMENT TRAY

CONTAINER - TRAY - 90126101004



Item No.	Description	Unit
236634	Trial Posterior Femoral Augment , 60*5, Aluminum	2
236641	Trial Posterior Femoral Augment , 60*10, Aluminum	2
236658	Trial Posterior Femoral Augment , 64*5, Aluminum	2
236665	Trial Posterior Femoral Augment , 64*10, Aluminum	2
236672	Trial Posterior Femoral Augment , 68*5, Aluminum	2
236689	Trial Posterior Femoral Augment , 68*10, Aluminum	2
236696	Trial Posterior Femoral Augment , 72*5, Aluminum	2
236702	Trial Posterior Femoral Augment , 72*10, Aluminum	2
236719	Trial Posterior Femoral Augment , 76*5, Aluminum	2
236726	Trial Posterior Femoral Augment , 76*10, Aluminum	2
236733	Trial Posterior Femoral Augment , 80*5, Aluminum	2
236740	Trial Posterior Femoral Augment , 80*10, Aluminum	2



Trial Femoral
Posterior Augment

Item No.	Description	Unit
236764	Trial Distal Femoral Augment, Right-Medial, 60*5, Aluminum	1
236801	Trial Distal Femoral Augment, Right-Medial, 60*10, Aluminum	1
236849	Trial Distal Femoral Augment, Right-Medial, 60*15, Aluminum	1
236887	Trial Distal Femoral Augment, Right-Medial, 64*5, Aluminum	1
236924	Trial Distal Femoral Augment, Right-Medial, 64*10, Aluminum	1
236962	Trial Distal Femoral Augment, Right-Medial, 64*15, Aluminum	1
236757	Trial Distal Femoral Augment, Right-Lateral, 60*5, Aluminum	1
236795	Trial Distal Femoral Augment, Right-Lateral, 60*10, Aluminum	1
236832	Trial Distal Femoral Augment, Right-Lateral, 60*15, Aluminum	1
236870	Trial Distal Femoral Augment, Right-Lateral, 64*5, Aluminum	1
236917	Trial Distal Femoral Augment, Right-Lateral, 64*10, Aluminum	1
236955	Trial Distal Femoral Augment, Right-Lateral, 64*15, Aluminum	1



Trial Femoral Distal
Augment

TPM REVISION KNEE SYSTEM TRIALS FEMORAL AUGMENT TRAY


Item No.	Description	Unit
236788	Trial Distal Femoral Augment, Left-Medial, 60*5, Aluminum	1
236825	Trial Distal Femoral Augment, Left-Medial, 60*10, Aluminum	1
236863	Trial Distal Femoral Augment, Left-Medial, 60*15, Aluminum	1
236900	Trial Distal Femoral Augment, Left-Medial, 64*5, Aluminum	1
236948	Trial Distal Femoral Augment, Left-Medial, 64*10, Aluminum	1
236986	Trial Distal Femoral Augment, Left-Medial, 64*15, Aluminum	1
236771	Trial Distal Femoral Augment, Left-Lateral, 60*5, Aluminum	1
236818	Trial Distal Femoral Augment, Left-Lateral, 60*10, Aluminum	1
236856	Trial Distal Femoral Augment, Left-Lateral, 60*15, Aluminum	1
236894	Trial Distal Femoral Augment, Left-Lateral, 64*5, Aluminum	1
236931	Trial Distal Femoral Augment, Left-Lateral, 64*10, Aluminum	1
236979	Trial Distal Femoral Augment, Left-Lateral, 64*15, Aluminum	1
237006	Trial Distal Femoral Augment, Right-Medial / Left-Lateral, 68*5, Aluminum	1
237020	Trial Distal Femoral Augment, Right-Medial / Left-Lateral, 68*10, Aluminum	1
237044	Trial Distal Femoral Augment, Right-Medial / Left-Lateral, 68*15, Aluminum	1
237068	Trial Distal Femoral Augment, Right-Medial / Left-Lateral, 72*5, Aluminum	1
237082	Trial Distal Femoral Augment, Right-Medial / Left-Lateral, 72*10, Aluminum	1
237105	Trial Distal Femoral Augment, Right-Medial / Left-Lateral, 72*15, Aluminum	1
237129	Trial Distal Femoral Augment, Right-Medial / Left-Lateral, 76*5, Aluminum	1
237143	Trial Distal Femoral Augment, Right-Medial / Left-Lateral, 76*10, Aluminum	1
237167	Trial Distal Femoral Augment, Right-Medial / Left-Lateral, 76*15, Aluminum	1
237181	Trial Distal Femoral Augment, Right-Medial / Left-Lateral, 80*5, Aluminum	1
237204	Trial Distal Femoral Augment, Right-Medial / Left-Lateral, 80*10, Aluminum	1
237228	Trial Distal Femoral Augment, Right-Medial / Left-Lateral, 80*15, Aluminum	1
236993	Trial Distal Femoral Augment, Right-Lateral / Left-Medial, 68*5, Aluminum	1
237013	Trial Distal Femoral Augment, Right-Lateral / Left-Medial, 68*10, Aluminum	1
237037	Trial Distal Femoral Augment, Right-Lateral / Left-Medial, 68*15, Aluminum	1
237051	Trial Distal Femoral Augment, Right-Lateral / Left-Medial, 72*5, Aluminum	1
237075	Trial Distal Femoral Augment, Right-Lateral / Left-Medial, 72*10, Aluminum	1
237099	Trial Distal Femoral Augment, Right-Lateral / Left-Medial, 72*15, Aluminum	1
237112	Trial Distal Femoral Augment, Right-Lateral / Left-Medial, 76*5, Aluminum	1
237136	Trial Distal Femoral Augment, Right-Lateral / Left-Medial, 76*10, Aluminum	1
237150	Trial Distal Femoral Augment, Right-Lateral / Left-Medial, 76*15, Aluminum	1
237174	Trial Distal Femoral Augment, Right-Lateral / Left-Medial, 80*5, Aluminum	1
237198	Trial Distal Femoral Augment, Right-Lateral / Left-Medial, 80*10, Aluminum	1
237211	Trial Distal Femoral Augment, Right-Lateral / Left-Medial, 80*15, Aluminum	1

TPM REVISION KNEE SYSTEM TIBIA BLOCK RESECTION TRAY

CONTAINER - TRAY - 90125101004



*
Tibial Offset
Correct Positioning
Guide



Item No.	Description	Unit
215073	TPM8 Knee Sys., External Guide, 6, Titanium	1
214946	TPM8 Knee Sys., Femoral Medullary Drill 9 mm, S Steel	1
215066	TPM8 Knee Sys., Medula Guide, 9, Long, Titanium	1
215059	TPM8 Knee Sys., Medula Guide, 9, Short, Titanium	1
350576	TPM Rev. Knee Sys., Intramedullary Tibia Proximal Cutting Guide Holder, S Steel	1
215141	TPM8 Knee Sys., Trial Tibia Insert Holder, Titanium	2
238140	TPM Rev. Knee Sys., Tibial Offset Reamer, 0, S Steel	1
238157	TPM Rev. Knee Sys., Tibial Offset Reamer, 2, S Steel	1
238164	TPM Rev. Knee Sys., Tibial Offset Reamer, 4, S Steel	1
238171	TPM Rev. Knee Sys., Tibial Offset Reamer, 6, S Steel	1
238188	TPM Rev. Knee Sys., Tibia Block Resection, S Steel	1
239581	TPM Rev. Knee Sys., Tibia Block Resection, Left, S Steel	1
239574	TPM Rev. Knee Sys., Tibia Block Resection, Right, S Steel	1



TPM REVISION KNEE SYSTEM TIBIA BLOCK RESECTION TRAY

Item No.	Description	Unit
239567	TPM Rev. Knee Sys., Tibia Augment Resection Block Holder, S Steel	1
350583	TPM Rev. Knee Sys., Tibial Augment Resection, Left, S Steel	1
350590	TPM Rev. Knee Sys., Tibial Augment Resection, Right, S Steel	1
239505	TPM Rev. Knee Sys., Tibia Sled Assembly, 60, S Steel	1
239512	TPM Rev. Knee Sys., Tibia Sled Assembly, 64, S Steel	1
239529	TPM Rev. Knee Sys., Tibia Sled Assembly, 68, S Steel	1
239536	TPM Rev. Knee Sys., Tibia Sled Assembly, 72, S Steel	1
239543	TPM Rev. Knee Sys., Tibia Sled Assembly, 76, S Steel	1
239550	TPM Rev. Knee Sys., Tibia Sled Assembly, 80, S Steel	1
238195	TPM Rev. Knee Sys., Tibial Offset Coin, 0, S Steel	1
238201	TPM Rev. Knee Sys., Tibial Offset Coin, 2, S Steel	1
238218	TPM Rev. Knee Sys., Tibial Offset Coin, 4, S Steel	1
238225	TPM Rev. Knee Sys., Tibial Offset Coin, 6, S Steel	1
238232	TPM Rev. Knee Sys., Tibial Boss Reamer Bushing, 0, S Steel	1
238249	TPM Rev. Knee Sys., Tibial Offset Reamer Bushing, 2, S Steel	1
238256	TPM Rev. Knee Sys., Tibial Offset Reamer Bushing, 4, S Steel	1
238263	TPM Rev. Knee Sys., Tibial Offset Reamer Bushing, 6, S Steel	1
238270	TPM Rev. Knee Sys., Tibial Punch, 60*64, S Steel	1
238287	TPM Rev. Knee Sys., Tibial Punch, 68*72, S Steel	1
238294	TPM Rev. Knee Sys., Tibial Punch, 76*80, S Steel	1
350569	TPM Rev. Knee Sys., Tibia Punch Impactor, S Steel	1
350620	TPM Rev. Knee Sys., Trial Offset Adapter, 0, S Steel	1
350637	TPM Rev. Knee Sys., Trial Offset Adapter, 2, S Steel	1
350644	TPM Rev. Knee Sys., Trial Offset Adapter, 4, S Steel	1
350651	TPM Rev. Knee Sys., Trial Offset Adapter, 6, S Steel	1



Tibial Augment
Resection Guide



TPM REVISION KNEE SYSTEM TRIAL TIBIAL AUGMENTS TRAY

CONTAINER - TRAY - 90125101002



Item No.	Description	Unit
237235	Trial Tibial Metal Back, 60, S Steel	1
237242	Trial Tibial Metal Back, 64, S Steel	1
237259	Trial Tibial Metal Back, 68, S Steel	1
237266	Trial Tibial Metal Back, 72, S Steel	1
237273	Trial Tibial Metal Back, 76, S Steel	1
237280	Trial Tibial Metal Back, 80, S Steel	1
350316	Gap Gauge, 10, Aluminum	1
350323	Gap Gauge, 12, Aluminum	1
350330	Gap Gauge, 14, Aluminum	1
350347	Gap Gauge, 16, Aluminum	1
350354	Gap Gauge, 18, Aluminum	1
350361	Gap Gauge, 20, Aluminum	1
350378	Gap Gauge, 22, Aluminum	1
350385	Gap Gauge, 24, Aluminum	1
237655	Tibial Impactor, Type 1	1
350392	Tibial Impactor, Type 2	1




TPM REVISION KNEE SYSTEM TRIAL TIBIAL AUGMENTS TRAY



Item No.	Description	Unit
237297	Trial Tibial Augment, Right-Lateral / Left-Medial, 60*5, Aluminum	1
237310	Trial Tibial Augment, Right-Lateral / Left-Medial, 60*10, Aluminum	1
237334	Trial Tibial Augment, Right-Lateral / Left-Medial, 60*15, Aluminum	1
237358	Trial Tibial Augment, Right-Lateral / Left-Medial, 64*5, Aluminum	1
237372	Trial Tibial Augment, Right-Lateral / Left-Medial, 64*10, Aluminum	1
237396	Trial Tibial Augment, Right-Lateral / Left-Medial, 64*15, Aluminum	1
237419	Trial Tibial Augment, Right-Lateral / Left-Medial, 68*5, Aluminum	1
237433	Trial Tibial Augment, Right-Lateral / Left-Medial, 68*10, Aluminum	1
237457	Trial Tibial Augment, Right-Lateral / Left-Medial, 68*15, Aluminum	1
237471	Trial Tibial Augment, Right-Lateral / Left-Medial, 72*5, Aluminum	1
237495	Trial Tibial Augment, Right-Lateral / Left-Medial, 72*10, Aluminum	1
237518	Trial Tibial Augment, Right-Lateral / Left-Medial, 72*15, Aluminum	1
237532	Trial Tibial Augment, Right-Lateral / Left-Medial, 76*5, Aluminum	1
237556	Trial Tibial Augment, Right-Lateral / Left-Medial, 76*10, Aluminum	1
237570	Trial Tibial Augment, Right-Lateral / Left-Medial, 76*15, Aluminum	1
237594	Trial Tibial Augment, Right-Lateral / Left-Medial, 80*5, Aluminum	1
237617	Trial Tibial Augment, Right-Lateral / Left-Medial, 80*10, Aluminum	1
237631	Trial Tibial Augment, Right-Lateral / Left-Medial, 80*15, Aluminum	1
237303	Trial Tibial Augment, Right-Medial / Left-Lateral, 60*5, Aluminum	1
237327	Trial Tibial Augment, Right-Medial / Left-Lateral, 60*10, Aluminum	1
237341	Trial Tibial Augment, Right-Medial / Left-Lateral, 60*15, Aluminum	1
237365	Trial Tibial Augment, Right-Medial / Left-Lateral, 64*5, Aluminum	1
237389	Trial Tibial Augment, Right-Medial / Left-Lateral, 64*10, Aluminum	1
237402	Trial Tibial Augment, Right-Medial / Left-Lateral, 64*15, Aluminum	1
237426	Trial Tibial Augment, Right-Medial / Left-Lateral, 68*5, Aluminum	1
237440	Trial Tibial Augment, Right-Medial / Left-Lateral, 68*10, Aluminum	1
237464	Trial Tibial Augment, Right-Medial / Left-Lateral, 68*15, Aluminum	1
237488	Trial Tibial Augment, Right-Medial / Left-Lateral, 72*5, Aluminum	1
237501	Trial Tibial Augment, Right-Medial / Left-Lateral, 72*10, Aluminum	1
237525	Trial Tibial Augment, Right-Medial / Left-Lateral, 72*15, Aluminum	1
237549	Trial Tibial Augment, Right-Medial / Left-Lateral, 76*5, Aluminum	1
237563	Trial Tibial Augment, Right-Medial / Left-Lateral, 76*10, Aluminum	1
237587	Trial Tibial Augment, Right-Medial / Left-Lateral, 76*15, Aluminum	1
237600	Trial Tibial Augment, Right-Medial / Left-Lateral, 80*5, Aluminum	1
237624	Trial Tibial Augment, Right-Medial / Left-Lateral, 80*10, Aluminum	1
237648	Trial Tibial Augment, Right-Medial / Left-Lateral, 80*15, Aluminum	1

TPM REVISION KNEE SYSTEM TRIALS TIBIAL INSERT TRAY

CONTAINER - TRAY - 90125101003




Item No.	Description	Unit
237662	Trial Tibial Insert, Posterior Stabilized, 10*60, Plastic	1
237679	Trial Tibial Insert, Posterior Stabilized, 12*60, Plastic	1
237686	Trial Tibial Insert, Posterior Stabilized, 14*60, Plastic	1
237693	Trial Tibial Insert, Posterior Stabilized, 16*60, Plastic	1
237709	Trial Tibial Insert, Posterior Stabilized, 18*60, Plastic	1
237716	Trial Tibial Insert, Posterior Stabilized, 20*60, Plastic	1
237723	Trial Tibial Insert, Posterior Stabilized, 22*60, Plastic	1
237730	Trial Tibial Insert, Posterior Stabilized, 24*60, Plastic	1
237747	Trial Tibial Insert, Posterior Stabilized, 10*64, Plastic	1
237754	Trial Tibial Insert, Posterior Stabilized, 12*64, Plastic	1
237761	Trial Tibial Insert, Posterior Stabilized, 14*64, Plastic	1
237778	Trial Tibial Insert, Posterior Stabilized, 16*64, Plastic	1
237785	Trial Tibial Insert, Posterior Stabilized, 18*64, Plastic	1
237792	Trial Tibial Insert, Posterior Stabilized, 20*64, Plastic	1
237808	Trial Tibial Insert, Posterior Stabilized, 22*64, Plastic	1
237815	Trial Tibial Insert, Posterior Stabilized, 24*64, Plastic	1
237822	Trial Tibial Insert, Posterior Stabilized, 10*68, Plastic	1
237839	Trial Tibial Insert, Posterior Stabilized, 12*68, Plastic	1
237846	Trial Tibial Insert, Posterior Stabilized, 14*68, Plastic	1
237853	Trial Tibial Insert, Posterior Stabilized, 16*68, Plastic	1
237860	Trial Tibial Insert, Posterior Stabilized, 18*68, Plastic	1
237877	Trial Tibial Insert, Posterior Stabilized, 20*68, Plastic	1
237884	Trial Tibial Insert, Posterior Stabilized, 22*68, Plastic	1
237891	Trial Tibial Insert, Posterior Stabilized, 24*68, Plastic	1



TPM REVISION KNEE SYSTEM TRIALS TIBIAL INSERT TRAY

St. Steel Container
Lid Perforated
530 x 240 x 60 mm



Item No.	Description	Unit
237907	Trial Tibial Insert, Posterior Stabilized, 10*72, Plastic	1
237914	Trial Tibial Insert, Posterior Stabilized, 12*72, Plastic	1
237921	Trial Tibial Insert, Posterior Stabilized, 14*72, Plastic	1
237938	Trial Tibial Insert, Posterior Stabilized, 16*72, Plastic	1
237945	Trial Tibial Insert, Posterior Stabilized, 18*72, Plastic	1
237952	Trial Tibial Insert, Posterior Stabilized, 20*72, Plastic	1
237969	Trial Tibial Insert, Posterior Stabilized, 22*72, Plastic	1
237976	Trial Tibial Insert, Posterior Stabilized, 24*72, Plastic	1
237983	Trial Tibial Insert, Posterior Stabilized, 10*76, Plastic	1
237990	Trial Tibial Insert, Posterior Stabilized, 12*76, Plastic	1
238003	Trial Tibial Insert, Posterior Stabilized, 14*76, Plastic	1
238010	Trial Tibial Insert, Posterior Stabilized, 16*76, Plastic	1
238027	Trial Tibial Insert, Posterior Stabilized, 18*76, Plastic	1
238034	Trial Tibial Insert, Posterior Stabilized, 20*76, Plastic	1
238041	Trial Tibial Insert, Posterior Stabilized, 22*76, Plastic	1
238058	Trial Tibial Insert, Posterior Stabilized, 24*76, Plastic	1
238065	Trial Tibial Insert, Posterior Stabilized, 10*80, Plastic	1
238072	Trial Tibial Insert, Posterior Stabilized, 12*80, Plastic	1
238089	Trial Tibial Insert, Posterior Stabilized, 14*80, Plastic	1
238096	Trial Tibial Insert, Posterior Stabilized, 16*80, Plastic	1
238102	Trial Tibial Insert, Posterior Stabilized, 18*80, Plastic	1
238119	Trial Tibial Insert, Posterior Stabilized, 20*80, Plastic	1
238126	Trial Tibial Insert, Posterior Stabilized, 22*80, Plastic	1
238133	Trial Tibial Insert, Posterior Stabilized, 24*80, Plastic	1



TPM REVISION KNEE SYSTEM PATELLA TRAY

CONTAINER - TRAY - 90127101001




Item No.	Description	Unit
215189	TPM8 Knee Sys., Pin With Head, 25, S Steel	4
215196	TPM8 Knee Sys., Pin With Head, 35, S Steel	4
215202	TPM8 Knee Sys., Pin With Head, 45, S Steel	4
215219	TPM8 Knee Sys., Pin With Head, 55, S Steel	4
215226	TPM8 Knee Sys., Pin With Head, 65, S Steel	4
215165	TPM8 Knee Sys., Pin Straight, 55, S Steel	4
215172	TPM8 Knee Sys., Pin Straight, 65, S Steel	4
359869	Pin with Drill, 65 mm, S Steel	4
359876	Pin with Drill, 80 mm, S Steel	4
215127	TPM8 Knee Sys., Pin Extractor, S Steel	1
215318	TPM8 Knee Sys., Patella Drill, 6, S Steel	1
214953	TPM8 Knee Sys., Bone Rasp	1
215325	TPM8 Knee Sys., Patella Drill Guide, 28-31, S Steel	1
215332	TPM8 Knee Sys., Patella Drill Guide, 34-37, S Steel	1



TPM REVISION KNEE SYSTEM PATELLA TRAY

St. Steel Container
Lid Perforated
530 x 240 x 65 mm



Item No.	Description	Unit
215257	TPM8 Knee Sys., Patella Holder Cutting Guide, S Steel	1
215271	TPM8 Knee Sys., Patella Trial, 28, Plastic	1
215288	TPM8 Knee Sys., Patella Trial, 31, Plastic	1
215295	TPM8 Knee Sys., Patella Trial, 34, Plastic	1
215301	TPM8 Knee Sys., Patella Trial, 37, Plastic	1
215424	TPM8 Knee Sys., Patellar Clamp	1
244103	Pin Adapter, S Steel	1
085065	Drill with Round Shaft, S Steel, 3.2*130	2
200888	Wire Plier, S Steel	1
244059	Trial Offset Adapter Screw For Femur, S Steel	4
244042	Trial Offset Adapter Screw For Tibia, S Steel	4
350149	TPM Revision Knee System, Posterior & Distal Femoral Augment Screw, Titanium	3
349730	TPM Revision Knee System, Tibial Augment Screw, Titanium	3
084068	TPM Revision Knee System, Tibial Augment Screw, Titanium, 5+15	5



TPM REVISION KNEE SYSTEM FEMORAL & TIBIAL TRIALS EXTENSION STEM TRAY

CONTAINER - TRAY - 90125101001

* For surgery double tray is provided




Item No.	Description	Unit
243571	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 40, Ø10, Titanium	2
243588	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 40, Ø11, Titanium	2
243595	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 40, Ø12, Titanium	2
243601	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 40, Ø13, Titanium	2
243618	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 40, Ø14, Titanium	2
243625	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 40, Ø15, Titanium	2
243632	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 40, Ø16, Titanium	2
243649	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 40, Ø17, Titanium	2
243656	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 40, Ø18, Titanium	1
243663	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 40, Ø19, Titanium	1
243670	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 40, Ø20, Titanium	1
243878	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 80, Ø10, Titanium	2
243885	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 80, Ø11, Titanium	2
243892	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 80, Ø12, Titanium	2
243908	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 80, Ø13, Titanium	2
243915	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 80, Ø14, Titanium	2
243922	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 80, Ø15, Titanium	2
243939	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 80, Ø16, Titanium	2
243946	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 80, Ø17, Titanium	2
243953	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 80, Ø18, Titanium	1
243960	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 80, Ø19, Titanium	1
243977	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 80, Ø20, Titanium	1



40 mm 80 mm

TPM REVISION KNEE SYSTEM FEMORAL & TIBIAL TRIALS EXTENSION STEM TRAY

St. Steel Container
Lid Perforated
530 x 240 x 55 mm



Item No.	Description	Unit
244172	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 120, Ø10, Titanium	2
244189	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 120, Ø11, Titanium	2
244196	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 120, Ø12, Titanium	2
244202	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 120, Ø13, Titanium	2
244219	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 120, Ø14, Titanium	2
244226	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 120, Ø15, Titanium	2
244233	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 120, Ø16, Titanium	2
244240	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 120, Ø17, Titanium	2
244257	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 120, Ø18, Titanium	1
244264	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 120, Ø19, Titanium	1
244271	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 120, Ø20, Titanium	1
244479	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 160, Ø10, Titanium	1
244486	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 160, Ø11, Titanium	1
244493	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 160, Ø12, Titanium	1
244509	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 160, Ø13, Titanium	1
244516	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 160, Ø14, Titanium	1
244523	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 160, Ø15, Titanium	1
244530	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 160, Ø16, Titanium	1
244547	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 160, Ø17, Titanium	1
244554	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 160, Ø18, Titanium	1
244561	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 160, Ø19, Titanium	1
244578	TPM Rev. Knee Sys., Trial Tibial & Femoral Stem Extension, 160, Ø20, Titanium	1



TPM REVISION KNEE SYSTEM FEMORAL & TIBIAL REAMERS TRAY

CONTAINER - TRAY - 90128101001



Item No.	Description	Unit
351153	Femoral & Tibial Reamer Ø10 mm	1
351160	Femoral & Tibial Reamer Ø11 mm	1
351177	Femoral & Tibial Reamer Ø12 mm	1
351184	Femoral & Tibial Reamer Ø13 mm	1
351191	Femoral & Tibial Reamer Ø14 mm	1
351207	Femoral & Tibial Reamer Ø15 mm	1
351214	Femoral & Tibial Reamer Ø16 mm	1
351221	Femoral & Tibial Reamer Ø17 mm	1
351238	Femoral & Tibial Reamer Ø18 mm	1
351245	Femoral & Tibial Reamer Ø19 mm	1
351252	Femoral & Tibial Reamer Ø20 mm	1



TPM REVISION KNEE SYSTEM FEMORAL & TIBIAL REAMERS TRAY

Item No.	Description	Unit
244110	T - Handle Type1, Short, S Steel	1
351269	Femoral & Tibial Reamer Adaptor	1
214991	Feeler Blade Long	1
214984	Feeler Blade Small	1
351078	Screw Driver A3 mm	2
351276	Screw Driver A5 mm	1
351283	Universal Driver A3 mm	1
308409	Pin-Punc Impactor - Extractor	1



Regain
your
mobility



RIISING SUN OF YOUR LIFE

You can feel the
experience, confidence
and quality in field
of arthroplasty,
traumatology and
spine.



Company Profile:

Tipmed is a prominent medical products manufacturer, specializing in implants and instruments in the field of orthopedics, traumatology and neurochirurgie. It was established in 2001 in Izmir by the industry veteran Sedat Onel to provide a path-breaking solution to the most demanding medical device manufacturing challenges.

Products:

Orthopedic instruments and implant systems such as Total Hip Replacement, Knee and Shoulder Replacement, Spine Systems, External Fixators and Trauma Systems like Locking Plating Systems, Intramedullary Nailing Systems, Screw Systems.

www.tipmed.com

*Regain
your
mobility*



TIPMED Medical Device Manufacturing Ltd. Co.

Kemalpaşa Mahallesi 7401/1 Sokak No:16
Bornova / İzmir - TÜRKİYE

T : +90 232 479 33 22

M : +90 543 479 33 22

F : +90 232 479 38 25

info@tipmed.com

marketing@tipmed.com

www.tipmed.com

