

*Regain
your
mobility*

S
S
R
M
T
P



TPM
Reverse
SHOULDER
SYSTEM

STANDARD HUMERAL STEM



Titanium Alloyed, Cemented & Uncemented, Porous Coated

Length : 122 mm
Diameter: 7 mm up to 20 mm with 1 mm increment

* 16-17-18-19-20mm are the optional sizes

MINI HUMERAL STEM



Titanium Alloyed, Cemented & Uncemented, Porous Coated

Length : 83 mm
Diameter: 7 mm up to 20 mm with 1 mm increment

* 16-17-18-19-20mm are the optional sizes

MICRO HUMERAL STEM



Titanium Alloyed, Cemented & Uncemented, Porous Coated

Length : 55 mm
Diameter: 7 mm up to 20 mm with 1 mm increment

Optional for all size

FRACTURE HUMERAL STEM



CoCrMo Alloyed, Cemented & Uncemented, Plasma Spray Coated

Length : 122 mm
Diameter: 6 mm up to 14 mm with 2 mm increment

REVISION HUMERAL STEM



Titanium Alloyed, Cemented & Uncemented, Porous Coated

Length : 160 mm, 194 mm
Diameter: 6 mm up to 12 mm with 2 mm increment

HUMERAL TRAY

Titanium Alloyed
Thickness : 0 - 5 - 10 mm, Size : 44 mm

MINI HUMERAL TRAY

Titanium Alloyed
Thickness : 0 - 5 - 10 mm, Size : 40 mm



HUMERAL CUP

UHMWPE Alloyed
Thickness : 0 - 3 mm, Size : 36 - 38 - 41 mm

MINI HUMERAL CUP

UHMWPE Alloyed
Thickness : 0 - 3 mm, Size : 36 - 38 - 41 mm



TPM Reverse SHOULDER SYSTEM

HUMERAL HEAD (ECCENTRIC)

Titanium & CoCrMo Alloyed
Size : 40 mm / 15 - 20 - 22 mm
42 - 46 mm / 18 - 21 - 24 mm
50 mm / 21 - 24 - 27 mm

Optional Size : 38
44 - 48
54 - 58



BASEPLATE (ECCENTRIC) & STANDARD BASEPLATE

Titanium Alloyed
Uncemented, Porous + H.A. Coated
Size : Ø25 mm



BASEPLATE WITH EXTENSION (ECCENTRIC)

Titanium Alloyed
Uncemented, Porous + H.A. Coated
Size : Ø25 mm



STANDARD BASEPLATE WITH EXTENSION

Titanium Alloyed
Uncemented, Porous + H.A. Coated
Size : Ø25 mm



4,75 MM SHOULDER SCREW

Titanium Alloyed
Locking & Non-Locking
Length : 15- 20 - 25 - 30 - 35 - 40 - 45 - 50 mm



6,5 MM SHOULDER SCREW

Titanium Alloyed
Length : 20 - 25 - 30 - 35 - 40 - 45 - 50 55 - 60 - 65 mm



STANDARD HUMERAL HEAD

Titanium & CoCrMo Alloyed
Size : 40 mm / 15 - 20 - 22 mm
42 - 46 mm / 18 - 21 - 24 mm
50 mm / 21 - 24 - 27 mm

Optional Size : 38
44 - 48
54 - 58



GLENOSPHERE (ECCENTRIC)

Titanium & CoCrMo Alloyed
Offset : 0 - 3 - 6 mm
Size : 36 - 38 - 41 mm



STANDARD GLENOSPHERE

Titanium & CoCrMo Alloyed
Offset : 0 - 3 - 6 mm
Size : 36 - 38 - 41 mm

