



*Regain  
your  
mobility*



STRAIGHT STEM  
**CRL TYPE HIP**  
PROSTHESIS

[www.tipmed.com](http://www.tipmed.com)

## STRAIGHT STEM, CRL TYPE HIP PROSTHESIS

Titanium Alloyed (Uncemented)  
St. Steel Alloyed (Cemented)

For uncemented prosthesis,  
two coating options;

- HA Coated
- Dual Coated (Titanium+HA Coated)

### TITANIUM + HA COATED




### CEMENTED



Stem Size	Stem Length
No.0	120 mm
No.1	125 mm
No.2	130 mm
No.3	135 mm
No.4	140 mm
No.5	145 mm
No.6	150 mm
No.7	155 mm
No.8	160 mm
No.9	165 mm
No.10	170 mm

### INSTRUMENT SET TRAY - 1



### INSTRUMENT SET TRAY - 2



STRAIGHT STEM  
**CRL TYPE HIP**  
PROSTHESIS