



*Regain  
your  
mobility*



**T I P M E D**  
**INTERTROCHANTERIC**  
ANTEGRADE NAILS

**TPM PFN SYSTEM,  
INTERTROCHANTERIC  
ANTEGRADE NAILS**

Titanium Alloyed

Length : 200 mm - 240 mm  
280 mm - 320 mm  
are the optional sizes

Diameter:  
10 mm - 11 mm  
12 mm - 13 mm



**TPM PFN SYSTEM,  
INTEGRATED  
INTERLOCKING  
SCREW**

Titanium Alloyed

Length :  
70 - 75 - 80 - 85 - 90 - 95 - 100 -  
105 - 110 - 115 - 120 mm

Diameter: 10,5 mm



**TPM PFN SYSTEM,  
INTEGRATED  
INTERLOCKING  
SCREW**

Titanium Alloyed

Length :  
65 - 70 - 75 - 80 - 85 - 90 - 95 -  
100 - 105 - 110 - 115 mm

Diameter: 6,5 mm



**5MM HEX.PULLER  
LOCKING SCREW  
SELF TAPPING FOR  
INTRAMEDULLARY  
NAIL, 4.8 mm**

Titanium Alloyed

Length :  
30 - 32 - 34 -  
36 - 38 - 40 -  
45 - 50 mm

Diameter: 4,8 mm



**TPM PFN SYSTEM,  
INTERTROCHANTERIC  
ANTEGRADE NAILS END CUP**

Titanium Alloyed  
0 - 5 - 10 - 15 mm



**TPM PFN SYSTEM,  
INTERTROCHANTERIC  
ANTEGRADE NAILS SET SCREW**

Titanium Alloyed




**TIPMED  
INTERTROCHANTERIC  
ANTEGRADE  
NAILS**



**INSTRUMENT  
SET TRAY - 1**



**INSTRUMENT  
SET TRAY - 2**



**INSTRUMENT  
SET TRAY - 3**



**INSTRUMENT  
SET TRAY - 4**



**INSTRUMENT  
SET TRAY - 5**

