



INTRAMEDULLARY NAILS

*Regain
your
mobility*

D

INTRAMEDULLARY NAILS IMPLANTS & INSTRUMENTS PAGE

CONTENTS

PAGE D4

ELASTIC NAIL SET



PAGE D6

TPM PFN SYSTEM,
INTERTROCHANTERIC
ANTEGRADE NAIL SET



PAGE D12

130°
PROXIMAL
FEMORAL
INTERLOCKING
NAIL SET
(TPM-PFN)



PAGE D20

CANNULATED
PFN PROXIMAL
FEMORAL
INTERLOCKING
NAIL SET



PAGE D26

MULTIAXIAL
FEMORAL
INTERLOCKING
NAIL SET



PAGE D34

FEMORAL
DYNAMIC
INTERLOCKING
NAIL SET



PAGE D40

FEMORAL
INTERLOCKING
NAIL SET



PAGE D46

SUPRACONDYLAR
INTERLOCKING
NAIL SET



PAGE D50

TIBIAL DYNAMIC
INTERLOCKING
NAIL SET



PAGE D56

TIBIAL
INTERLOCKING
NAIL SET



PAGE D61

HUMERAL
INTERLOCKING
NAIL SET




PAGE D66

FIBULA
INTERLOCKING
NAIL SET



ELASTIC NAIL SET
TITANIUM 90503022001



TRAY NO 1



290410075500169
Aluminium Container
500 x 169 x 75 mm




ORDER INFORMATION 8698791 ...


Length	Elastic Nail Plug		
	non-Sterile Titanium	Unit	Sterile Titanium
Ø2 - 2,5 mm	053361	4	498940
Ø3 - 3,5 - 4 mm	053583	4	498957

ORDER INFORMATION 8698791 ...

Specification	Elastic Nail		
	non-Sterile Titanium	Unit	Sterile Titanium
2.0 x 300 mm	123262	2	498964
2.0 x 450 mm	123279	2	498971
2.5 x 300 mm	123286	2	498988
2.5 x 450 mm	123293	2	498995
3.0 x 300 mm	123309	2	499466
3.0 x 450 mm	123316	2	499473
3.5 x 300 mm	123323	2	499480
3.5 x 450 mm	123330	2	499497
4.0 x 300 mm	123347	2	499503
4.0 x 450 mm	123354	2	499510




ELASTIC NAIL SET



Item No.	Description	Unit
210290	Elastic Pin Impactor	1
210238	Elastic Pin Trocar Straight	1
210245	Elastic Pin Trocar Curved	1
239956	Elastic Pin Plug Inserter	1
209256	Wire Cutter, S.Steel, 4mm	1
085065	Drill with Round Shaft, S.Steel, 3.2*130	1
085089	Drill with Round Shaft, S.Steel, 3.5*130	1
085102	Drill with Round Shaft, S.Steel, 4.0*130	1
200888	Wire Plier, S.Steel	1
210184	Key for Universal Chuck, S.Steel	1
457893	Elastic Pin T - Handle For Drill, S.Steel	1



**TPM PFN SYSTEM, INTERTROCHANTERIC ANTEGRADE NAIL SET
TITANIUM 90530082001**



TRAY NO 1



290410260580280
Aluminium Container
580 x 280 x 260 mm

ORDER INFORMATION 8698791 ...

**TPM PFN System,
Intertrochanteric Antegrade Nails**

Specification	non-Sterile Titanium	Unit	Sterile Titanium
10 x 200 mm	458630	1	506980
10 x 240 mm	458647	1	506997
10 x 280 mm *	458654	-	507130
10 x 320 mm *	458661	-	507147
11 x 200 mm	458678	1	507154
11 x 240 mm	458685	1	507161
11 x 280 mm *	458692	-	507178
11 x 320 mm *	458708	-	507185
12 x 200 mm	458715	1	507192
12 x 240 mm	458722	1	507338
12 x 280 mm *	458739	-	507345
12 x 320 mm *	458746	-	507352
13 x 200 mm	494218	1	507369
13 x 240 mm	494225	1	507376
13 x 280 mm *	494232	-	507383
13 x 320 mm *	494249	-	507390



TPM PFN SYSTEM, INTERTROCHANTERIC ANTEGRADE NAIL SET

ORDER INFORMATION 8698791 ...

TPM PFN System, Integrated Interlocking Screw				TPM PFN System, Integrated Interlocking Screw			
Length	non-Sterile Titanium	Unit	Sterile Titanium	Length	non-Sterile Titanium	Unit	Sterile Titanium
10.5 x 70 mm	458753	1	516583	6.5 x 65 mm	458869	1	516156
10.5 x 75 mm	458760	1	516590	6.5 x 70 mm	458876	1	516163
10.5 x 80 mm	458777	1	516606	6.5 x 75 mm	458883	1	516170
10.5 x 85 mm	458784	1	516613	6.5 x 80 mm	458890	1	516187
10.5 x 90 mm	458791	1	516620	6.5 x 85 mm	458906	1	516194
10.5 x 95 mm	458807	1	516637	6.5 x 90 mm	458913	1	516361
10.5 x 100 mm	458814	1	516644	6.5 x 95 mm	458920	1	516378
10.5 x 105 mm	458821	1	516811	6.5 x 100 mm	458937	1	516385
10.5 x 110 mm	458838	1	516828	6.5 x 105 mm	458944	1	516392
10.5 x 115 mm	458845	1	516835	6.5 x 110 mm	458951	1	516408
10.5 x 120 mm	458852	1	516842	6.5 x 115 mm	458968	1	516415

ORDER INFORMATION 8698791 ...

5mm Hex.Puller Locking Screw Self Tapping For Intramedullary Nail							
Length	non-Sterile Titanium	Unit	Sterile Titanium	Length	non-Sterile Titanium	Unit	Sterile Titanium
4,8 x 30 mm	489269	2	528999	4,8 x 55 mm *	489344	-	529071
4,8 x 32 mm	489276	2	529002	4,8 x 60 mm *	489351	-	529088
4,8 x 34 mm	489283	2	529019	4,8 x 65 mm *	489368	-	529095
4,8 x 36 mm	489290	2	529026	4,8 x 70 mm *	489375	-	529101
4,8 x 38 mm	489306	2	529033	4,8 x 75 mm *	489382	-	529118
4,8 x 40 mm	489313	2	529040	4,8 x 80 mm *	489399	-	529125
4,8 x 45 mm	489320	2	529057	4,8 x 85 mm *	489405	-	529132
4,8 x 50 mm	489337	2	529064	4,8 x 90 mm *	489412	-	529149

ORDER INFORMATION 8698791 ...

TPM PFN System, Intertrochanteric Antegrade Nails End Cap				TPM PFN System, Intertrochanteric Antegrade Nails Set Screw		
Size	non-Sterile Titanium	Unit	Sterile Titanium	non-Sterile Titanium	Unit	Sterile Titanium
0 mm	458975	1	507529	459019	2	507512
5 mm	458982	1	507536			
10 mm	458999	1	507543			
15 mm	459002	1	507550			


TPM PFN SYSTEM, INTERTROCHANTERIC ANTEGRADE NAIL SET



Item No.	Description	Unit
459040	TPM PFN System, Lag T Driver, S.Steel	1
459071	TPM PFN System, Hex Driver, S.Steel, 7	1
459309	TPM PFN System, Drill Bit, S.Steel Ø6.5 / Ø7.6	1
459262	TPM PFN System, Lag Drill, S.Steel, Ø10.5	1
459095	TPM PFN System, Stopper, S.Steel	1
473985	TPM PFN System, Canullated Cortex Reamer For Lag Screw, S.Steel	1
459064	TPM PFN System, Hex Driver, S.Steel, 5	1
459088	TPM PFN System, Hex Driver Connector Screw, S.Steel, 5	1
459408	TPM PFN System, Lag Screw Holder, S.Steel	2
204510	Allen Wrench, S.Steel, 3	1



TPM PFN SYSTEM, INTERTROCHANTERIC ANTEGRADE NAIL SET



Item No.	Description	Unit
459453	TPM PFN System, Sliding Hammer, S.Steel	1
459255	TPM PFN System, Flexible Driver, S.Steel, 5	1
459194	TPM PFN System, Entry Reamer, S.Steel	1
459323	TPM PFN System, Wrench, S.Steel, 11	1
459460	TPM PFN System, Nail Extractor Rod, S.Steel	1
459415	TPM PFN System, Thread Pin, S.Steel, Ø2,8*400	5
485612	Gripper For Guide Pin	1
457916	T - Handle For Drill, S.Steel	1
459101	TPM PFN System, Tissue Protector, S.Steel	1
210184	Key for Universal Chuck, S.Steel	1



TPM PFN SYSTEM, INTERTROCHANTERIC ANTEGRADE NAIL SET


Item No.	Description	Unit
457329	Drill Guide, S.Steel, 4.0*189	2
457336	Screw Guide, S.Steel, 8.0*169	2
459286	TPM PFN System, Drill Stop, S.Steel, Ø6.5 / Ø7.6	1
459293	TPM PFN System, Drill Stop, S.Steel, Long, Ø6.5 / Ø7.6	1
459446	TPM PFN System, Lag Depth Gauge, Ø3*400, S.Steel	1
200505	Flexible Reamer, S.Steel, 10	1
200512	Flexible Reamer, S.Steel, 11	1
200529	Flexible Reamer, S.Steel, 12	1
459279	TPM PFN System, Drill Stop, S.Steel, Ø4*300	1
085317	Drill with Round Shaft, S.Steel, 4.0*300	3
200635	Adaptor For Flexible Reamer, S.Steel	1
459163	TPM PFN System, Pin Sleeve, S.Steel	1
200628	Depth Gauge For Distal Screw, S.Steel, 120	1
200437	Trocar, Ø8*300	1




TPM PFN SYSTEM, INTERTROCHANTERIC ANTEGRADE NAIL SET



Item No.	Description	Unit
459026	TPM PFN System, Compression Nut, S.Steel	1
459132	TPM PFN System, Proximal Guide Sleeve, S.Steel	1
459149	TPM PFN System, Guide Handle, S.Steel	1
459170	TPM PFN System, Distal Dynamic Guide Arm, (200-240) Aluminum	1
459187	TPM PFN System, Proximal Guide Arm / Distal Static, PPSU	1
200277	Guide Wire, S.Steel, 2.5*690	2
459330	TPM PFN System, Handle-Distal Arm Connector Screw, S.Steel	2
459347	TPM PFN System, Handle Nail Connector Screw, S.Steel	2
213567	Cannulated AWL	1
459422	TPM PFN System, Impactor Connector, S.Steel	1



**130° PROXIMAL FEMORAL INTERLOCKING NAIL SET (TPM-PFN)
TITANIUM 90530072001**



TRAY NO 1



290410260580280
Aluminium Container
580 x 280 x 260 mm



ORDER INFORMATION 8698791 ...

Length	TPM-PFN End Cap		
	non-Sterile Titanium	Unit	Sterile Titanium
0 mm	526636	4	500223
5 mm	052784	1	500230
10 mm	052791	1	500247
15 mm	052807	1	500254

130° PROXIMAL FEMORAL INTERLOCKING NAIL SET (TPM-PFN)



ORDER INFORMATION 8698791 ...

130° TPM-PFN Proximal Femoral Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
9 x 200 mm *	525684	-	499527
9 x 240 mm *	525691	-	499534
9 x 280 mm *	525707	-	499541
9 x 320 mm *	525714	-	499558

130° TPM-PFN Proximal Femoral Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
10 x 200 mm	525745	1	499565
10 x 240 mm	525752	1	499572
10 x 280 mm	525769	1	499589
10 x 320 mm	525776	1	499596

130° TPM-PFN Proximal Femoral Interlocking Nail


Specification	non-Sterile Titanium	Unit	Sterile Titanium
11 x 200 mm	525806	1	499602
11 x 240 mm	525813	1	499619
11 x 280 mm	525820	1	499626
11 x 320 mm	525837	1	499633

130° TPM-PFN Proximal Femoral Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
12 x 200 mm	525868	1	499640
12 x 240 mm	525875	1	499657
12 x 280 mm	525882	1	499664
12 x 320 mm	525899	1	499671

* Optional size

130° PROXIMAL FEMORAL INTERLOCKING NAIL SET (TPM-PFN)



ORDER INFORMATION 8698791 ...

TPM-PFN Cannulated Compression Lag Screw Standard Threaded

Length	non-Sterile Titanium	Unit	Sterile Titanium
11 x 70 mm	526162	1	512387
11 x 75 mm	526179	1	512394
11 x 80 mm	526186	1	512554
11 x 85 mm	526193	1	512561
11 x 90 mm	526209	1	512578
11 x 95 mm	526216	1	512585
11 x 100 mm	526223	1	512592
11 x 105 mm	526230	1	512608
11 x 110 mm	526247	1	512615
11 x 115 mm	526254	1	512622
11 x 120 mm	526261	1	512639



ORDER INFORMATION 8698791 ...

TPM-PFN Cannulated Compression Lag Screw With Blade

Length	non-Sterile Titanium	Unit	Sterile Titanium
11 x 70 mm	526650	1	512646
11 x 75 mm	526667	1	512653
11 x 80 mm	526674	1	512660
11 x 85 mm	526681	1	512677
11 x 90 mm	526698	1	512684
11 x 95 mm	526704	1	512691
11 x 100 mm	526711	1	512851
11 x 105 mm	526728	1	512868
11 x 110 mm	526735	1	512875
11 x 115 mm	526742	1	512882
11 x 120 mm	526759	1	512899



130° PROXIMAL FEMORAL INTERLOCKING NAIL SET (TPM-PFN)

ORDER INFORMATION 8698791 ...

Locking Screw Self Tapping For Femoral Intramedullary Nail

Length	non-Sterile Titanium	Unit	Sterile Titanium
4,8x30 mm	069331	2	513223
4,8x35 mm	069348	2	513230
4,8x40 mm	069355	2	513247
4,8x45 mm	069362	2	513254
4,8x50 mm	069379	2	513261
4,8x55 mm	069386	2	513278
4,8x60 mm	069393	2	513285
4,8x65 mm	069409	2	513292
4,8x70 mm *	069416	-	513315
4,8x75 mm *	069423	-	513322
4,8x80 mm *	069430	-	513339
4,8x85 mm *	084129	-	513346
4,8x90 mm *	084136	-	513353



Item No.	Description	Unit
----------	-------------	------

457077	TPM-PFN Proximal & Distal Hole Device, Aluminum	1
--------	---	---



527626	TPM-PFN Nail Lock, S.Steel	1
--------	----------------------------	---



200239	Tissue Protector, S.Steel	1
--------	---------------------------	---




527602	TPM-PFN Directing Device, Right, Aluminum	1
--------	---	---



527619	TPM-PFN Directing Device, Left, Aluminum	1
--------	--	---


* Optional size

130° PROXIMAL FEMORAL INTERLOCKING NAIL SET (TPM-PFN)




TRAY NO 3

Item No.	Description	Unit
216544	Sliding Hammer Bar, S.Steel	1
458241	Slotted Hammer, 14	1
527534	Direction Indicator Device	1
527589	TPM-PFN Screw Driver & Holder for Spiral Blade Screw	1
200420	Pointed - Cortical Awl, S.Steel	1
216551	Supine Driver For Multiaxial Femur and TPM-PFN, S.Steel	1
210184	Key for Universal Chuck, S.Steel	1
457916	T - Handle For Drill, S.Steel	1
200628	Depth Gauge For Distal Screw, S.Steel, 120	1
527510	Wire Guide, S.Steel, 2.5*250	1
527503	Cannulated Screw Guide, S.Steel, 11	1
457329	Drill Guide, S.Steel, 4.0*189	2
457336	Screw Guide, S.Steel, 8.0*169	2
458128	Contra Wrench for Spiral Driver, S.Steel	2




130° PROXIMAL FEMORAL INTERLOCKING NAIL SET (TPM-PFN)




Item No.	Description	Unit
210375	Guide Wire, S.Steel, 2.5*423	2
085317	Drill with Round Shaft, S.Steel, 4.0*300	2
527572	TPM-PFN Screw Driver & Holder for Lag Screw, S.Steel	1
527640	Lag Screw Remover, S.Steel	1
210405	Tapped Reamer, S.Steel, 16	1
527558	Cannulated Screw Driver, (POMC), 6	1
527633	TPM-PFN Cannulated Threaded Drill with Scale for Spiral Blade Screw, S.Steel	1
527527	TPM-PFN Cannulated Threaded Drill with Scale for Lag Screw, S.Steel	1
308386	T - Holder for Set Screw, S.Steel, 3.5	1
210429	TPM-PFN Depth Gauge for Proximal Screw, S.Steel, 2.5	1
200437	Trocar, Ø8*300	1




130° PROXIMAL FEMORAL INTERLOCKING (TPM-PFN) NAIL SET



Item No.	Description	Unit
200505	Flexible Reamer, S.Steel, 10	1
200512	Flexible Reamer, S.Steel, 11	1
200529	Flexible Reamer, S.Steel, 12	1
457909	T - Handle For Flexible Reamer, S.Steel	1
200635	Adaptor For Flexible Reamer, S.Steel	1
130031	T - Reamer, S.Steel, 10	1
130048	T - Reamer, S.Steel, 11	1
130055	T - Reamer, S.Steel, 12	1
200390	Wrench, S.Steel, 16 - 17	2
239611	Cannulated First Reamer, S.Steel, 10	1
456995	Articulated Screwdriver, S.Steel, 8	1
200277	Guide Wire, S.Steel, 2.5*690	2



130° TPM-PFN PROXIMAL FEMORAL INTERLOCKING NAIL (RIGHT-LEFT)



TRAY *




ORDER INFORMATION 8698791 ...

130° TPM-PFN Proximal Femoral Interlocking Nail

Specification	Side	non-Sterile Titanium	Sterile Titanium	Side	non-Sterile Titanium	Sterile Titanium
9 x 340 mm *	Right	528234	499848	Left	528654	499688
9 x 360 mm *	Right	525721	499855	Left	528661	499695
9 x 380 mm *	Right	528241	499862	Left	528678	499701
9 x 400 mm *	Right	525738	499879	Left	528685	499718
10 x 340 mm *	Right	528173	499886	Left	528692	499725
10 x 360 mm *	Right	525783	499893	Left	528708	499732
10 x 380 mm *	Right	528180	499909	Left	528715	499749
10 x 400 mm *	Right	525790	499916	Left	528722	499756
11 x 340 mm *	Right	528197	499923	Left	528739	499763
11 x 360 mm *	Right	525844	499930	Left	528746	499770
11 x 380 mm *	Right	528203	499947	Left	528753	499787
11 x 400 mm *	Right	525851	499954	Left	528760	499794
12 x 340 mm *	Right	528210	499961	Left	528777	499800
12 x 360 mm *	Right	525905	499978	Left	528784	499817
12 x 380 mm *	Right	528227	499985	Left	528791	499824
12 x 400 mm *	Right	525912	499992	Left	528807	499831

* Optional for all size

CANNULATED PFN PROXIMAL FEMORAL INTERLOCKING NAIL SET
TITANIUM 90530062001



TRAY NO 1



290410150580280
 Aluminium Container
 580 x 280 x 150 mm

ORDER INFORMATION 8698791 ...



Nail End Cap For PFN

Length	non-Sterile Titanium	Unit	Sterile Titanium
0 mm	070313	4	508489
5 mm	052753 *	-	508496
10 mm	052760 *	-	508625
15 mm	052777 *	-	508632

ORDER INFORMATION 8698791 ...



Set Screw For PFN

non-Sterile Titanium	Unit	Sterile Titanium
070290	4	513834

CANNULATED PFN PROXIMAL FEMORAL INTERLOCKING NAIL SET

Ø10 mm
Yellow



Ø11 mm
Blue



Ø12 mm
Purple



ORDER INFORMATION 8698791 ...

Canullated PFN Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
10 x 200 mm	499008	1	508083
10 x 240 mm	499015	1	508090
10 x 280 mm	499022	1	508243
10 x 320 mm	499039	1	508250


Canullated PFN Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
11 x 200 mm	499046	1	508267
11 x 240 mm	499053	1	508274
11 x 280 mm	499060	1	508281
11 x 320 mm	499077	1	508298

Canullated PFN Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
12 x 200 mm	499084	1	508441
12 x 240 mm	499091	1	508458
12 x 280 mm	499107	1	508465
12 x 320 mm	499114	1	508472

CANNULATED PFN PROXIMAL FEMORAL INTERLOCKING NAIL SET



ORDER INFORMATION 8698791 ...

Cannulated Locking Screw For Canullated PFN				Locking Screw Self Tapping For Femoral Intramedullary Nail			
Length	non-Sterile Titanium	Unit	Sterile Titanium	Length	non-Sterile Titanium	Unit	Sterile Titanium
9 x 70 mm	499121	2	513483	4,8 x 30 mm	069331	2	513223
9 x 75 mm	499138	2	513490	4,8 x 35 mm	069348	2	513230
9 x 80 mm	499145	2	513742	4,8 x 40 mm	069355	2	513247
9 x 85 mm	499152	2	513759	4,8 x 45 mm	069362	2	513254
9 x 90 mm	499169	2	513766	4,8 x 50 mm	069379	2	513261
9 x 95 mm	499176	2	513773	4,8 x 55 mm	069386	2	513278
9 x 100 mm	499183	2	513780	4,8 x 60 mm	069393	2	513285
9 x 105 mm	499190	1	513797	4,8 x 65 mm	069409	2	513292
9 x 110 mm	499206	1	513803	4,8 x 70 mm	069416	2	513315
9 x 115 mm	499213	1	513810	4,8 x 75 mm	069423	2	513322
9 x 120 mm	499220	1	513827	4,8 x 80 mm	069430	2	513339

Item No.	Description	Unit
----------	-------------	------


210313	Cannulated PFN Main Device, Aluminum	1
--------	--------------------------------------	---

213550	Cannulated PFN Nail Lock, S.Steel	1
--------	-----------------------------------	---

200642	Elevator Adjustment T Wrench, S.Steel	1
--------	---------------------------------------	---



CANNULATED PFN PROXIMAL FEMORAL INTERLOCKING NAIL SET




Item No.	Description	Unit
200604	Sliding Hammer	1
200420	Pointed - Cortical Awl, S.Steel	1
210184	Key for Universal Chuck, S.Steel	1
457916	T - Handle For Drill, S.Steel	1
200611	Supine Driver, S.Steel	1
210351	Wire Guide, S.Steel, 2.5*235	2
210368	Screw Guide, S.Steel, 9.0	2
457329	Drill Guide, S.Steel, 4.0*189	2
457336	Screw Guide, S.Steel, 8.0*169	2
200628	Depth Gauge For Distal Screw, S.Steel, 120	1
085317	Drill with Round Shaft, S.Steel, 4.0*300	2
457350	Guide Wire, S.Steel, 2.5*453	2




CANNULATED PFN PROXIMAL FEMORAL INTERLOCKING NAIL SET



Item No.	Description	Unit
200505	Flexible Reamer, S.Steel, 10	1
200512	Flexible Reamer, S.Steel, 11	1
200529	Flexible Reamer, S.Steel, 12	1
457909	T - Handle For Flexible Reamer, S.Steel	1
200635	Adaptor For Flexible Reamer, S.Steel	1
130031	T - Reamer, S.Steel, 10	1
130048	T - Reamer, S.Steel, 11	1
130055	T - Reamer, S.Steel, 12	1
239611	Cannulated First Reamer, S.Steel, 10	1
200277	Guide Wire, S.Steel, 2.5*690	1



CANNULATED PFN PROXIMAL FEMORAL INTERLOCKING NAIL SET




Item No.	Description	Unit
200390	Wrench, S.Steel, 16 - 17	2
200437	Trocar, Ø8*300	1
309420	Cannulated PFN Screw Holder For Cannulated Locking Screw, S.Steel, 9	1
458180	Depth Gauge For Proximal Screw, S.Steel, 120	1
210412	T - Trocar, S.Steel, 9.0	1
210405	Tapped Reamer, S.Steel, 16	1
308386	T - Holder for Set Screw, S.Steel, 3.5	1
200673	T - Wrench For End Cap, S.Steel, 5.0	1
210436	Cannulated PFN Depth Gauge For Proximal Screw, S.Steel, 2.5	1
200666	T - Wrench For Locking Screw With Compression Screw, S.Steel, 3.5	1
308379	T - Wrench For End Cap With Threaded, S.Steel, 5.0 (Femoral)	1
210382	Cannulated Drill, S.Steel, 9.0	1
200239	Tissue Protector, S.Steel	1
308393	Puller Shaft For Femoral, Tibial, Humerus T - Wrench, S.Steel	2



MULTIAXIAL FEMORAL INTERLOCKING NAIL SET

TITANIUM 90550012001



290410260580280
Aluminium Container
for Instrument
580 x 280 x 260 mm




290410150580280
Aluminium Container
for Implant
580 x 280 x 150 mm

ORDER INFORMATION 8698791 ...

Multiaxial Femoral Locking Intramedullary Nail

Specification	Side	non-Sterile Titanium	Unit	Sterile Titanium	Side	non-Sterile Titanium	Unit	Sterile Titanium
9 x 280 mm	Right	400004	1	500261	Left	400554	1	501633
9 x 300 mm	Right	400011	1	500278	Left	400561	1	501640
9 x 320 mm	Right	400028	1	500285	Left	400578	1	501657
9 x 340 mm	Right	400035	1	500292	Left	400585	1	501664
9 x 360 mm	Right	400042	1	500414	Left	400592	1	501671
9 x 380 mm	Right	400059	1	500421	Left	400608	1	501688
9 x 400 mm	Right	400066	1	500438	Left	400615	1	501695


MULTIAXIAL FEMORAL INTERLOCKING NAIL SET



ORDER INFORMATION 8698791 ...

Multiaxial Femoral Locking Intramedullary Nail


Specification	Side	non-Sterile Titanium	Unit	Sterile Titanium	Side	non-Sterile Titanium	Unit	Sterile Titanium
10 x 300 mm	Right	400073	1	500445	Left	400622	1	501817
10 x 320 mm	Right	400080	1	500452	Left	400639	1	501824
10 x 340 mm	Right	400097	1	500469	Left	400646	1	501831
10 x 360 mm	Right	400103	1	500476	Left	400653	1	501848
10 x 380 mm	Right	400110	1	500483	Left	400660	1	501855
10 x 400 mm	Right	400127	1	500490	Left	400677	1	501862
10 x 420 mm	Right	400134	1	500612	Left	400684	1	501879
10 x 440 mm	Right	400141	1	500629	Left	400691	1	501886




ORDER INFORMATION 8698791 ...

Multiaxial Femoral Locking Intramedullary Nail

Specification	Side	non-Sterile Titanium	Unit	Sterile Titanium	Side	non-Sterile Titanium	Unit	Sterile Titanium
11 x 320 mm	Right	400158	1	500636	Left	400707	1	501893
11 x 340 mm	Right	400165	1	500643	Left	400714	1	502012
11 x 360 mm	Right	400172	1	500650	Left	400721	1	502029
11 x 380 mm	Right	400189	1	500667	Left	400738	1	502036
11 x 400 mm	Right	400196	1	500674	Left	400745	1	502043
11 x 420 mm	Right	400202	1	500681	Left	400752	1	502050
11 x 440 mm	Right	400219	1	500698	Left	400769	1	502067
11 x 460 mm	Right	400226	1	500728	Left	400776	1	502074

MULTIAXIAL FEMORAL INTERLOCKING NAIL SET

TRAY NO 4
ORDER INFORMATION 8698791 ...
Multiaxial Femoral Locking Intramedullary Nail

Specification	Side	non-Sterile Titanium	Unit	Sterile Titanium	Side	non-Sterile Titanium	Unit	Sterile Titanium
12 x 320 mm	Right	400233	1	500735	Left	400783	1	502081
12 x 340 mm	Right	400240	1	500742	Left	400790	1	502098
12 x 360 mm	Right	400257	1	500759	Left	400806	1	502210
12 x 380 mm	Right	400264	1	500766	Left	400813	1	502227
12 x 400 mm	Right	400271	1	500773	Left	400820	1	502234
12 x 420 mm	Right	400288	1	500780	Left	400837	1	502241
12 x 440 mm	Right	400295	1	500797	Left	400844	1	502258
12 x 460 mm	Right	400301	1	500919	Left	400851	1	502265


TRAY NO 5
ORDER INFORMATION 8698791 ...
Multiaxial Femoral Locking Intramedullary Nail

Specification	Side	non-Sterile Titanium	Unit	Sterile Titanium	Side	non-Sterile Titanium	Unit	Sterile Titanium
13 x 320 mm	Right	400318	1	500926	Left	400868	1	502272
13 x 340 mm	Right	400325	1	500933	Left	400875	1	502289
13 x 360 mm	Right	400332	1	500940	Left	400882	1	502296
13 x 380 mm	Right	400349	1	500957	Left	400899	1	502883
13 x 400 mm	Right	400356	1	500964	Left	400905	1	502890
13 x 420 mm	Right	400363	1	500971	Left	400912	1	502906
13 x 440 mm	Right	400370	1	500988	Left	400929	1	502913
13 x 460 mm	Right	400387	1	500995	Left	400936	1	502920

MULTIAXIAL FEMORAL INTERLOCKING NAIL SET *

ORDER INFORMATION 8698791 ...

Multiaxial Femoral Locking Intramedullary Nail

Specification	Side	non-Sterile Titanium	Sterile Titanium	Side	non-Sterile Titanium	Sterile Titanium
14 x 320 mm *	Right	400394	501473	Left	400943	502937
14 x 340 mm *	Right	400400	501480	Left	400950	502944
14 x 360 mm *	Right	400417	501497	Left	400967	502951
14 x 380 mm *	Right	400424	501503	Left	400974	502968
14 x 400 mm *	Right	400431	501510	Left	400981	502975
14 x 420 mm *	Right	400448	501527	Left	400998	502982
14 x 440 mm *	Right	400455	501534	Left	401001	502999
14 x 460 mm *	Right	400462	501541	Left	401018	503125

ORDER INFORMATION 8698791 ...

Multiaxial Femoral Locking Intramedullary Nail

Specification	Side	non-Sterile Titanium	Sterile Titanium	Side	non-Sterile Titanium	Sterile Titanium
15 x 320 mm *	Right	400479	501558	Left	401025	503132
15 x 340 mm *	Right	400486	501565	Left	401032	503149
15 x 360 mm *	Right	400493	501572	Left	401049	503156
15 x 380 mm *	Right	400509	501589	Left	401056	503163
15 x 400 mm *	Right	400516	501596	Left	401063	503170
15 x 420 mm *	Right	400523	501602	Left	401070	503187
15 x 440 mm *	Right	400530	501619	Left	401087	503194
15 x 460 mm *	Right	400547	501626	Left	401094	503781

* Optional size

MULTIAXIAL FEMORAL INTERLOCKING NAIL SET

TRAY NO 6
ORDER INFORMATION 8698791 ...
**Multiaxial Femoral Locking Intramedullary Nail
Canullated Prox. Locking Screw**


Length	non-Sterile Titanium	Unit	Sterile Titanium
6,4 x 50 mm *	402206	-	512912
6,4 x 55 mm *	402213	-	512929
6,4 x 60 mm	402220	2	512936
6,4 x 65 mm	402237	2	512943
6,4 x 70 mm	402244	2	512950
6,4 x 75 mm	402251	2	512967
6,4 x 80 mm	402268	2	512974
6,4 x 85 mm	402275	2	512981
6,4 x 90 mm	402282	2	512998
6,4 x 95 mm	402299	2	513162
6,4 x 100 mm	402305	2	513179
6,4 x 105 mm	402312	2	513186
6,4 x 110 mm	402329	2	513193
6,4 x 115 mm *	402336	-	513209
6,4 x 120 mm *	402343	-	513216

**Locking Screw Self Tapping For
Femoral Intramedullary Nail**


Length	non-Sterile Titanium	Unit	Sterile Titanium
4,8 x 30 mm	069331	2	513223
4,8 x 35 mm	069348	2	513230
4,8 x 40 mm	069355	2	513247
4,8 x 45 mm	069362	2	513254
4,8 x 50 mm	069379	2	513261
4,8 x 55 mm	069386	2	513278
4,8 x 60 mm	069393	2	513285
4,8 x 65 mm	069409	2	513292
4,8 x 70 mm	069416	2	513315
4,8 x 75 mm	069423	2	513322
4,8 x 80 mm	069430	2	513339
4,8 x 85 mm *	084129	-	513346
4,8 x 90 mm *	084136	-	513353


Nail Cannulated End Cap For Multiaxial Femoral Nail


Length	non-Sterile Titanium	Unit	Sterile Titanium
0 mm	052814	4	503798
5 mm	052821	1	503804
10 mm	052838	1	503811
15 mm	052845	1	503828

MULTIAXIAL FEMORAL INTERLOCKING NAIL SET

Item No.	Description	Unit
213567	Cannulated AWL	1
216582	Multiaxial Femur Cannulated Reamer, S.Steel, 135° - 4.5*300	1
085317	Drill with Round Shaft, S.Steel, 4.0*300	2
200659	Depth Gauge For Distal Screw, S.Steel, 148	1
216568	Multiaxial Femur Wire Guide, S.Steel, 135° - 2	1
216520	Multiaxial Femur Screw Guide, S.Steel, 135° - 8	1
200413	Drill Guide, S.Steel, 4.0*195	3
200406	Screw Guide, S.Steel, 8.0*180	3
216599	Canullated T Wrench for Locking Screw, 3.5	1
367772	Cannulated T Screw Driver For Proximal Locking Screw, S.Steel, 4.0	1
308393	Puller Shaft For Femoral, Tibial, Humerus T - Wrench, S.Steel	1
216575	Multiaxial Femur Depth Gauge for Proximal Screw, S.Steel, 2.5	1
308386	T - Holder for Set Screw, S.Steel, 3.5	1
457923	T Wrench For Connection Screw, S.Steel, Short, AA5	1



MULTIAXIAL FEMORAL INTERLOCKING NAIL SET



Item No.	Description	Unit
216537	Multiaxial Femur Main Device, Aluminum	1
308362	Main Device and Distal Bar Connection Screw, S.Steel	1
216513	Multiaxial Femur Distal Bar, Aluminum	1
216544	Sliding Hammer Bar, S.Steel	1
458241	Slotted Hammer, 14	1
210726	Guide Wire, S.Steel, 2.0*423	2
210184	Key for Universal Chuck, S.Steel	1
457916	T - Handle For Drill, S.Steel	1
200437	Trocar, Ø8*300	1
216551	Supine Driver For Multiaxial Femur and TPM-PFN, S.Steel	1
218890	Multiaxial Femur Nail Lock, S.Steel	1
200239	Tissue Protector, S.Steel	1
200697	Ruler, S.Steel, 500	1
200642	Elevator Adjustment T Wrench, S.Steel	1



MULTIAXIAL FEMORAL INTERLOCKING NAIL SET



Item No.	Description	Unit
200499	Flexible Reamer, S.Steel, 9	1
200505	Flexible Reamer, S.Steel, 10	1
200512	Flexible Reamer, S.Steel, 11	1
200529	Flexible Reamer, S.Steel, 12	1
200536	Flexible Reamer, S.Steel, 13	1
457909	T - Handle For Flexible Reamer, S.Steel	1
200635	Adaptor For Flexible Reamer, S.Steel	1
130024	T - Reamer, S.Steel, 9	1
130031	T - Reamer, S.Steel, 10	1
130048	T - Reamer, S.Steel, 11	1
130055	T - Reamer, S.Steel, 12	1
130062	T - Reamer, S.Steel, 13	1
200390	Wrench, S.Steel, 16 - 17	1
200468	Tapped Reamer, S.Steel, 13	1



FEMORAL DYNAMIC INTERLOCKING NAIL SET

TITANIUM 90530042002



290410150580280
Aluminium Container
580 x 280 x 150 mm

ORDER INFORMATION 8698791 ...

Femoral Dynamic Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
10 x 280 mm *	052326	-	506485
10 x 300 mm	108641	1	506492
10 x 320 mm	108665	1	506508
10 x 340 mm	108689	1	506515
10 x 360 mm	108702	1	506522
10 x 380 mm	108726	1	506539
10 x 400 mm	108740	1	506546
10 x 420 mm	108764	1	506553
10 x 440 mm	059134	1	506560

Femoral Dynamic Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
11 x 280 mm *	052302	-	506577
11 x 300 mm *	108795	-	506584
11 x 320 mm	108818	1	506591
11 x 340 mm	108832	1	506607
11 x 360 mm	108856	1	506614
11 x 380 mm	108870	1	506621
11 x 400 mm	108894	1	506638
11 x 420 mm	108917	1	506645
11 x 440 mm	059141	1	506652
11 x 460 mm	059158	1	506669

FEMORAL DYNAMIC INTERLOCKING NAIL SET



TRAY NO 2

ORDER INFORMATION 8698791 ...

Femoral Dynamic Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
9 x 280 mm	059066	1	506416
9 x 300 mm	059073	1	506423
9 x 320 mm	059080	1	506430
9 x 340 mm	059097	1	506447
9 x 360 mm	059103	1	506454
9 x 380 mm	059110	1	506461
9 x 400 mm	059127	1	506478

Femoral Dynamic Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
12 x 280 mm *	052333	-	506676
12 x 300 mm *	108948	-	506683
12 x 320 mm	108962	1	506690
12 x 340 mm	108986	1	506706
12 x 360 mm	109006	1	506713
12 x 380 mm	109020	1	506720
12 x 400 mm	109044	1	506737
12 x 420 mm	109068	1	506744
12 x 440 mm	059165	1	506751
12 x 460 mm	059172	1	506768



TRAY NO 3

ORDER INFORMATION 8698791 ...

Femoral Dynamic Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
8 x 280 mm	058991	1	506348
8 x 300 mm	059004	1	506355
8 x 320 mm	059011	1	506362
8 x 340 mm	059028	1	506379
8 x 360 mm	059035	1	506386
8 x 380 mm	059042	1	506393
8 x 400 mm	059059	1	506409

Femoral Dynamic Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
13 x 280 mm *	052319	-	506775
13 x 300 mm *	109099	-	506782
13 x 320 mm	109112	1	506799
13 x 340 mm	109136	1	506805
13 x 360 mm	109150	1	506812
13 x 380 mm	109174	1	506829
13 x 400 mm	109198	1	506836
13 x 420 mm	109211	1	506843
13 x 440 mm	059189	1	506850
13 x 460 mm	059196	1	506867

* Optional size

FEMORAL DYNAMIC INTERLOCKING NAIL SET *
ORDER INFORMATION 8698791 ...

Femoral Dynamic Interlocking Nail			Femoral Dynamic Interlocking Nail		
Specification	non-Sterile Titanium	Sterile Titanium	Specification	non-Sterile Titanium	Sterile Titanium
14 x 300 mm *	109242	506874	14 x 380 mm *	109327	506911
14 x 320 mm *	109266	506881	14 x 400 mm *	109341	506928
14 x 340 mm *	109280	506898	14 x 420 mm *	109365	506935
14 x 360 mm *	109303	506904			

* Optional for all size


ORDER INFORMATION 8698791 ...

Nail End Cap For Femoral Dynamic Nail				
Length	non-Sterile Titanium	Unit	Sterile Titanium	
0 mm	072447	4	506942	
5 mm *	052692	-	506959	
10 mm *	052708	-	506966	
15 mm *	052715	-	506973	

Compression Screw For Femoral Dynamic Nail			
non-Sterile Titanium	Unit	Sterile Titanium	
072423	4	513476	

ORDER INFORMATION 8698791 ...

Locking Screw Self Tapping For Femoral Intramedullary Nail			
Length	non-Sterile Titanium	Unit	Sterile Titanium
4,8x30 mm	069331	2	513223
4,8x35 mm	069348	2	513230
4,8x40 mm	069355	2	513247
4,8x45 mm	069362	2	513254
4,8x50 mm	069379	2	513261
4,8x55 mm	069386	2	513278
4,8x60 mm	069393	2	513285
4,8x65 mm	069409	2	513292
4,8x70 mm	069416	2	513315
4,8x75 mm	069423	2	513322
4,8x80 mm	069430	2	513339
4,8x85 mm *	084129	-	513346
4,8x90 mm *	084136	-	513353


Proximal Locking Screw Self Tapping For Femoral Intramedullary Nail			
Length	non-Sterile Titanium	Unit	Sterile Titanium
5x30 mm	072300	2	513360
5x35 mm	072317	2	513377
5x40 mm	072324	2	513384
5x45 mm	072331	2	513391
5x50 mm	072348	2	513407
5x55 mm	072355	2	513414
5x60 mm	072362	2	513421
5x65 mm	072379	2	513438
5x70 mm	072386	2	513445
5x75 mm	072393	2	513452
5x80 mm	072409	2	513469

FEMORAL DYNAMIC INTERLOCKING NAIL SET

Item No.	Description	Unit
200673	T - Wrench For End Cap, S.Steel, 5.0	1
308393	Puller Shaft For Femoral, Tibial, Humerus T - Wrench, S.Steel	2
085317	Drill with Round Shaft, S.Steel, 4.0*300	2
200666	T - Wrench For Locking Screw With Compression Screw, S.Steel, 3.5	1
200659	Depth Gauge For Distal Screw, S.Steel, 148	1
308379	T - Wrench For End Cap With Threaded, S.Steel, 5.0 (Femoral)	1
200406	Screw Guide, S.Steel, 8.0*180	4
200413	Drill Guide, S.Steel, 4.0*195	4
200437	Trocar, Ø8*300	1
200468	Tapped Reamer, S.Steel, 13	1
200420	Pointed - Cortical Awl, S.Steel	1
308386	T - Holder for Set Screw, S.Steel, 3.5	1




FEMORAL DYNAMIC INTERLOCKING NAIL SET




TRAY NO 5

Item No.	Description	Unit
200482	Flexible Reamer, S.Steel, 8	1
200499	Flexible Reamer, S.Steel, 9	1
200505	Flexible Reamer, S.Steel, 10	1
200512	Flexible Reamer, S.Steel, 11	1
200529	Flexible Reamer, S.Steel, 12	1
200536	Flexible Reamer, S.Steel, 13	1
457909	T - Handle For Flexible Reamer, S.Steel	1
200635	Adaptor For Flexible Reamer, S.Steel	1
130017	T - Reamer, S.Steel, 8	1
130024	T - Reamer, S.Steel, 9	1
130031	T - Reamer, S.Steel, 10	1
130048	T - Reamer, S.Steel, 11	1
130055	T - Reamer, S.Steel, 12	1
130062	T - Reamer, S.Steel, 13	1
200390	Wrench, S.Steel, 16 - 17	1
200680	Guide Wire, S.Steel, 3.0*690	1




FEMORAL DYNAMIC INTERLOCKING NAIL SET



Item No.	Description	Unit
216834	Femoral Dynamic Interlocking Nail Main Device, Alüminum	1
308362	Main Device and Distal Bar Connection Screw, S.Steel	1
216841	Femoral Dynamic Interlocking Nail Distal Bar, Alüminum	1
200611	Supine Driver, S.Steel	1
200604	Sliding Hammer	1
200239	Tissue Protector, S.Steel	1
210184	Key for Universal Chuck, S.Steel	1
457916	T - Handle For Drill, S.Steel	1
200567	Dynamic Femoral Nail Lock, S.Steel	1
200697	Ruler, S.Steel, 500	1
200642	Elevator Adjustment T Wrench, S.Steel	1




FEMORAL INTERLOCKING NAIL SET
TITANIUM 90530012001


290410150580280
Aluminium Container
580 x 280 x 150 mm

ORDER INFORMATION 8698791 ...

Femoral Interlocking Nail				Femoral Interlocking Nail			
Specification	non-Sterile Titanium	Unit	Sterile Titanium	Specification	non-Sterile Titanium	Unit	Sterile Titanium
10 x 300 mm	057222	1	504283	11 x 300 mm *	052340	-	504603
10 x 320 mm	057239	1	504290	11 x 320 mm	057307	1	504610
10 x 340 mm	057246	1	504412	11 x 340 mm	057314	1	504627
10 x 360 mm	057253	1	504429	11 x 360 mm	057321	1	504634
10 x 380 mm	057260	1	504436	11 x 380 mm	057338	1	504641
10 x 400 mm	057277	1	504443	11 x 400 mm	057345	1	505075
10 x 420 mm	057284	1	504450	11 x 420 mm	057352	1	505082
10 x 440 mm	057291	1	504597	11 x 440 mm	057369	1	505099
				11 x 460 mm	057376	1	505105


FEMORAL INTERLOCKING NAIL SET



TRAY NO 2

ORDER INFORMATION 8698791 ...

Femoral Interlocking Nail				Femoral Interlocking Nail			
Specification	non-Sterile Titanium	Unit	Sterile Titanium	Specification	non-Sterile Titanium	Unit	Sterile Titanium
9 x 280 mm	057154	1	504061	12 x 300 mm *	052357	-	505112
9 x 300 mm	057161	1	504078	12 x 320 mm	057383	1	505129
9 x 320 mm	057178	1	504085	12 x 340 mm	057390	1	505136
9 x 340 mm	057185	1	504092	12 x 360 mm	057406	1	505266
9 x 360 mm	057192	1	504245	12 x 380 mm	057413	1	505273
9 x 380 mm	057208	1	504252	12 x 400 mm	057420	1	505280
9 x 400 mm	057215	1	504269	12 x 420 mm	057437	1	505297
9 x 420 mm *	084099	-	504276	12 x 440 mm	057444	1	505723
				12 x 460 mm	057451	1	505730



TRAY NO 3

ORDER INFORMATION 8698791 ...

Femoral Interlocking Nail				Femoral Interlocking Nail			
Specification	non-Sterile Titanium	Unit	Sterile Titanium	Specification	non-Sterile Titanium	Unit	Sterile Titanium
8 x 280 mm	057086	1	503835	13 x 300 mm *	052364	-	505747
8 x 300 mm	057093	1	503842	13 x 320 mm	057468	1	505754
8 x 320 mm	057109	1	503859	13 x 340 mm	057475	1	505761
8 x 340 mm	057116	1	503866	13 x 360 mm	057482	1	505778
8 x 360 mm	057123	1	503873	13 x 380 mm	057499	1	505785
8 x 380 mm	057130	1	503880	13 x 400 mm	057505	1	505792
8 x 400 mm	057147	1	503897	13 x 420 mm	057512	1	505945
				13 x 440 mm	057529	1	505952
				13 x 460 mm	057536	1	505969

* Optional size


FEMORAL INTERLOCKING NAIL SET *
ORDER INFORMATION 8698791 ...

Femoral Interlocking Nail		
Specification	non-Sterile Titanium	Sterile Titanium
14 x 340 mm *	059486	505976
14 x 360 mm *	059493	505983
14 x 380 mm *	059509	505990
14 x 400 mm *	059516	506218
14 x 420 mm *	059523	506225
14 x 440 mm *	059530	506232

ORDER INFORMATION 8698791 ...

Femoral Interlocking Nail		
Specification	non-Sterile Titanium	Sterile Titanium
15 x 340 mm *	059547	506249
15 x 360 mm *	059554	506256
15 x 380 mm *	059561	506263
15 x 400 mm *	059578	506270
15 x 420 mm *	059585	506287
15 x 440 mm *	059592	506294

* Optional for all size


ORDER INFORMATION 8698791 ...
Locking Screw Self Tapping For Femoral Intramedullary Nail

Length	non-Sterile Titanium	Unit	Sterile Titanium
4,8x30 mm	069331	2	513223
4,8x35 mm	069348	2	513230
4,8x40 mm	069355	2	513247
4,8x45 mm	069362	2	513254
4,8x50 mm	069379	2	513261
4,8x55 mm	069386	2	513278
4,8x60 mm	069393	2	513285
4,8x65 mm	069409	2	513292
4,8x70 mm	069416	2	513315
4,8x75 mm	069423	2	513322
4,8x80 mm	069430	2	513339
4,8x85 mm *	084129	-	513346
4,8x90 mm *	084136	-	513353


ORDER INFORMATION 8698791 ...
Nail End Cap For Femoral Nail


Length	non-Sterile Titanium	Unit	Sterile Titanium
0 mm	069454	4	506300
5 mm *	052630	-	506317
10 mm *	409298	-	506324
15 mm *	409304	-	506331

FEMORAL INTERLOCKING NAIL SET

Item No.	Description	Unit
200673	T - Wrench For End Cap, S.Steel, 5.0	1
085317	Drill with Round Shaft, S.Steel, 4.0*300	2
308386	T - Holder for Set Screw, S.Steel, 3.5	1
200659	Depth Gauge For Distal Screw, S.Steel, 148	1
308379	T - Wrench For End Cap With Threaded, S.Steel, 5.0 (Femoral)	2
200406	Screw Guide, S.Steel, 8.0*180	4
200413	Drill Guide, S.Steel, 4.0*195	4
200468	Tapped Reamer, S.Steel, 13	1
200420	Pointed - Cortical Awl, S.Steel	1
308393	Puller Shaft For Femoral, Tibial, Humerus T - Wrench, S.Steel	2




FEMORAL INTERLOCKING NAIL SET



Item No.	Description	Unit
200482	Flexible Reamer, S.Steel, 8	1
200499	Flexible Reamer, S.Steel, 9	1
200505	Flexible Reamer, S.Steel, 10	1
200512	Flexible Reamer, S.Steel, 11	1
200529	Flexible Reamer, S.Steel, 12	1
200536	Flexible Reamer, S.Steel, 13	1
457909	T - Handle For Flexible Reamer, S.Steel	1
200635	Adaptor For Flexible Reamer, S.Steel	1
130017	T - Reamer, S.Steel, 8	1
130024	T - Reamer, S.Steel, 9	1
130031	T - Reamer, S.Steel, 10	1
130048	T - Reamer, S.Steel, 11	1
130055	T - Reamer, S.Steel, 12	1
130062	T - Reamer, S.Steel, 13	1
200390	Wrench, S.Steel, 16 - 17	2
200680	Guide Wire, S.Steel, 3.0*690	2



FEMORAL INTERLOCKING NAIL SET




Item No.	Description	Unit
216810	Femoral Interlocking Nail Main Device, Aluminum	1
308362	Main Device and Distal Bar Connection Screw, S.Steel	1
200611	Supine Driver, S.Steel	1
200604	Sliding Hammer	1
200239	Tissue Protector, S.Steel	1
216827	Femoral Interlocking Nail Distal Bar, Aluminum	1
210184	Key for Universal Chuck, S.Steel	1
457916	T - Handle For Drill, S.Steel	1
200437	Trocar, Ø8*300	1
200376	Femoral Nail Lock, S.Steel	1
200697	Ruler, S.Steel, 500	1
200642	Elevator Adjustment T Wrench, S.Steel	1



SUPRACONDYLAR INTERLOCKING NAIL SET

TITANIUM 90530022001



TRAY NO 1



290410200580280
Aluminium Container
for Instrument
580 x 280 x 200 mm

ORDER INFORMATION 8698791 ...

Supracondylar Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
9 x 200 mm *	052265	-	507567
9 x 240 mm *	052272	-	507574
9 x 280 mm *	052289	-	507581
9 x 320 mm *	052296	-	507598

Supracondylar Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
10 x 200 mm	057826	1	507710
10 x 240 mm	057833	1	507727
10 x 280 mm	057840	1	507734
10 x 320 mm	057857	1	507741

Supracondylar Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
11 x 200 mm	057864	1	507758
11 x 240 mm	057871	1	507765
11 x 280 mm	057888	1	507987
11 x 320 mm	057895	1	507994


Supracondylar Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
12 x 200 mm	057901	1	508007
12 x 240 mm	057918	1	508014
12 x 280 mm	057925	1	508021
12 x 320 mm	057932	1	508038

Supracondylar Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
13 x 200 mm	057949	1	508045
13 x 240 mm	057956	1	508052
13 x 280 mm	057963	1	508069
13 x 320 mm	057970	1	508076

SUPRACONDYLAR INTERLOCKING NAIL SET



ORDER INFORMATION 8698791 ...

Locking Screw Self Tapping For Femoral Intramedullary Nail

Length	non-Sterile Titanium		Sterile Titanium		Length	non-Sterile Titanium		Sterile Titanium	
	Unit	Unit	Unit	Unit		Unit	Unit	Unit	Unit
4,8x30 mm	069331	2	513223		4,8x60 mm	069393	2	513285	
4,8x35 mm	069348	2	513230		4,8x65 mm	069409	2	513292	
4,8x40 mm	069355	2	513247		4,8x70 mm	069416	2	513315	
4,8x45 mm	069362	2	513254		4,8x75 mm	069423	2	513322	
4,8x50 mm	069379	2	513261		4,8x80 mm	069430	2	513339	
4,8x55 mm	069386	2	513278		4,8x85 mm	084129	2	513346	
					4,8x90 mm	084136	2	513353	


Nail End Cap For Femoral Nail

non-Sterile Titanium	Unit
069454	4
Sterile Titanium	Unit
506300	

Item No.	Description	Unit
200673	T - Wrench For End Cap, S.Steel, 5.0	1
085317	Drill with Round Shaft, S.Steel, 4.0*300	2
308386	T - Holder for Set Screw, S.Steel, 3.5	1
200659	Depth Gauge For Distal Screw, S.Steel, 148	1
308379	T - Wrench For End Cap With Threaded, S.Steel, 5.0 (Femoral)	1
200406	Screw Guide, S.Steel, 8.0*180	4
200413	Drill Guide, S.Steel, 4.0*195	4
200468	Tapped Reamer, S.Steel, 13	1
200420	Pointed - Cortical Awl, S.Steel	1
308393	Puller Shaft For Femoral, Tibial, Humerus T - Wrench, S.Steel	1




SUPRACONDYLAR INTERLOCKING NAIL SET



Item No.	Description	Unit
200505	Flexible Reamer, S.Steel, 10	1
200512	Flexible Reamer, S.Steel, 11	1
200529	Flexible Reamer, S.Steel, 12	1
200536	Flexible Reamer, S.Steel, 13	1
457909	T - Handle For Flexible Reamer, S.Steel	1
200635	Adaptor For Flexible Reamer, S.Steel	1
130031	T - Reamer, S.Steel, 10	1
130048	T - Reamer, S.Steel, 11	1
130055	T - Reamer, S.Steel, 12	1
130062	T - Reamer, S.Steel, 13	1
200390	Wrench, S.Steel, 16 - 17	1
200680	Guide Wire, S.Steel, 3.0*690	2



SUPRACONDYLAR INTERLOCKING NAIL SET




Item No.	Description	Unit
441588	Supracondylar Interlocking Nail Main Device, Aluminum	1
308362	Main Device and Distal Bar Connection Screw, S.Steel	1
200611	Supine Driver, S.Steel	1
200604	Sliding Hammer	1
200239	Tissue Protector, S.Steel	1
200246	Supracondylar Nail Distal Bar, Aluminum	1
210184	Key for Universal Chuck, S.Steel	1
457916	T - Handle For Drill, S.Steel	1
200437	Trocar, Ø8*300	1
200376	Femoral Nail Lock, S.Steel	1
200697	Ruler, S.Steel, 500	1
200642	Elevator Adjustment T Wrench, S.Steel	1



TIBIAL DYNAMIC INTERLOCKING NAIL SET

TITANIUM 90531022001



TRAY NO 1



290410150580280
Aluminium Container
580 x 280 x 150 mm

ORDER INFORMATION 8698791 ...


Tibial Dynamic Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
8 x 280 mm	107231	1	510154
8 x 300 mm	107255	1	510161
8 x 320 mm	107279	1	510178
8 x 340 mm	107293	1	510185
8 x 360 mm	107316	1	510192
8 x 380 mm	107330	1	510208
8 x 400 mm	107354	1	510215

Tibial Dynamic Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
9 x 280 mm	107361	1	510369
9 x 300 mm	107385	1	510376
9 x 320 mm	107408	1	510383
9 x 340 mm	107422	1	510390
9 x 360 mm	107446	1	510406
9 x 380 mm	107460	1	510413
9 x 400 mm	107484	1	510420

TIBIAL DYNAMIC INTERLOCKING NAIL SET



TRAY NO 2


ORDER INFORMATION 8698791 ...

Tibial Dynamic Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
10 x 280 mm	107491	1	510437
10 x 300 mm	107514	1	510581
10 x 320 mm	107538	1	510598
10 x 340 mm	107552	1	510604
10 x 360 mm	107576	1	510611
10 x 380 mm	107590	1	510628
10 x 400 mm	107613	1	510635

Tibial Dynamic Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
11 x 280 mm	107620	1	510727
11 x 300 mm	107644	1	510734
11 x 320 mm	107668	1	510741
11 x 340 mm	107682	1	510758
11 x 360 mm	107705	1	510765
11 x 380 mm	107729	1	510772
11 x 400 mm	107743	1	510789



TRAY NO 3

ORDER INFORMATION 8698791 ...

Tibial Dynamic Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
12 x 280 mm	107750	1	510796
12 x 300 mm	107774	1	510963
12 x 320 mm	107798	1	510970
12 x 340 mm	107811	1	510987
12 x 360 mm	107835	1	510994
12 x 380 mm	107859	1	511007
12 x 400 mm	107873	1	511014

Tibial Dynamic Interlocking Nail

Specification	non-Sterile Titanium	Unit	Sterile Titanium
13 x 280 mm	059271	1	511021
13 x 300 mm	059288	1	511038
13 x 320 mm	059295	1	511120
13 x 340 mm	059301	1	511137
13 x 360 mm	059318	1	511144
13 x 380 mm	059325	1	511151
13 x 400 mm	059332	1	511168

TIBIAL DYNAMIC INTERLOCKING NAIL SET

ORDER INFORMATION 8698791 ...

Nail End Cap For Tibial Dynamic Nail				
Length	non-Sterile Titanium	Unit	Sterile Titanium	
0 mm	072706	4	511199	
5 mm *	052722 *	-	511205	
10 mm *	052739 *	-	511489	
15 mm *	052746 *	-	511496	

Compression Screw For Tibial Dynamic Nail			
non-Sterile Titanium	Unit	Sterile Titanium	
072683	4	514787	

ORDER INFORMATION 8698791 ...

Locking Screw Self Tapping For Tibia Nail				
Length	non-Sterile Titanium	Unit	Sterile Titanium	
4,5x30 mm	361916	2	514084	
4,5x35 mm	361923	2	514091	
4,5x40 mm	361930	2	514107	
4,5x45 mm	361947	2	514114	
4,5x50 mm	361954	2	514121	
4,5x55 mm	361961	2	514138	
4,5x60 mm	361978	2	514145	
4,5x65 mm	361985	2	514398	
4,5x70 mm	361992	2	514404	
4,5x75 mm	362005	2	514411	
4,5x80 mm	362012	2	514428	


Proximal Locking Screw Self Tapping For Tibial Dynamic Nail			
Length	non-Sterile Titanium	Unit	Sterile Titanium
4,5x30 mm	072560	2	514435
4,5x35 mm	072577	2	514442
4,5x40 mm	072584	2	514459
4,5x45 mm	072591	2	514701
4,5x50 mm	072607	2	514718
4,5x55 mm	072614	2	514725
4,5x60 mm	072621	2	514732
4,5x65 mm	072638	2	514749
4,5x70 mm	072645	2	514756
4,5x75 mm	072652	2	514763
4,5x80 mm	072669	2	514770

TIBIAL DYNAMIC INTERLOCKING NAIL SET

Item No.	Description	Unit
200673	T - Wrench For End Cap, S.Steel, 5.0	1
085362	Drill with Round Shaft for Nail, S.Steel, 3.5*300	2
200666	T - Wrench For Locking Screw With Compression Screw, S.Steel, 3.5	1
457336	Screw Guide, S.Steel, 8.0*169	4
308386	T - Holder for Set Screw, S.Steel, 3.5	1
200598	Drill Guide, S.Steel, 3.5*190	4
200628	Depth Gauge For Distal Screw, S.Steel, 120	1
200772	Tapped Reamer, S.Steel, 12	1
200420	Pointed - Cortical Awl, S.Steel	1
308393	Puller Shaft For Femoral, Tibial, Humerus T - Wrench, S.Steel	1



TIBIAL DYNAMIC INTERLOCKING NAIL SET




TRAY NO 5

Item No.	Description	Unit
200482	Flexible Reamer, S.Steel, 8	1
200499	Flexible Reamer, S.Steel, 9	1
200505	Flexible Reamer, S.Steel, 10	1
200512	Flexible Reamer, S.Steel, 11	1
200529	Flexible Reamer, S.Steel, 12	1
200536	Flexible Reamer, S.Steel, 13	1
457909	T - Handle For Flexible Reamer, S.Steel	1
200635	Adaptor For Flexible Reamer, S.Steel	1
130017	T - Reamer, S.Steel, 8	1
130024	T - Reamer, S.Steel, 9	1
130031	T - Reamer, S.Steel, 10	1
130048	T - Reamer, S.Steel, 11	1
130055	T - Reamer, S.Steel, 12	1
130062	T - Reamer, S.Steel, 13	1
200390	Wrench, S.Steel, 16 - 17	2
457343	Guide Wire, S.Steel, 2.0*690	1



TIBIAL DYNAMIC INTERLOCKING NAIL SET



Item No.	Description	Unit
216872	Tibial Dynamic Interlocking Nail Main Device, Aluminum	1
308362	Main Device and Distal Bar Connection Screw, S.Steel	1
200611	Supine Driver, S.Steel	1
200604	Sliding Hammer	1
200239	Tissue Protector, S.Steel	1
216889	Tibial Dynamic Interlocking Nail Distal Bar, Aluminum	1
210184	Key for Universal Chuck	1
457916	T - Handle For Drill, S.Steel	1
200437	Trocar, Ø8*300	1
200574	Tibial Dynamic Interlocking Nail Lock, S.Steel	1
200697	Ruler, S.Steel, 500	1
200642	Elevator Adjustment T Wrench, S.Steel	1



TIBIAL INTERLOCKING NAIL SET

TITANIUM 90531012001



TRAY NO 1



290410150580280
Aluminium Container
580 x 280 x 150 mm

ORDER INFORMATION 8698791 ...


Tibial Interlocking Nail				Tibial Interlocking Nail			
Specification	non-Sterile Titanium	Unit	Sterile Titanium	Specification	non-Sterile Titanium	Unit	Sterile Titanium
8 x 280 mm	056201	1	508649	9 x 280 mm	056270	1	508830
8 x 300 mm	056218	1	508656	9 x 300 mm	056287	1	508847
8 x 320 mm	056225	1	508663	9 x 320 mm	056294	1	508977
8 x 340 mm	056232	1	508793	9 x 340 mm	056300	1	508984
8 x 360 mm	056249	1	508809	9 x 360 mm	056317	1	508991
8 x 380 mm	056256	1	508816	9 x 380 mm	056324	1	509004
8 x 400 mm	056263	1	508823	9 x 400 mm	056331	1	509011

TIBIAL INTERLOCKING NAIL SET




ORDER INFORMATION 8698791 ...

Tibial Interlocking Nail				Tibial Interlocking Nail			
Specification	non-Sterile Titanium	Unit	Sterile Titanium	Specification	non-Sterile Titanium	Unit	Sterile Titanium
10 x 280 mm	056348	1	509141	11 x 280 mm	056416	1	509264
10 x 300 mm	056355	1	509158	11 x 300 mm	056423	1	509271
10 x 320 mm	056362	1	509165	11 x 320 mm	056430	1	509288
10 x 340 mm	056379	1	509172	11 x 340 mm	056447	1	509295
10 x 360 mm	056386	1	509189	11 x 360 mm	056454	1	509479
10 x 380 mm	056393	1	509240	11 x 380 mm	056461	1	509486
10 x 400 mm	056409	1	509257	11 x 400 mm	056478	1	509493



ORDER INFORMATION 8698791 ...

Tibial Interlocking Nail				Tibial Interlocking Nail			
Specification	non-Sterile Titanium	Unit	Sterile Titanium	Specification	non-Sterile Titanium	Unit	Sterile Titanium
12 x 280 mm	056485	1	509592	13 x 280 mm	056553	1	509769
12 x 300 mm	056492	1	509608	13 x 300 mm	056560	1	509912
12 x 320 mm	056508	1	509615	13 x 320 mm	056577	1	509929
12 x 340 mm	056515	1	509721	13 x 340 mm	056584	1	509936
12 x 360 mm	056522	1	509738	13 x 360 mm	056591	1	509943
12 x 380 mm	056539	1	509745	13 x 380 mm	056607	1	509950
12 x 400 mm	056546	1	509752	13 x 400 mm	056614	1	509967

TIBIAL INTERLOCKING NAIL SET

TRAY NO 4
ORDER INFORMATION 8698791 ...
Nail End Cap For Tibia Nail


Length	non-Sterile Titanium	Unit	Sterile Titanium
0 mm	069690	4	509974
5 mm *	491330 *	-	509981
10 mm *	491347 *	-	509998
15 mm *	491354 *	-	510147

Locking Screw Self Tapping For Tibia Nail


Length	non-Sterile Titanium	Unit	Sterile Titanium
4,5x30 mm	361916	2	514084
4,5x35 mm	361923	2	514091
4,5x40 mm	361930	2	514107
4,5x45 mm	361947	2	514114
4,5x50 mm	361954	2	514121


Locking Screw Self Tapping For Tibia Nail

Length	non-Sterile Titanium	Unit	Sterile Titanium
4,5x55 mm	361961	2	514138
4,5x60 mm	361978	2	514145
4,5x65 mm	361985	2	514398
4,5x70 mm	361992	2	514404
4,5x75 mm	362005	2	514411
4,5x80 mm	362012	2	514428

Item No.	Description	Unit
200666	T - Wrench For Locking Screw With Compression Screw, S.Steel, 3,5	1
200673	T - Wrench For End Cap, S.Steel, 5,0	1
085362	Drill with Round Shaft, S.Steel, 3,5*300	2
308386	T - Holder for Set Screw, S.Steel, 3,5	1
457336	Screw Guide, S.Steel, 8,0*169	4
200598	Drill Guide, S.Steel, 3,5*190	4
200628	Depth Gauge For Distal Screw, S.Steel, 120	1
200772	Tapped Reamer, S.Steel, 12	1
200420	Pointed - Cortical Awl, S.Steel	1
308393	Puller Shaft For Femoral, Tibial, Humerus T - Wrench, S.Steel	1




TIBIAL INTERLOCKING NAIL SET



Item No.	Description	Unit
200482	Flexible Reamer, S.Steel, 8	1
200499	Flexible Reamer, S.Steel, 9	1
200505	Flexible Reamer, S.Steel, 10	1
200512	Flexible Reamer, S.Steel, 11	1
200529	Flexible Reamer, S.Steel, 12	1
200536	Flexible Reamer, S.Steel, 13	1
457909	T - Handle For Flexible Reamer, S.Steel	1
200635	Adaptor For Flexible Reamer, S.Steel	1
130017	T - Reamer, S.Steel, 8	1
130024	T - Reamer, S.Steel, 9	1
130031	T - Reamer, S.Steel, 10	1
130048	T - Reamer, S.Steel, 11	1
130055	T - Reamer, S.Steel, 12	1
130062	T - Reamer, S.Steel, 13	1
200390	Wrench, S.Steel, 16 - 17	2
457343	Guide Wire, S.Steel, 2.0*690	2



TIBIAL INTERLOCKING NAIL SET




TRAY NO 6

Item No.	Description	Unit
216858	Tibial Interlocking Nail Main Device, Aluminum	1
308362	Main Device and Distal Bar Connection Screw, S.Steel	1
200611	Supine Driver, S.Steel	1
200604	Sliding Hammer	1
200239	Tissue Protector, S.Steel	1
216865	Tibial Interlocking Nail Distal Bar, Aluminum	1
210184	Key for Universal Chuck, S.Steel	1
457916	T - Handle For Drill, S.Steel	1
200437	Trocar, Ø8*300	1
200260	Tibia Interlocking Nail Lock, S.Steel	1
200697	Ruler, S.Steel, 500	1
200642	Elevator Adjustment T Wrench, S.Steel	1



HUMERAL INTERLOCKING NAIL SET
TITANIUM 90525012001




290410200580280
 Aluminium Container
 580 x 280 x 200 mm

ORDER INFORMATION 8698791 ...

Humeral Interlocking Nail, Titanium

Specification	non-Sterile Titanium	Unit	Sterile Titanium
7 x 180 mm	055433	1	511557
7 x 200 mm	055440	1	511564
7 x 220 mm	055457	1	511571
7 x 240 mm	055464	1	511588
7 x 260 mm	055471	1	511595
7 x 280 mm	055488	1	511878
7 x 300 mm	055495	1	511885
8 x 180 mm	055570	1	511892
8 x 200 mm	055587	1	511908
8 x 220 mm	055594	1	511915
8 x 240 mm	055600	1	511922
8 x 260 mm	055617	1	511939
8 x 280 mm	055624	1	511946
8 x 300 mm	055631	1	511953
9 x 180 mm	055716	1	511960
9 x 200 mm	055723	1	511977
9 x 220 mm	055730	1	511984
9 x 240 mm	055747	1	511991
9 x 260 mm	055754	1	512158
9 x 280 mm	055761	1	512165
9 x 300 mm	055778	1	512172

HUMERAL INTERLOCKING NAIL SET



ORDER INFORMATION 8698791 ...



Nail End Cap For Humeral Nail

non-Sterile Titanium	Unit	Sterile Titanium
068228	4	512189

ORDER INFORMATION 8698791 ...



Locking Screw Self Tapping For Humeral Nail

Length	non-Sterile Titanium	Unit	Sterile Titanium
4 x 20 mm	067603	2	515753
4 x 22 mm	067610	2	515760
4 x 24 mm	067627	2	515777
4 x 26 mm	067634	2	515784
4 x 28 mm	067641	2	515791
4 x 30 mm	067658	2	515807
4 x 32 mm	067665	2	515906

Locking Screw Self Tapping For Humeral Nail

Length	non-Sterile Titanium	Unit	Sterile Titanium
4 x 34 mm	067672	2	515913
4 x 36 mm	067689	2	515920
4 x 38 mm	067696	2	515937
4 x 40 mm	067702	2	515944
4 x 42 mm	067719	2	516118
4 x 44 mm	067726	2	516125

Item No.	Description	Unit
210184	Key for Universal Chuck, S.Steel	1
457916	T - Handle For Drill, S.Steel	1
458203	Depth Gauge For Distal Screw, S.Steel, 110	1




HUMERAL INTERLOCKING NAIL SET

Item No.	Description	Unit
200666	T - Wrench For Locking Screw With Compression Screw, S.Steel, 3.5	1
457923	T Wrench For Connection Screw, S.Steel, Short, AA5	1
203445	Tapped Reamer, S.Steel, 11	1
308386	T - Holder for Set Screw, S.Steel, 3.5	1
308393	Puller Shaft For Femoral, Tibial, Humerus T - Wrench, S.Steel	1
203476	Screw Guide, S.Steel, 8.0*159	4
203407	Drill Guide, S.Steel, 2.7*179	4
200420	Pointed - Cortical Awl, S.Steel	1




HUMERAL INTERLOCKING NAIL SET



Item No.	Description	Unit
203360	Humeral Nail Main Device, Aluminum	1
244158	Humeral Nail Main Device and Distal Bar Connection Screw, S.Steel	1
203353	Humeral Nail Distal Bar, Aluminum	1
200611	Supine Driver, S.Steel	1
200604	Sliding Hammer	1
200239	Tissue Protector, S.Steel	1
200437	Trocar, Ø8*300	1
203438	Humeral Nail Lock, S.Steel	1
200642	Elevator Adjustment T Wrench, S.Steel	1



HUMERAL INTERLOCKING NAIL SET



Item No.	Description	Unit
203452	Ruler, S.Steel, 300	1
457602	T - Reamer, S.Steel, 7*350	1
457619	T - Reamer, S.Steel, 8*350	1
457626	T - Reamer, S.Steel, 9*350	1
457633	Flexible Reamer, S.Steel, 7*380	1
457640	Flexible Reamer, S.Steel, 8*380	1
457657	Flexible Reamer, S.Steel, 9*380	1
200390	Wrench, S.Steel, 16 - 17	2
457909	T - Handle For Flexible Reamer, S.Steel	1
200635	Adaptor For Flexible Reamer, S.Steel	1
085270	Drill with Round Shaft, S.Steel, 2.7*300	2
200253	Guide Wire, S.Steel, 2.0*500	2



FIBULA INTERLOCKING NAIL SET
TITANIUM - 90532012001

ORDER INFORMATION 8698791 ...
Fibula Intramedullary Nail

Length	Size	non-Sterile Titanium	Unit	Sterile Titanium
6.5 x 90 mm	Short	052555	1	511502
6.5 x 150 mm	Long	052562	1	511519
6.5 x 200 mm	X-Long	083320	1	511526
6.5 - 3 x 90 mm *	Short	083313	-	511533

ORDER INFORMATION 8698791 ...
Fibula Intramedullary Nail End Cap

non-Sterile Titanium	Unit	Sterile Titanium
052579	4	511540

ORDER INFORMATION 8698791 ...
Fibula Intramedullar Nail Locking Screw


Length	non-Sterile Titanium	Unit	Sterile Titanium
3.5 x 10 mm *	068075	-	514794
3.5 x 12 mm *	068082	-	515043
3.5 x 14 mm	068099	2	515050
3.5 x 16 mm	068105	2	515067
3.5 x 18 mm	068112	2	515074
3.5 x 20 mm	068129	2	515081

Fibula Intramedullar Nail Locking Screw

Length	non-Sterile Titanium	Unit	Sterile Titanium
3.5 x 22 mm	068136	2	515098
3.5 x 24 mm *	068143	-	515104
3.5 x 26 mm	068150	2	515111
3.5 x 28 mm *	068167	-	515128
3.5 x 30 mm	068174	2	515135
3.5 x 35 mm	083146	2	515142
3.5 x 40 mm	083153	2	515395

FIBULA INTERLOCKING NAIL SET

ORDER INFORMATION 8698791 ...

Fibula Intramedullar Nail Compression Screw				Fibula Intramedullar Nail Locking Screw Full Threaded		
Size	non-Sterile Titanium	Unit	Sterile Titanium	Length	non-Sterile Titanium	Sterile Titanium
M3 x 50 mm	083252	1	515494	3.5 x 14 mm *	083160	515401
M3 x 55 mm	083269	1	515548	3.5 x 16 mm *	083177	515418
M3 x 60 mm	083276	1	515555	3.5 x 18 mm *	083184	515425
M3 x 65 mm	083283	1	515562	3.5 x 20 mm *	083191	515432
				3.5 x 22 mm *	083207	515449
				3.5 x 26 mm *	083214	515456
				3.5 x 30 mm *	083221	515463
				3.5 x 35 mm *	083238	515470
				3.5 x 40 mm *	083245	515487

Fibula Intramedullar Nail Compression Nut, M3			
non-Sterile Titanium	Unit	Sterile Titanium	
083290	4	516132	

ORDER INFORMATION 8698791 ...

Fibula Intramedullar Nail Compression Washer for Screw				Fibula Intramedullar Nail Compression Washer for Nut			
Size	non-Sterile Titanium	Unit	Sterile Titanium	Size	non-Sterile Titanium	Unit	Sterile Titanium
14 x 8 mm	083306	4	515579	14 x 10 mm	083511	4	516149

Item No.	Description	Unit
200321	Fibula Nail Main Device, S.Steel	1
080268	1.8 Trocar Point Kirschner Wire, S.Steel, 290	2
457046	Screw Driver for Screw Holder AA2.5	1
219491	Screw Driver for Screw Holder AA3.5	1
457091	Fibular Nail Lock, S.Steel	1



FIBULA INTERLOCKING NAIL SET

Item No.	Description	Unit	
457480	Drill Guide, S.Steel, 2.7*77	1	
457503	Drill Guide, S.Steel, 3.2*68	1	
457510	Drill Guide, S.Steel, 4.0*68	1	
457527	Sleeve Guide, S.Steel, 6.0*65	1	
457534	Sleeve Guide, S.Steel, 6.0*70	1	
457541	Canullated Drill Guide, S.Steel, 6.5	1	
085041	Drill with Round Shaft, S.Steel, 2.7*130	1	
085065	Drill with Round Shaft, S.Steel, 3.2*130	1	
085317	Drill with Round Shaft, S.Steel, 4.0*300	1	
457824	Canullated Drill, S.Steel, 6.5	1	
458173	Fibula Nail Impactor, S.Steel	1	
458197	Fibula Nail Depth Gauge For Distal Screw, S.Steel, 45	1	

Regain
your
mobility

www.tipmed.com



Tipmed Tıbbi Medikal Ürünler San. Tic. Ltd. Şti.

Kemalpaşa Mahallesi 7401/1 Sokak No:16 Bornova / İzmir - TÜRKİYE

T : +90 232 479 33 22

M : +90 543 479 33 22

F : +90 232 479 38 25

info@tipmed.com

marketing@tipmed.com

