



*Regain
your
mobility*



TPM
Revision
KNEE

www.tipmed.com

FEMORAL PS COMPONENT

Co.Cr.Mo Alloyed

Size : 60 - 64 - 66* - 68 - 70* - 72 - 76 - 80 mm (Right & Left)



POSTERIOR FEMORAL AUGMENT

Titanium Alloyed

Size : 60 - 64 - 66* - 68 - 70* - 72 - 76 - 80 mm

Thickness : 5 - 10 mm



Fixation Screw to Femoral Component

DISTAL FEMORAL AUGMENT


Titanium Alloyed

Right Lateral & Left Medial

Right Medial & Left Lateral

Size : 60 - 64 - 66* - 68 - 70* - 72 - 76 - 80 mm

Thickness : 5 - 10 - 15 mm



Fixation Screw to Femoral Component

PATELLA COMPONENT

UHMWPE Alloyed

3 Peg Flat Length :

28mm - 31mm

34mm - 37mm



TPM
Revision
KNEE



EXTENSION STEM

Titanium Alloyed

To be used for both Femur and Tibia

Diameter: 10 mm up to 20 mm with 1 mm increment

Length : 40 - 80 - 120 - 160 mm



TIBIAL AUGMENT

Titanium Alloyed

Right Lateral / Left Medial

Right Medial / Left Lateral

Thickness : 5 - 10 - 15 mm

Size : 60 - 64 - 66* - 68 - 70* - 72 - 76 - 80 mm



Fixation Screw to Tibial Metal Back

OFFSET ADAPTER

Titanium Alloyed

Offset : 0 - 2 - 4 - 6 mm



Offset Stem Adapter Fixation Screw to Tibial Metal Back

Offset Stem Adapter Fixation Screw to Femoral Component

TIBIAL METAL BACK

Co.Cr.Mo Alloyed

Cemented & Uncemented

Size : 60 - 64 - 66* - 68 - 70* - 72 - 76 - 80 mm



UHMWPE POSTERIOR STABILIZED (PS) TIBIAL INSERT

UHMWPE & High Crosslink PE Alloyed

Size : 60 - 64 - 66* - 68 - 70* - 72 - 76 - 80 mm

Diameter: 10 mm up to 24 mm with 2 mm increment



Regain
your
mobility

 **Tipmed**

* On demand